

XISL: X-ray Imaging Software Library

Generated by Doxygen 1.7.6.1

Mon Nov 21 2016 14:43:42

Contents

1	Module Index	1
1.1	Modules	1
2	Class Index	3
2.1	Class List	3
3	Module Documentation	5
3.1	Init Functions for XRpad	5
3.1.1	Function Documentation	6
3.1.1.1	Acquisition_GetVersionInfo	6
3.1.1.2	Acquisition_Set_OnboardOptions	6
3.1.1.3	Acquisition_Set_OnboardOptionsPostOffset	7
3.1.1.4	Acquisition_Set_OnboardOptionsPostOffsetEx	8
3.1.1.5	Acquisition_wpe_ActivateNetworkConfig	9
3.1.1.6	Acquisition_wpe_getAvailableSystems	9
3.1.1.7	Acquisition_wpe_GetNetworkConfigs	10
3.1.1.8	Acquisition_wpe_ReadCameraRegisters	10
3.2	Special XRpad Functions	12
3.2.1	Function Documentation	16
3.2.1.1	Acq_wpe_LoadCorrectionImageToBuffer	16
3.2.1.2	Acq_wpe_SetImageTransferInterface	17
3.2.1.3	Acq_wpe_SystemControl	17
3.2.1.4	Acquisition_AcknowledgeImage	18
3.2.1.5	Acquisition_AckSDCardForceFsck	18
3.2.1.6	Acquisition_AckSDCardForceFsckError	18
3.2.1.7	Acquisition_CloseFile	19

3.2.1.8	Acquisition_CreateFakeShockCriticalLevel	19
3.2.1.9	Acquisition_CreateFakeShockWarningLevel	19
3.2.1.10	Acquisition_DisableEventCallback	20
3.2.1.11	Acquisition_Enable_EMI_Data_Readout	20
3.2.1.12	Acquisition_FactoryResetShock	20
3.2.1.13	Acquisition_FreeFTPFileBuffer	21
3.2.1.14	Acquisition_FTP_CloseSession	21
3.2.1.15	Acquisition_FTP_InitSession	21
3.2.1.16	Acquisition_Get_Current_Voltage	22
3.2.1.17	Acquisition_GetAutoPowerOnLocations	22
3.2.1.18	Acquisition_GetBatteryStatus	22
3.2.1.19	Acquisition_GetChargeMode	23
3.2.1.20	Acquisition_GetChargeModePAcqDesc	23
3.2.1.21	Acquisition_GetFTPFile	24
3.2.1.22	Acquisition_GetGridSensorStatus	24
3.2.1.23	Acquisition_GetLocation	24
3.2.1.24	Acquisition_GetMissedImageCount	25
3.2.1.25	Acquisition_GetNetwork	25
3.2.1.26	Acquisition_GetPowerstate	25
3.2.1.27	Acquisition_GetProvidedEnhancedFeatures	26
3.2.1.28	Acquisition_GetSDCardInfo	26
3.2.1.29	Acquisition_GetSDCardTimeout	27
3.2.1.30	Acquisition_GetTemperatureThresholds	27
3.2.1.31	Acquisition_IdentifyDevice	27
3.2.1.32	Acquisition_IsPreviewImage	28
3.2.1.33	Acquisition_OpenMissedImage	28
3.2.1.34	Acquisition_Resend_All_Messages	28
3.2.1.35	Acquisition_Reset_OnboardOptions	29
3.2.1.36	Acquisition_ResetOnboardShockEvent	29
3.2.1.37	Acquisition_ResetTemperatureTimeout	29
3.2.1.38	Acquisition_Set_FPGA_Power_Mode	30
3.2.1.39	Acquisition_Set_OnboardOffsetImageAcquisition	30
3.2.1.40	Acquisition_Set_OnboardOptionPreview	30
3.2.1.41	Acquisition_SetAutoPowerOnLocations	31

3.2.1.42	Acquisition_SetChargeMode	31
3.2.1.43	Acquisition_SetChargeModePAcqDesc	32
3.2.1.44	Acquisition_SetEventCallback	32
3.2.1.45	Acquisition_SetFakeTemperature	33
3.2.1.46	Acquisition_SetFTPFile	33
3.2.1.47	Acquisition_SetIdleTimeout	33
3.2.1.48	Acquisition_SetNetworkSpeed	34
3.2.1.49	Acquisition_SetPhototimedParams	34
3.2.1.50	Acquisition_SetPrivateKey	35
3.2.1.51	Acquisition_SetSDCardForceFsck	35
3.2.1.52	Acquisition_SetSDCardTimeout	36
3.2.1.53	Acquisition_SetSystemTime	36
3.2.1.54	Acquisition_SetTailTimeforTriggerMode	37
3.2.1.55	Acquisition_SetTemperatureThresholds	37
3.2.1.56	Acquisition_SetTemperatureTimeout	38
3.2.1.57	Acquisition_VerifyGenuineness	38
3.2.1.58	Acquisition_wpe_ChangeNetworkConfig	39
3.2.1.59	Acquisition_wpe_FillDefaultNetworkConfiguration	39
3.2.1.60	Acquisition_wpe_ForceIP	39
3.2.1.61	Acquisition_wpe_GetExamFlag	40
3.2.1.62	Acquisition_wpe_SetMaxOnboardCorrValue	40
3.2.1.63	Acquisition_wpe_SetUniqueImageTag	40
3.3	Common Init functions	42
3.3.1	Function Documentation	43
3.3.1.1	Acquisition_Close	43
3.3.1.2	Acquisition_CloseAll	43
3.3.1.3	Acquisition_EnableLogging	44
3.3.1.4	Acquisition_EnumSensors	44
3.3.1.5	Acquisition_GetCommChannel	45
3.3.1.6	Acquisition_GetConfiguration	45
3.3.1.7	Acquisition_GetHwHeaderInfo	46
3.3.1.8	Acquisition_GetHwHeaderInfoEx	47
3.3.1.9	Acquisition_GetIntTimes	47
3.3.1.10	Acquisition_GetNextSensor	48

3.3.1.11	Acquisition_Init	48
3.3.1.12	Acquisition_SetCallbacksAndMessages	50
3.3.1.13	Acquisition_SetLogLevel	51
3.3.1.14	Acquisition_SetLogOutput	51
3.3.1.15	Acquisition_TogglePerformanceLogging	52
3.4	Init Functions for GbIF	53
3.4.1	Function Documentation	55
3.4.1.1	Acquisition_GbIF_CheckNetworkSpeed	55
3.4.1.2	Acquisition_GbIF_DiscoverDetectors	55
3.4.1.3	Acquisition_GbIF_DiscoveredDetectorByIndex	55
3.4.1.4	Acquisition_GbIF_DiscoveredDetectorCount	56
3.4.1.5	Acquisition_GbIF_ForceIP	56
3.4.1.6	Acquisition_GbIF_GetConnectionSettings	57
3.4.1.7	Acquisition_GbIF_GetDetectorProperties	57
3.4.1.8	Acquisition_GbIF_GetDevice	58
3.4.1.9	Acquisition_GbIF_GetDeviceCnt	58
3.4.1.10	Acquisition_GbIF_GetDeviceList	59
3.4.1.11	Acquisition_GbIF_GetDeviceParams	59
3.4.1.12	Acquisition_GbIF_GetFilterDrvState	59
3.4.1.13	Acquisition_GbIF_GetPacketDelay	60
3.4.1.14	Acquisition_GbIF_GetVersion	60
3.4.1.15	Acquisition_GbIF_Init	61
3.4.1.16	Acquisition_GbIF_SetConnectionSettings	62
3.4.1.17	Acquisition_GbIF_SetDiscoveryTimeout	62
3.4.1.18	Acquisition_GbIF_SetPacketDelay	63
3.4.1.19	Acquisition_GbIF_SetPortRange	63
3.4.1.20	Acquisition_wpe_GetVersion	64
3.5	Acquire and correction functions	65
3.5.1	Function Documentation	68
3.5.1.1	Acquisition_Abort	68
3.5.1.2	Acquisition_AbortCurrentFrame	69
3.5.1.3	Acquisition_Acquire_GainImage	69
3.5.1.4	Acquisition_Acquire_GainImage_Ex	70
3.5.1.5	Acquisition_Acquire_GainImage_Ex_ROI	71

3.5.1.6	Acquisition_Acquire_GainImage_Ex_ROI_PreloadCorr	72
3.5.1.7	Acquisition_Acquire_GainImage_PreloadCorr	72
3.5.1.8	Acquisition_Acquire_Image	73
3.5.1.9	Acquisition_Acquire_Image_Ex	74
3.5.1.10	Acquisition_Acquire_Image_PreloadCorr	75
3.5.1.11	Acquisition_Acquire_OffsetImage	75
3.5.1.12	Acquisition_Acquire_OffsetImage_Ex	76
3.5.1.13	Acquisition_Acquire_OffsetImage_PreloadCorr	77
3.5.1.14	Acquisition_CreateGainMap	77
3.5.1.15	Acquisition_CreateGainMap32	78
3.5.1.16	Acquisition_CreateOnboardPixelMaskFrom16Bit-PixelMask	78
3.5.1.17	Acquisition_CreatePixelMap	79
3.5.1.18	Acquisition_DefineDestBuffers	79
3.5.1.19	Acquisition_DoOffsetCorrection	80
3.5.1.20	Acquisition_DoOffsetGainCorrection	80
3.5.1.21	Acquisition_DoOffsetGainCorrection32	81
3.5.1.22	Acquisition_DoOffsetGainCorrection_Ex	81
3.5.1.23	Acquisition_DoPixelCorrection	82
3.5.1.24	Acquisition_GetCorrData	82
3.5.1.25	Acquisition_GetCorrData_Ex	83
3.5.1.26	Acquisition_GetLatestFrameHeader	83
3.5.1.27	Acquisition_SetCorrData	84
3.5.1.28	Acquisition_SetCorrData_Ex	84
3.6	Setting up the Device	86
3.6.1	Function Documentation	87
3.6.1.1	Acquisition_GetCameraBinningMode	87
3.6.1.2	Acquisition_GetCameraFOVMode	87
3.6.1.3	Acquisition_GetCameraROI	88
3.6.1.4	Acquisition_GetCameraTriggerMode	88
3.6.1.5	Acquisition_GetIpAddress	89
3.6.1.6	Acquisition_ResetFrameCnt	89
3.6.1.7	Acquisition_SetCameraBinningMode	89
3.6.1.8	Acquisition_SetCameraFOVMode	90

3.6.1.9	Acquisition_SetCameraGain	90
3.6.1.10	Acquisition_SetCameraMode	91
3.6.1.11	Acquisition_SetCameraROI	91
3.6.1.12	Acquisition_SetCameraTriggerMode	92
3.6.1.13	Acquisition_SetDACOffset	92
3.6.1.14	Acquisition_SetDACOffsetBinningFPS	93
3.6.1.15	Acquisition_SetFrameSyncMode	93
3.6.1.16	Acquisition_SetFrameSyncTimeMode	94
3.6.1.17	Acquisition_SetTimerSync	95
3.6.1.18	Acquisition_SetTriggerOutSignalOptions	95
3.7	Callback function	97
3.7.1	Function Documentation	97
3.7.1.1	Acquisition_GetAcqData	97
3.7.1.2	Acquisition_GetActFrame	98
3.7.1.3	Acquisition_GetReady	98
3.7.1.4	Acquisition_GetWinHandle	98
3.7.1.5	Acquisition_SetAcqData	99
3.7.1.6	Acquisition_SetReady	99
3.8	Common error handling	101
3.8.1	Function Documentation	101
3.8.1.1	Acquisition_GetErrorCode	101
3.8.1.2	Acquisition_wpe_GetErrorCode	101
3.8.1.3	Acquisition_wpe_GetErrorCodeEx	101
3.9	functions provided by Acq.h	103
3.9.1	Function Documentation	105
3.9.1.1	Acq_wpe_GetSystemInformation	105
3.9.1.2	Acq_WPE_Init	105
3.9.1.3	Acquisition_ActivateServiceMode	105
3.9.1.4	Acquisition_CreateXISFileInMemory	106
3.9.1.5	Acquisition_DeleteFile	106
3.9.1.6	Acquisition_DoOffsetCorrection32	106
3.9.1.7	Acquisition_DoOffsetGainCorrection_Ex32	107
3.9.1.8	Acquisition_GetConnectionStatus	107
3.9.1.9	Acquisition_GetDACOffsetFloorValueFromFlash	107

3.9.1.10	Acquisition_GetDetectorProperties	108
3.9.1.11	Acquisition_GetFileInfo	108
3.9.1.12	Acquisition_GetRotationAngle	108
3.9.1.13	Acquisition_GetTriggerOutStatus	109
3.9.1.14	Acquisition_GetVersion	109
3.9.1.15	Acquisition_GetXISFileBufferSize	109
3.9.1.16	Acquisition_IsAcquiringData	110
3.9.1.17	Acquisition_LoadFile	110
3.9.1.18	Acquisition_LoadXISFileToMemory	111
3.9.1.19	Acquisition_SaveFile	111
3.9.1.20	Acquisition_SaveRawData	111
3.9.1.21	Acquisition_SetConsoleLogging	112
3.9.1.22	Acquisition_SetDACOffsetFloorValueByMode	112
3.9.1.23	Acquisition_SetDACOffsetFloorValueInFlash	113
3.9.1.24	Acquisition_SetDACOffsetFloorValueInFlashInternal	113
3.9.1.25	Acquisition_SetFileLogging	113
3.9.1.26	Acquisition_SetFPGACameraMode	114
3.9.1.27	Acquisition_SetRotationAngle	114
3.9.1.28	Acquisition_wpe_GetVersionNEW	115
3.10	enumerations provided by Acq.h	116
3.10.1	Define Documentation	120
3.10.1.1	HIS_ALL_OK	121
3.10.1.2	HIS_ERROR_ABORT	121
3.10.1.3	HIS_ERROR_ABORTCURRFRAME	121
3.10.1.4	HIS_ERROR_ACKNOWLEDGE_IMAGE	121
3.10.1.5	HIS_ERROR_ACQ	121
3.10.1.6	HIS_ERROR_ACQ_ALREADY_RUNNING	121
3.10.1.7	HIS_ERROR_ACQABORT	121
3.10.1.8	HIS_ERROR_ACQUISITION	121
3.10.1.9	HIS_ERROR_ALREADY_EXISTS	121
3.10.1.10	HIS_ERROR_AVERAGED_LOST	121
3.10.1.11	HIS_ERROR_BAD_SORTING_PARAM	122
3.10.1.12	HIS_ERROR_BOARDINIT	122
3.10.1.13	HIS_ERROR_BUFFERSPACE_NOT_SUFF	122

3.10.1.14 HIS_ERROR_CORRBUFFER_INCOMPATIBLE . . .	122
3.10.1.15 HIS_ERROR_CREATE_MEMORYMAPPING	122
3.10.1.16 HIS_ERROR_CREATE_MUTEX	122
3.10.1.17 HIS_ERROR_CURL	122
3.10.1.18 HIS_ERROR_DESC_NOT_LOCAL	122
3.10.1.19 HIS_ERROR_DOES_NOT_EXIST	122
3.10.1.20 HIS_ERROR_EMI_NOT_SET	122
3.10.1.21 HIS_ERROR_ENABLE_INTERRUPTS	123
3.10.1.22 HIS_ERROR_ENABLE_ONBOARD_GAINOFFSET .	123
3.10.1.23 HIS_ERROR_ENABLE_ONBOARD_MEAN	123
3.10.1.24 HIS_ERROR_ENABLE_ONBOARD_OFFSET	123
3.10.1.25 HIS_ERROR_ENABLE_ONBOARD_PREVIEW . . .	123
3.10.1.26 HIS_ERROR_FRAME_INV	123
3.10.1.27 HIS_ERROR_FUNC_NOTIMPL	123
3.10.1.28 HIS_ERROR_GET_CHARGE_MODE	123
3.10.1.29 HIS_ERROR_GET_NUM_BOARDS	123
3.10.1.30 HIS_ERROR_GET_ONBOARD_OFFSET	123
3.10.1.31 HIS_ERROR_GETHWHEADERINFO	124
3.10.1.32 HIS_ERROR_GETOSVERSION	124
3.10.1.33 HIS_ERROR_HEADER_TIMEOUT	124
3.10.1.34 HIS_ERROR_HW_ALREADY_OPEN_BY_ANOTH- ER_PROCESS	124
3.10.1.35 HIS_ERROR_HW_BOARD_CHANNEL_ALREADY- _USED	124
3.10.1.36 HIS_ERROR_HWHEADER_INV	124
3.10.1.37 HIS_ERROR_ILLEGAL_INDEX	124
3.10.1.38 HIS_ERROR_INVALID_FILENAME	124
3.10.1.39 HIS_ERROR_INVALID_FUNC_CALL	124
3.10.1.40 HIS_ERROR_INVALID_HANDLE	124
3.10.1.41 HIS_ERROR_INVALID_PARAM	125
3.10.1.42 HIS_ERROR_INVALIDACQDESC	125
3.10.1.43 HIS_ERROR_INVALIDBUFFERNR	125
3.10.1.44 HIS_ERROR_LOAD_COORECTIONIMAGETOBUF- FER	125
3.10.1.45 HIS_ERROR_LOADDRIVER	125

3.10.1.46 HIS_ERROR_MEMORY	125
3.10.1.47 HIS_ERROR_MEMORY_MAPPING	125
3.10.1.48 HIS_ERROR_MISSING_VERSION_INFORMATION	125
3.10.1.49 HIS_ERROR_NO_BOARD_IN_SUBNET	125
3.10.1.50 HIS_ERROR_NO_FPGA_ACK	125
3.10.1.51 HIS_ERROR_NOCAMERA	126
3.10.1.52 HIS_ERROR_NODESC_AVAILABLE	126
3.10.1.53 HIS_ERROR_NOT_DISCOVERED	126
3.10.1.54 HIS_ERROR_NOT_INITIALIZED	126
3.10.1.55 HIS_ERROR_NR_OF_BOARDS_CHANGED	126
3.10.1.56 HIS_ERROR_ONBOARDVGFFAILED	126
3.10.1.57 HIS_ERROR_OPEN_FILE	126
3.10.1.58 HIS_ERROR_READ_DATA	126
3.10.1.59 HIS_ERROR_RETRIEVE_ENHANCED_HEADER	126
3.10.1.60 HIS_ERROR_SET_CHARGE_MODE	126
3.10.1.61 HIS_ERROR_SET_IDLE_TIMEOUT	127
3.10.1.62 HIS_ERROR_SET_IMAGE_TAG	127
3.10.1.63 HIS_ERROR_SET_IMAGE_TAG_LENGTH	127
3.10.1.64 HIS_ERROR_SET_ONBOARD_BINNING	127
3.10.1.65 HIS_ERROR_SET_PROC_SCRIPT	127
3.10.1.66 HIS_ERROR_SETBAUDRATE	127
3.10.1.67 HIS_ERROR_SETCAMERAMODE	127
3.10.1.68 HIS_ERROR_SETDISCOVERYTIMEOUT	127
3.10.1.69 HIS_ERROR_SETEXAMFLAG	127
3.10.1.70 HIS_ERROR_SETFRMSYNC	127
3.10.1.71 HIS_ERROR_SETFRMSYNCMODE	128
3.10.1.72 HIS_ERROR_SETLINETRIG_MODE	128
3.10.1.73 HIS_ERROR_SETTIMERSYNC	128
3.10.1.74 HIS_ERROR_SLOW_SYSTEM	128
3.10.1.75 HIS_ERROR_TIMEOUT	128
3.10.1.76 HIS_ERROR_UNABLE_TO_ACCESS_DETECTOR- _FLASH	128
3.10.1.77 HIS_ERROR_UNABLE_TO_CLOSE_BOARD	128
3.10.1.78 HIS_ERROR_UNABLE_TO_OPEN_BOARD	128

3.10.1.79 HIS_ERROR_UNKNOWN_IP_MAC_NAME	128
3.10.1.80 HIS_ERROR_VXD_REGISTER_DMA_ADDRESS	128
3.10.1.81 HIS_ERROR_VXD_REGISTER_IRQ	129
3.10.1.82 HIS_ERROR_VXD_REGISTER_STAT_ADDR	129
3.10.1.83 HIS_ERROR_VXD_REGISTER_STATADR	129
3.10.1.84 HIS_ERROR_VXD_UNMASK_IRQ	129
3.10.1.85 HIS_ERROR_VXDGETDMAADR	129
3.10.1.86 HIS_ERROR_VXDNOTFOUND	129
3.10.1.87 HIS_ERROR_VXDNOTOPEN	129
3.10.1.88 HIS_ERROR_VXDUNKNOWNERROR	129
3.10.1.89 HIS_ERROR_WRITE_DATA	129
3.10.1.90 HIS_ERROR_WRONG_CAMERA_MODE	129
3.10.1.91 HIS_ERROR_WRONGBOARDSELECT	130
3.10.1.92 HIS_ERROR_XRPD_CONNECT	130
3.10.1.93 HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EV- ENT	130
3.10.1.94 HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EV- ENT_CRIT	130
3.10.1.95 HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EV- ENT_WARN	130
3.10.1.96 HIS_ERROR_XRPD_FACTORY_RESET_SHOCK_- EVENT	130
3.10.1.97 HIS_ERROR_XRPD_GET_AUTOPOWERONLOCA- TIONS	130
3.10.1.98 HIS_ERROR_XRPD_GET_CURRENT_VOLTAGE	130
3.10.1.99 HIS_ERROR_XRPD_GET_SDCARD_INFO	130
3.10.1.100 HIS_ERROR_XRPD_GET_SDCARD_TIMEOUT	130
3.10.1.101 HIS_ERROR_XRPD_GET_TEMPERATURE_THR- ESHOLDS	131
3.10.1.102 HIS_ERROR_XRPD_NO_EVENT_INTERFACE	131
3.10.1.103 HIS_ERROR_XRPD_NO_EVENTCALLBACK_DEF- INED	131
3.10.1.104 HIS_ERROR_XRPD_NO_LOCATION	131
3.10.1.105 HIS_ERROR_XRPD_NO_NETWORK	131
3.10.1.106 HIS_ERROR_XRPD_NOT_CONNECTED	131
3.10.1.107 HIS_ERROR_XRPD_REQUEST_POWERSTATE	131

3.10.1.108	HIS_ERROR_XRPD_RESEND_ALL_MSG	131
3.10.1.109	HIS_ERROR_XRPD_RESET_SHOCK	131
3.10.1.110	HIS_ERROR_XRPD_RESET_TEMPERATURE_TIMEOUT	131
3.10.1.111	HIS_ERROR_XRPD_SDCARDPERFORMANCE	132
3.10.1.112	HIS_ERROR_XRPD_SESSION_ERROR	132
3.10.1.113	HIS_ERROR_XRPD_SET_AUTOPOWERONLOCATIONS	132
3.10.1.114	HIS_ERROR_XRPD_SET_CPUFREQ_GOVERNOR	132
3.10.1.115	HIS_ERROR_XRPD_SET_DATE_TIME	132
3.10.1.116	HIS_ERROR_XRPD_SET_EVENT	132
3.10.1.117	HIS_ERROR_XRPD_SET_FORCE_FSCK	132
3.10.1.118	HIS_ERROR_XRPD_SET_NETWORK	132
3.10.1.119	HIS_ERROR_XRPD_SET_PRIVATE_KEY	132
3.10.1.120	HIS_ERROR_XRPD_SET_SDCARD_TIMEOUT	132
3.10.1.121	HIS_ERROR_XRPD_SET_TEMP_FAKE_MODE	133
3.10.1.122	HIS_ERROR_XRPD_SET_TEMPERATURE_THRESHOLDS	133
3.10.1.123	HIS_ERROR_XRPD_SET_TEMPERATURE_TIMEOUT	133
3.10.1.124	HIS_ERROR_XRPD_VERIFY_GENUINENESS	133
3.10.2	Typedef Documentation	133
3.10.2.1	ACQDESCPOS	133
3.10.2.2	HACQDESC	133
3.10.2.3	OnboardBinningMode	133
3.10.2.4	ProcScriptOperation	133
3.10.2.5	XIS_Acquisition_Event	133
3.10.2.6	XIS_Battery_Event	133
3.10.2.7	XIS_Detector_Event	133
3.10.2.8	XIS_Detector_TRIGOUT_SignalMode	133
3.10.2.9	XIS_DetectorTriggerMode	133
3.10.2.10	XIS_Event	134
3.10.2.11	XIS_FileType	134
3.10.2.12	XIS_Library_Event	134
3.10.2.13	XIS_Sensor_Event	134

3.10.2.14	XislFileEntryType	134
3.10.2.15	XislFileHandle	134
3.10.2.16	XislFileStorageLocation	134
3.10.2.17	XislFtpSession	134
3.10.2.18	XRpad_BatteryHealth	134
3.10.2.19	XRpad_BatteryPresence	134
3.10.2.20	XRpad_BatteryStatus	134
3.10.2.21	XRpad_ChargeMode	134
3.10.2.22	XRpad_DataInterfaceControlEnum	134
3.10.2.23	XRpad_ShockEvent	134
3.10.2.24	XRpad_ShockSensorReport	134
3.10.2.25	XRpad_TempSensor	134
3.10.2.26	XRpad_TempSensorReport	134
3.10.2.27	XRpad_VersionInfo	134
3.10.3	Enumeration Type Documentation	134
3.10.3.1	OnboardBinningMode	134
3.10.3.2	ProcScriptOperation	135
3.10.3.3	XIS_Acquisition_Event	135
3.10.3.4	XIS_Battery_Event	135
3.10.3.5	XIS_Detector_Event	135
3.10.3.6	XIS_Detector_TRIGOUT_SignalMode	136
3.10.3.7	XIS_DetectorTriggerMode	136
3.10.3.8	XIS_Event	136
3.10.3.9	XIS_FileType	137
3.10.3.10	XIS_Library_Event	137
3.10.3.11	XIS_Sensor_Event	137
3.10.3.12	XislFileEntryType	137
3.10.3.13	XislFileStorageLocation	138
3.10.3.14	XislLoggingLevels	138
3.10.3.15	XRpad_BatteryHealth	138
3.10.3.16	XRpad_BatteryPresence	138
3.10.3.17	XRpad_ChargeMode	139
3.10.3.18	XRpad_DataInterfaceControlEnum	139
3.10.3.19	XRpad_SystemControlEnum	139

4	Class Documentation	141
4.1	CHwHeaderInfo Struct Reference	141
4.1.1	Member Data Documentation	142
4.1.1.1	bAddRow	142
4.1.1.2	bAddRow	142
4.1.1.3	bPwrSave	142
4.1.1.4	bPwrSave	142
4.1.1.5	bSyncMode	142
4.1.1.6	bSyncMode	142
4.1.1.7	dwAccess	142
4.1.1.8	dwAccess	142
4.1.1.9	dwAcqMode	142
4.1.1.10	dwAcqMode	142
4.1.1.11	dwBias	142
4.1.1.12	dwBias	143
4.1.1.13	dwDataSorting	143
4.1.1.14	dwDataSorting	143
4.1.1.15	dwDataType	143
4.1.1.16	dwDataType	143
4.1.1.17	dwFrmFillRowIntervalls	143
4.1.1.18	dwFrmFillRowIntervalls	143
4.1.1.19	dwFrmNrRows	143
4.1.1.20	dwFrmNrRows	143
4.1.1.21	dwFrmRowType	143
4.1.1.22	dwFrmRowType	143
4.1.1.23	dwGain	143
4.1.1.24	dwGain	143
4.1.1.25	dwHeaderID	143
4.1.1.26	dwHeaderID	143
4.1.1.27	dwLeakRows	143
4.1.1.28	dwLeakRows	143
4.1.1.29	dwNrColumns	143
4.1.1.30	dwNrColumns	143
4.1.1.31	dwNrOfFillingRows	143

4.1.1.32	dwNrOfFillingRows	143
4.1.1.33	dwNrRows	143
4.1.1.34	dwNrRows	143
4.1.1.35	dwOffset	143
4.1.1.36	dwOffset	144
4.1.1.37	dwPROMID	144
4.1.1.38	dwPROMID	144
4.1.1.39	dwTiming	144
4.1.1.40	dwTiming	144
4.1.1.41	dwZoomBRColumn	144
4.1.1.42	dwZoomBRColumn	144
4.1.1.43	dwZoomBRRow	144
4.1.1.44	dwZoomBRRow	144
4.1.1.45	dwZoomULColumn	144
4.1.1.46	dwZoomULColumn	144
4.1.1.47	dwZoomULRow	144
4.1.1.48	dwZoomULRow	144
4.2	CHwHeaderInfoEx Struct Reference	144
4.2.1	Member Data Documentation	146
4.2.1.1	wAccess	146
4.2.1.2	wAccess	146
4.2.1.3	wBias	146
4.2.1.4	wBias	146
4.2.1.5	wBinningMode	146
4.2.1.6	wBinningMode	146
4.2.1.7	wCameratype	146
4.2.1.8	wCameratype	146
4.2.1.9	wClock	146
4.2.1.10	wClock	146
4.2.1.11	wCommand1	146
4.2.1.12	wCommand1	146
4.2.1.13	wCommand2	146
4.2.1.14	wCommand2	146
4.2.1.15	wCommand3	146

4.2.1.16	wCommand3	146
4.2.1.17	wCommand4	146
4.2.1.18	wCommand4	146
4.2.1.19	wDataSorting	146
4.2.1.20	wDataSorting	146
4.2.1.21	wDummy	146
4.2.1.22	wDummy	146
4.2.1.23	wFrameCnt	146
4.2.1.24	wFrameCnt	147
4.2.1.25	wFrmNrRows	147
4.2.1.26	wFrmNrRows	147
4.2.1.27	wFrmRowType	147
4.2.1.28	wFrmRowType	147
4.2.1.29	wGain	147
4.2.1.30	wGain	147
4.2.1.31	wHeaderID	147
4.2.1.32	wHeaderID	147
4.2.1.33	wLeakRows	147
4.2.1.34	wLeakRows	147
4.2.1.35	wNrColumns	147
4.2.1.36	wNrColumns	147
4.2.1.37	wNrRows	147
4.2.1.38	wNrRows	147
4.2.1.39	wPROMID	147
4.2.1.40	wPROMID	147
4.2.1.41	wRealInttime_microSec	147
4.2.1.42	wRealInttime_microSec	147
4.2.1.43	wRealInttime_milliSec	147
4.2.1.44	wRealInttime_milliSec	147
4.2.1.45	wResolutionX	147
4.2.1.46	wResolutionX	147
4.2.1.47	wResolutionY	147
4.2.1.48	wResolutionY	148
4.2.1.49	wRowTime	148

4.2.1.50	wRowTime	148
4.2.1.51	wStatus	148
4.2.1.52	wStatus	148
4.2.1.53	wTiming	148
4.2.1.54	wTiming	148
4.2.1.55	wUgComp	148
4.2.1.56	wUgComp	148
4.2.1.57	wZoomBRColumn	148
4.2.1.58	wZoomBRColumn	148
4.2.1.59	wZoomBRRow	148
4.2.1.60	wZoomBRRow	148
4.2.1.61	wZoomULColumn	148
4.2.1.62	wZoomULColumn	148
4.2.1.63	wZoomULRow	148
4.2.1.64	wZoomULRow	148
4.3	DETECTOR_BATTERY Struct Reference	148
4.3.1	Member Data Documentation	149
4.3.1.1	capacity	149
4.3.1.2	charge	149
4.3.1.3	current	149
4.3.1.4	cycle_count	149
4.3.1.5	energy	149
4.3.1.6	serial_no	149
4.3.1.7	status	149
4.3.1.8	temperature	149
4.3.1.9	voltage	149
4.4	DETECTOR_CURRENT_VOLTAGE Struct Reference	149
4.4.1	Member Data Documentation	149
4.4.1.1	imA1	149
4.4.1.2	imA2	149
4.4.1.3	imA3	149
4.4.1.4	iV1	149
4.4.1.5	iV2	149
4.4.1.6	iV3	149

4.5	deviceInfo Struct Reference	150
4.5.1	Detailed Description	150
4.5.2	Member Data Documentation	150
4.5.2.1	device_version	150
4.5.2.2	manufacturer_name	150
4.5.2.3	manufacturer_specific	150
4.5.2.4	model_name	150
4.5.2.5	serial_number	151
4.5.2.6	spec_version	151
4.5.2.7	user_name	151
4.6	discoveryReply Struct Reference	151
4.6.1	Detailed Description	151
4.6.2	Member Data Documentation	151
4.6.2.1	deviceInfo	151
4.6.2.2	gvcp_ip	151
4.6.2.3	lanInfo	151
4.6.2.4	wlanInfo	152
4.7	discoveryReplyEx Struct Reference	152
4.7.1	Detailed Description	152
4.7.2	Member Data Documentation	152
4.7.2.1	deviceInfo	152
4.7.2.2	gvcp_ip	152
4.7.2.3	lanInfo	152
4.7.2.4	messageCount	152
4.7.2.5	messages	153
4.7.2.6	wlanInfo	153
4.8	discoveryReplyMsg Struct Reference	153
4.8.1	Detailed Description	153
4.8.2	Member Data Documentation	153
4.8.2.1	local_ip	153
4.8.2.2	remote_ip	153
4.9	EPC_REGISTER Struct Reference	153
4.9.1	Member Data Documentation	154
4.9.1.1	active_network_config	154

4.9.1.2	battery	154
4.9.1.3	channel	154
4.9.1.4	exam_flag	154
4.9.1.5	lan_status_register	154
4.9.1.6	power_state	154
4.9.1.7	rtc_value	154
4.9.1.8	sdcard_state	154
4.9.1.9	sdcard_usage	154
4.9.1.10	signal_strength	154
4.9.1.11	spartan_id	154
4.9.1.12	spartan_register	154
4.9.1.13	temperature_error_level	154
4.9.1.14	temperature_value	154
4.9.1.15	temperature_warning_level	154
4.9.1.16	version	154
4.9.1.17	wlan_status_register	154
4.10	FPGAType Struct Reference	154
4.10.1	Member Data Documentation	155
4.10.1.1	wTiming	155
4.10.1.2	wValue0	155
4.10.1.3	wValue1	155
4.10.1.4	wValue2	155
4.10.1.5	wValue3	155
4.10.1.6	wValue4	155
4.10.1.7	wValue5	155
4.10.1.8	wValue6	155
4.11	GBIF_Detector_Properties Struct Reference	155
4.11.1	Member Data Documentation	155
4.11.1.1	cDetectorType	155
4.11.1.2	cDeviceIdentifier	155
4.11.1.3	cDummy	155
4.11.1.4	cManufacturingDate	155
4.11.1.5	cPlaceOfManufacture	156
4.11.1.6	cUniqueDeviceIdentifier	156

4.12 GBIF_DEVICE_PARAM Struct Reference	156
4.12.1 Member Data Documentation	156
4.12.1.1 cDeviceName	156
4.12.1.2 cDeviceName	156
4.12.1.3 cGBIFFirmwareVersion	156
4.12.1.4 cGBIFFirmwareVersion	156
4.12.1.5 cManufacturerName	156
4.12.1.6 cManufacturerName	156
4.12.1.7 cModelName	156
4.12.1.8 cModelName	157
4.12.1.9 dwIPCurrentBootOptions	157
4.12.1.10 dwIPCurrentBootOptions	157
4.12.1.11 ucAdapterIP	157
4.12.1.12 ucAdapterMask	157
4.12.1.13 ucGateway	157
4.12.1.14 ucIP	157
4.12.1.15 ucMacAddress	157
4.12.1.16 ucSubnetMask	157
4.13 networkAdapterConfiguration Struct Reference	157
4.13.1 Detailed Description	158
4.13.2 Member Data Documentation	158
4.13.2.1 bridged	158
4.13.2.2 dns	158
4.13.2.3 enabled	158
4.13.2.4 gateway	158
4.13.2.5 hw_accel	158
4.13.2.6 ifname	158
4.13.2.7 ipaddr	158
4.13.2.8 macaddr	158
4.13.2.9 netmask	158
4.13.2.10 not_used	159
4.13.2.11 proto	159
4.14 networkConfiguration Struct Reference	159
4.14.1 Detailed Description	159

4.14.2	Member Data Documentation	159
4.14.2.1	gbif_enabled	159
4.14.2.2	hostname	160
4.14.2.3	lan	160
4.14.2.4	name	160
4.14.2.5	notUsed	160
4.14.2.6	path	160
4.14.2.7	readonly	160
4.14.2.8	sshd_enabled	160
4.14.2.9	wifi	160
4.14.2.10	wlan	160
4.15	networkInfo Struct Reference	160
4.15.1	Detailed Description	161
4.15.2	Member Data Documentation	161
4.15.2.1	broadcast	161
4.15.2.2	ip	161
4.15.2.3	mac	161
4.15.2.4	mask	161
4.16	RTC_STRUCT Struct Reference	161
4.16.1	Member Data Documentation	162
4.16.1.1	day	162
4.16.1.2	hour	162
4.16.1.3	minute	162
4.16.1.4	month	162
4.16.1.5	second	162
4.16.1.6	year	162
4.17	wifiConfiguration Struct Reference	162
4.17.1	Detailed Description	162
4.17.2	Member Data Documentation	162
4.17.2.1	agmode	162
4.17.2.2	channel	162
4.17.2.3	description	163
4.17.2.4	mode	163
4.17.2.5	ssid	163

4.18	wifiConfigurationEx Struct Reference	163
4.18.1	Detailed Description	163
4.18.2	Member Data Documentation	163
4.18.2.1	agmode	163
4.18.2.2	channel	164
4.18.2.3	description	164
4.18.2.4	mode	164
4.18.2.5	passphrase	164
4.18.2.6	scan_ssid	164
4.18.2.7	ssid	164
4.19	WinHeaderType Struct Reference	164
4.19.1	Member Data Documentation	165
4.19.1.1	BRX	165
4.19.1.2	BRX	165
4.19.1.3	BRY	165
4.19.1.4	BRY	165
4.19.1.5	Correction	165
4.19.1.6	Correction	165
4.19.1.7	FileSize	165
4.19.1.8	FileSize	165
4.19.1.9	FileType	165
4.19.1.10	FileType	165
4.19.1.11	HeaderSize	165
4.19.1.12	HeaderSize	165
4.19.1.13	HeaderVersion	165
4.19.1.14	HeaderVersion	165
4.19.1.15	ImageHeaderSize	165
4.19.1.16	ImageHeaderSize	165
4.19.1.17	IntegrationTime	166
4.19.1.18	NrOfFrames	166
4.19.1.19	NrOfFrames	166
4.19.1.20	TypeOfNumbers	166
4.19.1.21	TypeOfNumbers	166
4.19.1.22	ULX	166

4.19.1.23 ULX	166
4.19.1.24 ULY	166
4.19.1.25 ULY	166
4.19.1.26 x	166
4.19.1.27 x	166
4.20 WinHeaderType101 Struct Reference	166
4.20.1 Member Data Documentation	166
4.20.1.1 BRX	166
4.20.1.2 BRY	167
4.20.1.3 Correction	167
4.20.1.4 FileSize	167
4.20.1.5 FileType	167
4.20.1.6 HeaderSize	167
4.20.1.7 HeaderVersion	167
4.20.1.8 ImageHeaderSize	167
4.20.1.9 IntegrationTime	167
4.20.1.10 NrOfFrames	167
4.20.1.11 TypeOfNumbers	167
4.20.1.12 ULX	167
4.20.1.13 ULY	167
4.20.1.14 wMedianValue	167
4.20.1.15 x	167
4.21 WinImageHeaderType Struct Reference	167
4.21.1 Member Data Documentation	167
4.21.1.1 dwPROMID	167
4.21.1.2 dwPROMID	168
4.21.1.3 fAmpere	168
4.21.1.4 fKVolt	168
4.21.1.5 n_avframes	168
4.21.1.6 n_avframes	168
4.21.1.7 strPrefilter	168
4.21.1.8 strProject	168
4.21.1.9 strSystemused	168
4.22 XislLoggingErrorHandler Class Reference	168

4.22.1	Constructor & Destructor Documentation	168
4.22.1.1	XislLoggingErrorHandler	168
4.22.1.2	~XislLoggingErrorHandler	168
4.22.2	Member Function Documentation	168
4.22.2.1	error	168
4.22.2.2	reset	168
4.23	XRpad_BatteryStatus Struct Reference	168
4.23.1	Member Data Documentation	169
4.23.1.1	authenticated	169
4.23.1.2	charge_mode	169
4.23.1.3	charge_state	169
4.23.1.4	cycle_count	169
4.23.1.5	design_capacity	169
4.23.1.6	health	169
4.23.1.7	presence	169
4.23.1.8	remaining_capacity	169
4.23.1.9	temperature	169
4.24	XRpad_ShockEvent Struct Reference	169
4.24.1	Member Data Documentation	169
4.24.1.1	critical_sensor1	169
4.24.1.2	critical_sensor2	169
4.24.1.3	critical_sensor3	170
4.24.1.4	timestamp	170
4.24.1.5	warning_sensor1	170
4.24.1.6	warning_sensor2	170
4.24.1.7	warning_sensor3	170
4.25	XRpad_ShockSensorReport Struct Reference	170
4.25.1	Member Data Documentation	170
4.25.1.1	largest	170
4.25.1.2	latest	170
4.26	XRpad_TempSensor Struct Reference	170
4.26.1	Member Data Documentation	170
4.26.1.1	index	170
4.26.1.2	is_virtual	170

4.26.1.3	name	170
4.26.1.4	temperature	170
4.26.1.5	warn_level	170
4.27	XRpad_TempSensorReport Struct Reference	171
4.27.1	Member Data Documentation	171
4.27.1.1	sensor_count	171
4.27.1.2	sensors	171
4.27.1.3	shutdown_time	171
4.27.1.4	system_warn_level	171
4.28	XRpad_VersionInfo Struct Reference	171
4.28.1	Member Data Documentation	171
4.28.1.1	hwdriver	171
4.28.1.2	linux_kernel	171
4.28.1.3	mcp_firmware	171
4.28.1.4	pld_firmware	171
4.28.1.5	software	172
4.28.1.6	spartan_firmware	172
4.28.1.7	subversion	172
4.28.1.8	wlan	172
4.28.1.9	xrpd	172
4.28.1.10	zynq_firmware	172

Chapter 1

Module Index

1.1 Modules

Here is a list of all modules:

Acquire and correction functions	65
Callback function	97
Common Init functions	42
Common error handling	101
Init Functions for GbIF	53
Init Functions for XRpad	5
Setting up the Device	86
Special XRpad Functions	12
enumerations provided by Acq.h	116
functions provided by Acq.h	103

Chapter 2

Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

CHwHeaderInfo	141
CHwHeaderInfoEx	144
DETECTOR_BATTERY	148
DETECTOR_CURRENT_VOLTAGE	149
deviceInfo	150
discoveryReply	151
discoveryReplyEx	152
discoveryReplyMsg	153
EPC_REGISTER	153
FPGAType	154
GBIF_Detector_Properties	155
GBIF_DEVICE_PARAM	156
networkAdapterConfiguration	157
networkConfiguration	159
networkInfo	160
RTC_STRUCT	161
wifiConfiguration	162
wifiConfigurationEx	163
WinHeaderType	164
WinHeaderType101	166
WinImageHeaderType	167
XislLoggingErrorHandler	168
XRpad_BatteryStatus	168
XRpad_ShockEvent	169
XRpad_ShockSensorReport	170
XRpad_TempSensor	170
XRpad_TempSensorReport	171
XRpad_VersionInfo	171

Chapter 3

Module Documentation

3.1 Init Functions for XRpad

Functions

- HIS_RETURN [Acquisition_GetVersionInfo](#) (HACQDESC hAcqDesc, [XRpad_VersionInfo](#) *versionInfo)
Acquisition_GetVersionInfo Retrieves version information of the detector.
- HIS_RETURN [Acquisition_Set_OnboardOptions](#) (HACQDESC hAcqDesc, BOOL bEnableAck, BOOL bOffset, BOOL bGain, BOOL bPixel)
Activate image acquisition options for onboard offset, gain, mean correction and acknowledgement of frames.
- HIS_RETURN [Acquisition_Set_OnboardOptionsPostOffset](#) (HACQDESC hAcqDesc, BOOL bNoOnboardCorr, BOOL bSendPreviewFrist, BOOL bSendFULLFirst, BOOL bEnableAckFirst, BOOL bEnableAckSecond, BOOL bEnableOffsetFirst, BOOL bEnablePostOffsetCorr, BOOL bGain, BOOL bPixel)
This function defines the image acquisition options for PhotoTimed mode including preview on/off offset/gain/pixel correction on/off and which image will be send to the client.
- HIS_RETURN [Acquisition_Set_OnboardOptionsPostOffsetEx](#) (HACQDESC hAcqDesc, BOOL bNoOnboardCorr, BOOL bSendPreviewFrist, BOOL bSendFULLFirst, BOOL bEnableAckFirst, BOOL bEnableAckSecond, BOOL bEnableOffsetFirst, BOOL bEnablePostOffsetCorr, BOOL bGain, BOOL bPixel, BOOL bStoreOffsetToSD)
This function can be used to verify the pgototimed mode. It defines the image acquisition options for PhotoTimed mode including preview on/off offset/gain/pixel correction on/off and which image will be send to the client.
- HIS_RETURN [Acquisition_wpe_ActivateNetworkConfig](#) (const char *ipAddress, int configIndex)
This function activates the selected network configuration of the detector.
- HIS_RETURN [Acquisition_wpe_getAvailableSystems](#) (struct [discoveryReply](#) *reply, int *numDevices, int timeout, int port)

This function sends a network broadcast and retrieves all available XRpad detectors in the Network.

- HIS_RETURN [Acquisition_wpe_GetNetworkConfigs](#) (const char *ipAddress, struct [networkConfiguration](#) *configs, int *arrayLength, int *activeConfig)

This function retrieves the Network settings of a specific Detector defined by its ip address.

- HIS_RETURN [Acquisition_wpe_ReadCameraRegisters](#) (const char *ipAddress, unsigned long *buffer)

This function retrieves the Status register of a specific Detector defined by its ip address.

3.1.1 Function Documentation

3.1.1.1 HIS_RETURN Acquisition_GetVersionInfo (HACQDESC hAcqDesc, XRpad_VersionInfo * versionInfo)

Acquisition_GetVersionInfo Retrieves version information of the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>versionInfo</i>	Pointer to XRpad_VersionInfo structure, which will be filled with the version information.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.1.1.2 HIS_RETURN Acquisition_Set_OnboardOptions (HACQDESC hAcqDesc, BOOL bEnableAck, BOOL bOffset, BOOL bGain, BOOL bPixel)

Activate image acquisition options for onboard offset, gain, mean correction and acknowledgement of frames.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>bEnableAck</i>	Enables or disables the Ack mechanism. Image will be stored on sd card for redundancy/safety when not acknowledged
<i>bOffset</i>	Enables or disables onboard offset corr.
<i>bGain</i>	Enables or disables onboard gain corr
<i>bPixel</i>	Enables or disables onboard pixel corr

Note

Acq_wpe_LoadCorrectionImageToBuffer must be used to upload the correction data from sd-card to memory

Acquisition_wpe_SetUniqueImageTag must be used to tag the image for acknowledgement

Acquisition_AcknowledgeImage must be used to acknowledge the image if activated. otherwise it will be stored on the sdcard

Returns

An ErroCode is returned when unsuccessful otherwise HIS_ALL_OK

3.1.1.3 HIS_RETURN Acquisition_Set_OnboardOptionsPostOffset (HACQDESC *hAcqDesc*, BOOL *bNoOnboardCorr*, BOOL *bSendPreviewFrist*, BOOL *bSendFULLFirst*, BOOL *bEnableAckFirst*, BOOL *bEnableAckSecond*, BOOL *bEnableOffsetFirst*, BOOL *bEnablePostOffsetCorr*, BOOL *bGain*, BOOL *bPixel*)

This function defines the image acquisition options for PhotoTimed mode including preview on/off offset/gain/pixel correction on/off and which image will be send to the client.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descripto
	<i>bNo-Onboard-Corr</i>	if this parameters is set all onboard corrections and preview generation are disabled
	<i>bSend-PreviewFrist</i>	if this parameters is set all a preview for the first image is generated and will be send to the client
	<i>bSendFULL-First</i>	if this parameters is set the first image will be send to the client in full size
	<i>bEnableAck-First</i>	if this parameters is set first full size image will be stored to sd-card if not acknowledged
	<i>bEnableAck-Second</i>	if this parameters is set second full size image will be stored to sd-card if not acknowledged
	<i>bEnable-OffsetFirst</i>	if this parameters is set the first image and the preview image (if selected) will be offset corrected with a predefined offset image
	<i>bEnable-PostOffset-Corr</i>	if this parameters is set the first image will be offset corrected with the second image
	<i>bGain</i>	if this parameters is set the images image will be gain corrected with the predefined gain image
	<i>bPixel</i>	if this parameters is set the images image will be pixel corrected with the predefined pixel mask image

Returns

An `ErroCode` is returned when unsuccessful otherwise `HIS_ALL_OK`

*note the clients image buffer defined using `Acquisition_DefineDestBuffers` must have the size to retrieve all selected images (max 3xFull size for Preview, Bright, Offset/-Brightoc)

3.1.1.4 HIS_RETURN Acquisition_Set_OnboardOptionsPostOffsetEx (HACQDESC *hAcqDesc*, BOOL *bNoOnboardCorr*, BOOL *bSendPreviewFrist*, BOOL *bSendFULLFirst*, BOOL *bEnableAckFirst*, BOOL *bEnableAckSecond*, BOOL *bEnableOffsetFirst*, BOOL *bEnablePostOffsetCorr*, BOOL *bGain*, BOOL *bPixel*, BOOL *bStoreOffsetToSD*)

This function can be used to verify the pgototimed mode. It defines the image acquisition options for PhotoTimed mode including preview on/off offset/gain/pixel correction on/off and which image will be send to the client.

Parameters

<code>in</code>	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
	<i>bNoOnboardCorr</i>	if this parameters is set all onboard corrections and preview generation are disabled
	<i>bSendPreviewFrist</i>	if this parameters is set all a preview for the first image is generated and will be send to the client
	<i>bSendFULLFirst</i>	if this parameters is set the first image will be send to the client in full size
	<i>bEnableAckFirst</i>	if this parameters is set first full size image will be stored to sd-card if not acknowledged
	<i>bEnableAckSecond</i>	if this parameters is set second full size image will be stored to sd-card if not acknowledged
	<i>bEnableOffsetFirst</i>	if this parameters is set the first image and the preview image (if selected) will be offset corrected with a predefined offset image
	<i>bEnablePostOffsetCorr</i>	if this parameters is set the first image will be offset corrected with the second image
	<i>bGain</i>	if this parameters is set the images image will be gain corrected with the predefined gain image
	<i>bPixel</i>	if this parameters is set the images image will be pixel corrected with the predefined pixel mask image
	<i>bStoreOffsetToSD</i>	if this parameters is set second full size image will be stored to sd-card (for debug/test only)

Returns

An `ErroCode` is returned when unsuccessful otherwise `HIS_ALL_OK`

*note the clients image buffer defined using `Acquisition_DefineDestBuffers` must have the size to retrieve all selected images (max 3xFull size for Preview, Bright, Offset/-

Brightoc) *note this function is for test/debug only. If bStoreOffsetToSD is enabled the other Acknowledge/StoreToSD functions will be disabled

3.1.1.5 HIS_RETURN Acquisition_wpe_ActivateNetworkConfig (const char * *ipAddress*, int *configIndex*)

This function activates the selected network configuration of the detector.

Parameters

<i>ipAddress</i>	IP adress of the detectors network interface LAN/WLAN
<i>configIndex</i>	config index to activate

Returns

HIS_ALL_OK function successfull otherwise an error code

3.1.1.6 HIS_RETURN Acquisition_wpe_getAvailableSystems (struct [discoveryReply](#) * *reply*, int * *numDevices*, int *timeout*, int *port*)

This function sends a network broadcast and retrieves all available XRpad detectors in the Network.

Parameters

<i>reply</i>	array of discoveryReply structures
<i>numDevices</i>	pointer to a value containing the number of allocated structures when called
<i>timeout</i>	The timeout shall be greater or erqual to 0 (zero). The physical unit of timeout is ms. The value 0 means to make use of the WPE default timeout of about 2000ms.
<i>port</i>	The port value shall be 0 or 57635. The value 0 means to make use of the default port. The value 57635 is the default port. Actually the PKI detectors do not support other ports than the default port.

The parameter is reserved for future use and exists to kepp the interface unchanged.

Note

Please refer to [Acq.h](#) for the parameter definitions. The user has to allocate the [discoveryReply](#) array before calling the function.

- numDevices has to contain the number of allocated structures and will contain the retrieved number of devices after the call. An array length of 10 is recommended
- Timeout and port must be 0 to use the predefined default values

Returns

HIS_ALL_OK function successful otherwise an error code

3.1.1.7 HIS_RETURN Acquisition_wpe_GetNetworkConfigs (const char * *ipAddress*, struct networkConfiguration * *configs*, int * *arrayLength*, int * *activeConfig*)

This function retrieves the Network settings of a specific Detector defined by its ip address.

Parameters

<i>ipAddress</i>	IP adress of the detectors LAN/WLAN interface
<i>configs</i>	array of networkConfiguration structures
<i>arrayLength</i>	number of allocated strucures
<i>activeConfig</i>	currently activated configuration

Note

Please refer to [Acq.h](#) for the parameter definitions.

- The user has to allocate the [networkConfiguration](#) array before calling the function.
- arrayLength has to contain the number of allocated structures and will contain the retrieved number of devices after the call. An array length of 20 is recommended
- activeConfig will contain the current activated configuration

Returns

HIS_ALL_OK function successful otherwise an error code

3.1.1.8 HIS_RETURN Acquisition_wpe_ReadCameraRegisters (const char * *ipAddress*, unsigned long * *buffer*)

This function retrieves the Status register of a specific Detector defined by its ip address.

Parameters

<i>ipAddress</i>	IP ddress of the detector LAN/WLAN
<i>buffer</i>	Buffer to retrieve the contend of the camera register

Note

Please refer to [acq.h](#) for the parameter definitions.

- The user has to allocate the buffer for the register first
- Buffer size has to be EPC_REGISTER_LENGTH

- Buffer will contain Status Register as struct [EPC_REGISTER](#)

Returns

HIS_ALL_OK status register could be retrieved otherwise an error code

3.2 Special XRpad Functions

Functions

- HIS_RETURN [Acq_wpe_LoadCorrectionImageToBuffer](#) (HACQDESC hAcqDesc, const char *pccCorrectionFilePath, [ProcScriptOperation](#) Operation)
Load a correction file from detector SDCARD to a specific correction buffer.
- HIS_RETURN [Acq_wpe_SetImageTransferInterface](#) (const char *ipAddress, [XRpad_DataInterfaceControlEnum](#) eDataInterface)
This function is used to set the desired detector interface for image transfer.
- HIS_RETURN [Acq_wpe_SystemControl](#) (const char *ipAddress, [XRpad_SystemControlEnum](#) eAction)
This function is used to control the status of the detector. It can be used e.g. for reboot, shutdown.
- HIS_RETURN [Acquisition_AcknowledgeImage](#) (HACQDESC hAcqDesc, const char *tag)
Acquisition_AcknowledgeImage Acknowledges the successful receipt of an image.
- HIS_RETURN [Acquisition_AckSDCardForceFsck](#) (HACQDESC hAcqDesc)
If the filesystem was checked and no errors were found, a message will be send to the client. The reception of the message can be acknowledged/removed with this method. No further message will be send.
- HIS_RETURN [Acquisition_AckSDCardForceFsckError](#) (HACQDESC hAcqDesc)
In the case check of filesystem found errors and fixed them, there will be files like FS-CK0000.REC on the SD card and a message will be send to the client. This method won't touch these files. => They have to be handled "manually" through ftp. The reception of the message can be acknowledged with this method. No further message will be send.
- HIS_RETURN [Acquisition_CloseFile](#) (XislFileHandle fileHandle)
Acquisition_CloseFile Closes a file and releases allocated memory.
- HIS_RETURN [Acquisition_CreateFakeShockCriticalLevel](#) (HACQDESC hAcqDesc)
Simulates a shock event at critical level.
- HIS_RETURN [Acquisition_CreateFakeShockWarningLevel](#) (HACQDESC hAcqDesc)
Simulates a shock event at warning level.
- HIS_RETURN [Acquisition_DisableEventCallback](#) (HACQDESC hAcqDesc)
Disconnects the event callback mechanism from from the detector.
- HIS_RETURN [Acquisition_Enable_EMI_Data_Readout](#) (HACQDESC hAcqDesc, unsigned int uiOnOff)
Enables/Disables the EMI data readout and transfer.
- HIS_RETURN [Acquisition_FactoryResetShock](#) (HACQDESC hAcqDesc)
Resets all shock events.
- HIS_RETURN [Acquisition_FreeFTPFileBuffer](#) (void *pdatatbuffer)
This function to frees a filebuffer allocated by Acquisition_GetFTPFile.
- HIS_RETURN [Acquisition_FTP_CloseSession](#) (XislFtpSession session)

Acquisition_FTP_CloseSession Closes an FTP session.

- HIS_RETURN [Acquisition_FTP_InitSession](#) (HACQDESC hAcqDesc, [XisIFtpSession](#) *session)

Acquisition_FTP_InitSession Initializes a FTP session between client and detector (if supported).

- HIS_RETURN [Acquisition_Get_Current_Voltage](#) (HACQDESC hAcqDesc, [DETECTOR_CURRENT_VOLTAGE](#) *pstructCurrentVoltage)

Retrieves the internal Currents and Voltages from the detector.

- HIS_RETURN [Acquisition_GetAutoPowerOnLocations](#) (HACQDESC hAcqDesc, unsigned int *autopoweronlocations)

Retrieves the current locations from which the detector switches on automatically when the location changed.

- HIS_RETURN [Acquisition_GetBatteryStatus](#) (HACQDESC hAcqDesc, [XRpad_BatteryStatus](#) *batteryStatus)

Retrieves the battery status.

- HIS_RETURN [Acquisition_GetChargeMode](#) (HACQDESC hAcqDesc, unsigned char *charge_mode_req, unsigned char *charge_mode_charger)

Retrieves the battery charge mode.

- HIS_RETURN [Acquisition_GetChargeModePAcqDesc](#) (PAcquisitionDesc pAcqDesc, unsigned char *charge_mode_req, unsigned char *charge_mode_charger)

Retrieves the battery charge mode.

- HIS_RETURN [Acquisition_GetFTPFile](#) (const char *ipAddress, const char *filename, void **databuffer, long *filesize)

This function retrieves data from a file from the detector into a memory buffer.

- HIS_RETURN [Acquisition_GetGridSensorStatus](#) (HACQDESC hAcqDesc, unsigned int *uiStatus)

This function retrieves the status of the Grid Sensors.

- HIS_RETURN [Acquisition_GetLocation](#) (HACQDESC hAcqDesc, unsigned int *location)

Retrieves the location info of the detector.

- HIS_RETURN [Acquisition_GetMissedImageCount](#) ([XisIFtpSession](#) session, UINT *count)

Acquisition_FTP_GetMissedFileCount Retrieves the count of missed images stored on the detector.

- HIS_RETURN [Acquisition_GetNetwork](#) (HACQDESC hAcqDesc, unsigned int *network)

Retrieves the network link speed of the detector.

- HIS_RETURN [Acquisition_GetPowerstate](#) (HACQDESC hAcqDesc, unsigned int *powerstate)

Retrieves the current powerstate info of the detector.

- HIS_RETURN [Acquisition_GetProvidedEnhancedFeatures](#) (HACQDESC hAcqDesc, unsigned int *puiEnhancesFeatures)

Acquisition_GetProvidedEnhancedFeatures retrieves whether the Detector Provides and Enhanced features.

- HIS_RETURN [Acquisition_GetSDCardInfo](#) (HACQDESC hAcqDesc, unsigned int *total, unsigned int *avail)
Retrieves the sdcard info about available storage space on the sd-card from the detector.
- HIS_RETURN [Acquisition_GetSDCardTimeout](#) (HACQDESC hAcqDesc, unsigned short *sdcard_timeout)
Retrieves the SD card timeout value of the detector.
- HIS_RETURN [Acquisition_GetTemperatureThresholds](#) (HACQDESC hAcqDesc, unsigned int *threshold_warning, unsigned int *threshold_critical)
Retrieves the warning level and the critical temperature thresholds of the detector.
- HIS_RETURN [Acquisition_IdentifyDevice](#) (HACQDESC hAcqDesc)
trigger device to identify itself by flashing the display
- HIS_RETURN [Acquisition_IsPreviewImage](#) (HACQDESC hAcqDesc, unsigned int *uilsPreview)
This function retrieves the information whether the latest received frame is a preview image.
- HIS_RETURN [Acquisition_OpenMissedImage](#) (XislFtpSession session, UINT index, XislFileHandle *fileHandle)
Acquisition_FTP_OpenMissedImage Retrieves a handle of a missed image file.
- HIS_RETURN [Acquisition_Resend_All_Messages](#) (HACQDESC hAcqDesc)
resets all messages to be resend
- HIS_RETURN [Acquisition_Reset_OnboardOptions](#) (HACQDESC hAcqDesc)
Reset all onboard features like preview and corrections.
- HIS_RETURN [Acquisition_ResetOnboardShockEvent](#) (HACQDESC hAcqDesc, unsigned int latestShock_Timestamp)
Resets the onboard shockevent. Log files will not be deleted.
- HIS_RETURN [Acquisition_ResetTemperatureTimeout](#) (HACQDESC hAcqDesc)
Reset the timeout the detector waits before a thermal shutdown.
- HIS_RETURN [Acquisition_Set_FPGA_Power_Mode](#) (HACQDESC hAcqDesc, unsigned int uiMode)
switches the Analog control FPGA On (IDLE or off DEEP_SLEEP) and waits for the device to be ready
- HIS_RETURN [Acquisition_Set_OnboardOffsetImageAcquisition](#) (HACQDESC hAcqDesc, BOOL bEnable, BOOL bSend, BOOL bStoreSD)
Activate onboard offset acquisition to acquire an image into the onboard offset correction buffer.
- HIS_RETURN [Acquisition_Set_OnboardOptionPreview](#) (HACQDESC hAcqDesc, BOOL bEnablePreview, BOOL bPreviewOptionSendFull, [OnboardBinning-Mode](#) eMode, unsigned int uiSelectedScript)
This functions Enables/Disables Onboard Preview generation for single shot acquisition.
- HIS_RETURN [Acquisition_SetAutoPowerOnLocations](#) (HACQDESC hAcqDesc, unsigned int autopoweronlocations)
sets the locations from which the detector switches on automatically when location changed

- HIS_RETURN [Acquisition_SetChargeMode](#) (HACQDESC hAcqDesc, unsigned char charge_mode)
Sets the battery charge mode.
- HIS_RETURN [Acquisition_SetChargeModePAcqDesc](#) (PAcquisitionDesc pAcqDesc, unsigned char charge_mode)
Sets the battery charge mode.
- HIS_RETURN [Acquisition_SetEventCallback](#) (HACQDESC hAcqDesc, XIS_EventCallback EventCallback, void *userData)
Sets a callback function to retrieve events.
- HIS_RETURN [Acquisition_SetFakeTemperature](#) (HACQDESC hAcqDesc, BOOL bEnableFakeMode, int iFakeTemperature)
Acquisition_SetFakeTemperature Enables or disables a fake temperature mode of the virtual temperature sensor on XRpad detectors.
- HIS_RETURN [Acquisition_SetFTPFile](#) (const char *ipAddress, const char *filename, void *databuffer, long filesize)
This function stores a databuffer on the sd card folder /mnt/sdcard.
- HIS_RETURN [Acquisition_SetIdleTimeout](#) (HACQDESC hAcqDesc, unsigned short timeout)
This function sets the idle timeout period.
- HIS_RETURN [Acquisition_SetNetworkSpeed](#) (HACQDESC hAcqDesc, unsigned int network)
sets the max speed of the network of the detector
- HIS_RETURN [Acquisition_SetPhototimedParams](#) (HACQDESC hAcqDesc, unsigned short usNrOfScrubs, unsigned short usMaxDelay)
This function configures detector Parameters for Phototimed mode.
- HIS_RETURN [Acquisition_SetPrivateKey](#) (HACQDESC hAcqDesc, unsigned char(*key_old)[64], unsigned char(*key_new)[64])
Updates the private key stored on the device.
- HIS_RETURN [Acquisition_SetSDCardForceFsck](#) (HACQDESC hAcqDesc)
Set flag to force check of filesystem on SD card. After setting the flag the detector has to be rebooted to execute the fsck.
- HIS_RETURN [Acquisition_SetSDCardTimeout](#) (HACQDESC hAcqDesc, unsigned short sdcard_timeout)
Sets the timeout after which an unacknowledged image is stored to the SD card.
- HIS_RETURN [Acquisition_SetSystemTime](#) (HACQDESC hAcqDesc, char *c-DateTime)
sets the system time of the detector.
- HIS_RETURN [Acquisition_SetTailTimeforTriggerMode](#) (HACQDESC hAcqDesc, unsigned short usTailTime, XIS_DetectorTriggerMode eTriggerMode)
set the tail time in front of the readout in msec
- HIS_RETURN [Acquisition_SetTemperatureThresholds](#) (HACQDESC hAcqDesc, unsigned int threshold_warning, unsigned int threshold_critical)
Sets the warning level and the critical threshold for the surface temperature of the detector.
- HIS_RETURN [Acquisition_SetTemperatureTimeout](#) (HACQDESC hAcqDesc, unsigned short timeout)

This function sets the timeout period the detector waits before the thermal shutdown.

- HIS_RETURN [Acquisition_VerifyGenuineness](#) (HACQDESC hAcqDesc, char(*msg)[128], size_t *msg_len, unsigned char(*md)[20])

Acquisition_VerifyGenuineness Verifies if the device is genuine.

- HIS_RETURN [Acquisition_wpe_ChangeNetworkConfig](#) (const char *ipAddress, int configIndex, struct [networkConfiguration](#) *config)

Change the network configuration via API call.

- HIS_RETURN [Acquisition_wpe_FillDefaultNetworkConfiguration](#) (struct [networkConfiguration](#) *config)

this function will fill a structure with default network settings

- HIS_RETURN [Acquisition_wpe_ForceIP](#) (const char *macAddress, struct [networkConfiguration](#) *config, int port, int *isAnswered)

Send a device a "force IP" request as broadcast to be used a temporary IP-adress.

- HIS_RETURN [Acquisition_wpe_GetExamFlag](#) (const char *ipAddress, unsigned long *pExamFlag)

This function retrieved the status of the exam flag which is set during a running acquisition.

- HIS_RETURN [Acquisition_wpe_SetMaxOnboardCorrValue](#) (HACQDESC hAcqDesc, unsigned short usMax, unsigned short usReplace)

Set Max Value for onboard corrections.

- HIS_RETURN [Acquisition_wpe_SetUniqueImageTag](#) (HACQDESC hAcqDesc, const char *imageTag)

Sets an Image Tag via wpe200 library.

3.2.1 Function Documentation

3.2.1.1 HIS_RETURN Acq_wpe_LoadCorrectionImageToBuffer (HACQDESC hAcqDesc, const char * pccCorrectionFilePath, ProcScriptOperation Operation)

Load a correction file from detector SDCARD to a specific correction buffer.

Parameters

<i>hAcqDesc</i>	handle to the detector
<i>pcc-Correction-FilePath</i>	Full path of the correction file on the SDCARD including /mnt/sdcard
<i>Operation</i>	0-OFFSET 1-GAIN 2-MEAN

Note

The correction will not be enabled with this function.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.2.1.2 HIS_RETURN Acq_wpe_SetImageTransferInterface (const char * *ipAddress*, XRpad_DataInterfaceControlEnum *eDataInterface*)

This function is used to set the desired detector interface for image transfer.

Parameters

<i>ipAddress</i>	IP-Address of the device to control
<i>eData-Interface</i>	Interface to use (0 - LAN 1 - WLAN)

Note

This function may return -10004 since the device will not be reachable in case of shutdown/reboot

This function overwrites the current settings for the interface temporarily

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.2.1.3 HIS_RETURN Acq_wpe_SystemControl (const char * *ipAddress*, XRpad_SystemControlEnum *eAction*)

This function is used to control the status of the detector. It can be used e.g. for reboot, shutdown.

Parameters

<i>ipAddress</i>	IP-Address of the device to control
<i>eAction</i>	Action to execute

Note

This function may return -10004 since the device will not be reachable in case of shutdown/reboot

WPE_SYSTEM_CONTROL_REBOOT - requires AcquisitonClose and re init
XRpad_SYSTEM_CONTROL_SET_DEEP_SLEEP / XRpad_SYSTEM_CONTROL_SET_IDLE If switching to deep sleep and switching to idle is used it has to be checked whether the device is back in IDLE stated by using Acquisition_GetHWHeader(..) and the image tag must be set to a value <>"". Use - Acquisition_Set_FPGA_Power_Mode to controll DEEPSLEEP<> IDLE instead to ensure device is back to idle.

The function will return HIS_ERROR_ACQ_ALREADY_RUNNING when called during a running acquisition.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.2.1.4 HIS_RETURN Acquisition_AcknowledgeImage (HACQDESC *hAcqDesc*, const char * *tag*)

Acquisition_AcknowledgeImage Acknowledges the successful receipt of an image.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>tag</i>	The tag of the image to be acknowledged.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.5 HIS_RETURN Acquisition_AckSDCardForceFsck (HACQDESC *hAcqDesc*)

If the filesystem was checked and no errors were found, a message will be send to the client. The reception of the message can be acknowledged/removed with this method. No further message will be send.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.6 HIS_RETURN Acquisition_AckSDCardForceFsckError (HACQDESC *hAcqDesc*)

In the case check of filesystem found errors and fixed them, there will be files like FSC-K0000.REC on the SD card and a message will be send to the client. This method won't touch these files. => They have to be handled "manually" through ftp. The reception of the message can be acknowledged with this method. No further message will be send.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.7 HIS_RETURN Acquisition_CloseFile (XislFileHandle fileHandle)

Acquisition_CloseFile Closes a file and releases allocated memory.

Parameters

in	<i>fileHandle</i>	A valid file handle.
----	-------------------	----------------------

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Warning

The fileHandle is no longer valid after this operation. Further usage of this handle may lead to undefined behavior.

3.2.1.8 HIS_RETURN Acquisition_CreateFakeShockCriticalLevel (HACQDESC hAcqDesc)

Simulates a shock event at critical level.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.9 HIS_RETURN Acquisition_CreateFakeShockWarningLevel (HACQDESC hAcqDesc)

Simulates a shock event at warning level.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.10 HIS_RETURN Acquisition_DisableEventCallback (HACQDESC *hAcqDesc*)

Disconnects the event callback mechanism from from the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Note

Beware that the detector will shut down after the configured timeout after this disconnect. To re-enable the callback-mechanism, call Acquisition_SetEventCallback again.

Returns

Returns always HIS_ALL_OK.

3.2.1.11 HIS_RETURN Acquisition_Enable_EMI_Data_Readout (HACQDESC *hAcqDesc*, unsigned int *uiOnOff*)

Enables/Disables the EMI data readout and transfer.

Note

Use Acquisition_GbIF_DiscoverDevices in a preceding call to discover the network.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>uiOnOff</i>	Defines whether the EMI readout shall be enabled 1 or disabled 0

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.2.1.12 HIS_RETURN Acquisition_FactoryResetShock (HACQDESC *hAcqDesc*)

Resets all shock events.

Parameters

<i>in</i>	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
-----------	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

Log file will not be deleted

3.2.1.13 HIS_RETURN Acquisition_FreeFTPFileBuffer (void * *pdatabuffer*)

This function to frees a filebuffer allocated by Acquisition_GetFTPFile.

Parameters

<i>pdatabuffer</i>	Pointer to buffer that should be freed up.
--------------------	--

Note

customer/client app has to set the pointer to null

Returns

Always HIS_ALL_OK

3.2.1.14 HIS_RETURN Acquisition_FTP_CloseSession (XislFtpSession *session*)

Acquisition_FTP_CloseSession Closes an FTP session.

Parameters

<i>in</i>	<i>session</i>	A valid FTP session descriptor.
-----------	----------------	---------------------------------

Returns

Always return HIS_ALL_OK-

3.2.1.15 HIS_RETURN Acquisition_FTP_InitSession (HACQDESC *hAcqDesc*, XislFtpSession * *session*)

Acquisition_FTP_InitSession Initializes a FTP session between client and detector (if supported).

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>session</i>	If the function succeeds, this parameter retrieves a descriptor of the current FTP session.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.16 HIS_RETURN Acquisition_Get_Current_Voltage (HACQDESC *hAcqDesc*, DETECTOR_CURRENT_VOLTAGE * *pstructCurrentVoltage*)

Retrieves the internal Currents and Voltages from the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>pstruct-Current-Voltage</i>	Will retrieve internal detector Voltages and Currents

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.17 HIS_RETURN Acquisition_GetAutoPowerOnLocations (HACQDESC *hAcqDesc*, unsigned int * *autopoweronlocations*)

Retrieves the current locations from which the detector switches on automatically when the location changed.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>autopoweron-locations</i>	Will retrieve the bitwise coded location for auto power on from the detector

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.18 HIS_RETURN Acquisition_GetBatteryStatus (HACQDESC *hAcqDesc*, XRpad_BatteryStatus * *batteryStatus*)

Retrieves the battery status.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>battery-Status</i>	Pointer to an XRpad_BatteryStatus structure to retrieve the battery status.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.19 HIS_RETURN Acquisition_GetChargeMode (HACQDESC *hAcqDesc*, unsigned char * *charge_mode_req*, unsigned char * *charge_mode_charger*)

Retrieves the battery charge mode.

The *charge_mode_req* requested by the software defines the max charge mode The *charge_mode_charger* defines the max charge allowed by the charger.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>charge_mode_req</i>	The value of charge mode requested by the software, 0..3 0 no charging, 3 max charging
out	<i>charge_mode_charger</i>	The value of charge mode set by the charger, 0..3 0 no charging, 3 max charging

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.20 HIS_RETURN Acquisition_GetChargeModePAcqDesc (PAcquisitionDesc *pAcqDesc*, unsigned char * *charge_mode_req*, unsigned char * *charge_mode_charger*)

Retrieves the battery charge mode.

The *charge_mode_req* requested by the software defines the max charge mode The *charge_mode_charger* defines the max charge allowed by the charger.

Parameters

in	<i>pAcqDesc</i>	pointer to a valid Acquisition Descriptor.
out	<i>charge_mode_req</i>	The value of charge mode requested by the software, 0..3 0 no charging, 3 max charging
out	<i>charge_mode_charger</i>	The value of charge mode set by the charger, 0..3 0 no charging, 3 max charging

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.21 HIS_RETURN Acquisition_GetFTPFile (const char * *ipAddress*, const char * *filename*, void ** *databuffer*, long * *filesize*)

This function retrieves data from a file from the detector into a memory buffer.

Parameters

<i>ipAddress</i>	Pointer to char array of ip address
<i>filename</i>	Pointer to value to retrieve the actual Field Of View mode.
<i>databuffer</i>	Pointer to databuffer
<i>filesize</i>	Pointer to long value of retrieved filesize

Note

default directory is /mnt/sdcard/

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.2.1.22 HIS_RETURN Acquisition_GetGridSensorStatus (HACQDESC *hAcqDesc*, unsigned int * *uiStatus*)

This function retrieves the status of the Grid Sensors.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>uiStatus</i>	pointer to a unsigned int value to retrieve the current status

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

0 - No Grid 1 - Grid Type1 2 - Grid Type2 3 - Grid Type3

3.2.1.23 HIS_RETURN Acquisition_GetLocation (HACQDESC *hAcqDesc*, unsigned int * *location*)

Retrieves the location info of the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>location</i>	Will retrieve location of detector

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

**3.2.1.24 HIS_RETURN Acquisition_GetMissedImageCount (XisIFtpSession *session*,
UINT * *count*)**

Acquisition_FTP_GetMissedFileCount Retrieves the count of missed images stored on the detector.

Parameters

in	<i>session</i>	A valid FTP session descriptor.
out	<i>count</i>	If the function succeeds, this parameter retrieves the count of missed images stored on the detector.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

**3.2.1.25 HIS_RETURN Acquisition_GetNetwork (HACQDESC *hAcqDesc*, unsigned int *
network)**

Retrieves the network link speed of the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>network</i>	will retrieve network link speed of detector

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

**3.2.1.26 HIS_RETURN Acquisition_GetPowerstate (HACQDESC *hAcqDesc*, unsigned
int * *powerstate*)**

Retrieves the current powerstate info of the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>powerstate</i>	Will retrieve powerstate of the detector

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

powerstate 4 - IDLE; 3 - DEEP SLEEP

3.2.1.27 HIS_RETURN Acquisition_GetProvidedEnhancedFeatures (HACQDESC *hAcqDesc*, unsigned int * *puiEnhancesFeatures*)

Acquisition_GetProvidedEnhancedFeatures retrieves whether the Detector Provides and Enhanced features.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>pui-Enhances-Features</i>	pointer to a vaule retrieving the whether enhanced modes are provided

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.28 HIS_RETURN Acquisition_GetSDCardInfo (HACQDESC *hAcqDesc*, unsigned int * *total*, unsigned int * *avail*)

Retrieves the sdcard info about available storage space on the sd-card from the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>total</i>	Will retrieve total storage in MB
out	<i>avail</i>	Will retrieve available storage in MB

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.29 HIS.RETURN Acquisition_GetSDCardTimeout (HACQDESC *hAcqDesc*, unsigned short * *sdcard_timeout*)

Retrieves the SD card timeout value of the detector.

See Acquisition_SetSDCardTimeout for more information.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>sdcard_timeout</i>	Pointer to the variable to retrieve SD card timeout value in seconds.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.30 HIS.RETURN Acquisition_GetTemperatureThresholds (HACQDESC *hAcqDesc*, unsigned int * *threshold_warning*, unsigned int * *threshold_critical*)

Retrieves the warning level and the critical temperature thresholds of the detector.

See Acquisition_SetTemperatureThresholds for more information.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>threshold_warning</i>	Pointer to the variable to retrieve the warning level temperature in millidegree Celsius.
out	<i>threshold_critical</i>	Pointer to the variable to retrieve the critical temperature threshold in millidegree Celsius.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.31 HIS.RETURN Acquisition_IdentifyDevice (HACQDESC *hAcqDesc*)

trigger device to identify itself by flashing the display

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.32 HIS_RETURN Acquisition_IsPreviewImage (HACQDESC *hAcqDesc*, unsigned int * *uilsPreview*)

This function retrieves the information whether the latest received frame is a preview image.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>uilsPreview</i>	Flag whether the image was a preview - 1 otherwise - 0

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

This function must be called in the EndframeCallback since the information whether the image is a preview will be overwritten when the next image was retrieved

3.2.1.33 HIS_RETURN Acquisition_OpenMissedImage (XisIFtpSession *session*, UINT *index*, XisIFileHandle * *fileHandle*)

Acquisition_FTP_OpenMissedImage Retrieves a handle of a missed image file.

Parameters

in	<i>session</i>	A valid FTP session descriptor.
in	<i>index</i>	Index between 0 and (count - 1), where count is retrieved by Acquisition_FTP_GetMissedImageCount.
out	<i>fileHandle</i>	If the function succeeds, this parameter retrieves the handle of the missed image file.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Warning

This function allocates a memory buffer to carry the whole data of the file. Call Acquisition_CloseFile to release this memory. Otherwise it will leak.

3.2.1.34 HIS_RETURN Acquisition_Resend_All_Messages (HACQDESC *hAcqDesc*)

resets all messages to be resend

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.35 HIS_RETURN Acquisition_Reset_OnboardOptions (HACQDESC *hAcqDesc*)

Reset all onboard features like preview and corrections.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.36 HIS_RETURN Acquisition_ResetOnboardShockEvent (HACQDESC *hAcqDesc*, unsigned int *latestShock_Timestamp*)

Resets the onboard shockevent. Log files will not be deleted.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
	<i>latestShock-Timestamp</i>	has to be the timestamp of the latest shock that is going to be resetted. to avoid a race condition

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.37 HIS_RETURN Acquisition_ResetTemperatureTimeout (HACQDESC *hAcqDesc*)

Reset the timeout the detector waits before a thermal shutdown.

If the temperature of the detector is above the critical threshold, it waits for a predefined time before it shuts down. This shutdown timeout is reset with this function. You can use Acquisition_SetTemperatureTimeout to set the timeout value.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.38 HIS_RETURN Acquisition_Set_FPGA_Power_Mode (HACQDESC *hAcqDesc*, unsigned int *uiMode*)

switches the Analog control FPGA On (IDLE or off DEEP_SLEEP) and waits for the device to be ready

Parameters

<i>in</i>	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
<i>in</i>	<i>uiMode</i>	0 - Off 1 - On

Returns

Returns HIS_ALL_OK on success or an appropriate error code including HIS_ERROR_NO_FPGA_ACK otherwise.

Note

Currently only valid for network connected detectors

3.2.1.39 HIS_RETURN Acquisition_Set_OnboardOffsetImageAcquisition (HACQDESC *hAcqDesc*, BOOL *bEnable*, BOOL *bSend*, BOOL *bStoreSD*)

Activate onboard offset acquisition to acquire an image into the onboard offset correction buffer.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>bEnable</i>	Enable Store next image to onboard offset correction buffer
<i>bSend</i>	Send out the offset image to client
<i>bStoreSD</i>	Store the image to SDcard when not acknowledged

Returns

An ErrorCode is returned when unsuccessful

3.2.1.40 HIS_RETURN Acquisition_Set_OnboardOptionPreview (HACQDESC *hAcqDesc*, BOOL *bEnablePreview*, BOOL *bPreviewOptionSendFull*, OnboardBinningMode *eMode*, unsigned int *uiSelectedScript*)

This functions Enables/Disables Onboard Preview generation for single shot acquisition.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>bEnable-Preview</i>	Enables or disables the preview option. Binned Image will be send in front of the full size image for preview.
<i>bPreview-OptionSend-Full</i>	Enables or disables the send Fullsize Image option.
<i>eMode</i>	Selected Binning mode for Preview
<i>uiSelected-Script</i>	Must be 0 in current implementation

Returns

An ErroCode is returned when unsuccessful

Note

If Preview is enabled the destination buffer must have the size of 2 full size images

3.2.1.41 HIS_RETURN Acquisition_SetAutoPowerOnLocations (HACQDESC *hAcqDesc*, unsigned int *autopoweronlocations*)

sets the locations from which the detecor switches on automatically when location changed

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>autopoweron-locations</i>	Bitmask to define the locations

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

locations are ored location0 - 1, location1 - 2, location3 - 4, .. , location7 - 128; all on 255

3.2.1.42 HIS_RETURN Acquisition_SetChargeMode (HACQDESC *hAcqDesc*, unsigned char *charge_mode*)

Sets the battery charge mode.

To avoid charging influences on the image quality and to reduce heating up of the device the charge mode can be modified.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>charge_mode</i>	The value of charge mode, 0..3 0 no charging, 3 max charging

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.43 HIS_RETURN Acquisition_SetChargeModePAcqDesc (PAcquisitionDesc *pAcqDesc*, unsigned char *charge_mode*)

Sets the battery charge mode.

To avoid charging influences on the image quality and to reduce heating up of the device the charge mode can be modified.

Parameters

in	<i>pAcqDesc</i>	pointer to a valid Acquisition Descriptor.
in	<i>charge_mode</i>	The value of charge mode, 0..3 0 no charging, 3 max charging

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.44 HIS_RETURN Acquisition_SetEventCallback (HACQDESC *hAcqDesc*, XIS_EventCallback *EventCallback*, void * *userData*)

Sets a callback function to retrieve events.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>EventCallback</i>	Pointer to the callback function which retrieves the events.
	<i>userData</i>	pointer to the user data

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.45 HIS_RETURN Acquisition_SetFakeTemperature (HACQDESC *hAcqDesc*, BOOL *bEnableFakeMode*, int *iFakeTemperature*)

Acquisition_SetFakeTemperature Enables or disables a fake temperature mode of the virtual temperature sensor on XRpad detectors.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>bEnable-FakeMode</i>	TRUE enables the fake mode, FALSE disables it.
in	<i>iFake-Temperature</i>	Temperature to set in 1/1000 degree Celsius. (42500 == 42.5 C). This parameter is ignored if bEnableFakeMode is FALSE.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.46 HIS_RETURN Acquisition_SetFTPFile (const char * *ipAddress*, const char * *filename*, void * *databuffer*, long *filesize*)

This function stores a databuffer on the sd card folder /mnt/sdcard.

Parameters

<i>ipAddress</i>	Pointer to char array of ip address
<i>filename</i>	Pointer to value to retrieve the actual Field Of View mode.
<i>databuffer</i>	Pointer to databuffer
<i>filesize</i>	Pointer to long value of uploaded filesize

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.2.1.47 HIS_RETURN Acquisition_SetIdleTimeout (HACQDESC *hAcqDesc*, unsigned short *timeout*)

This function sets the idle timeout period.

To reduce power consumption and extend battery lifetime, the XRpad detectors shut down automatically after 10 minutes in idle state.

Use this function to adjust the idle-timeout period, after which the automatic shutdown is initiated.

Note

The shutdown timeout is NOT permanently stored and must be (re-)set after every (re-)boot of the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>timeout</i>	The value of the idle timeout period in seconds. The timeout value must be at least 20 seconds.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.48 HIS_RETURN Acquisition_SetNetworkSpeed (HACQDESC *hAcqDesc*, unsigned int *network*)

sets the max speed of the network of the detector

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>network</i>	Will provide network of detector

Note

this may take several seconds to settle

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.49 HIS_RETURN Acquisition_SetPhototimedParams (HACQDESC *hAcqDesc*, unsigned short *usNrOfScrubs*, unsigned short *usMaxDelay*)

This function configures detector Parameters for Phototimed mode.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>usNrOfScrubs</i>	after bright image readout before offset readout (default is 4 min 1)
in	<i>usMaxDelay</i>	Max Exposure Delay in milliseconds (default is 5000)

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

currently only available fro XRPad2

3.2.1.50 HIS_RETURN Acquisition_SetPrivateKey (HACQDESC hAcqDesc, unsigned char(*) key_old[64], unsigned char(*) key_new[64])

Updates the private key stored on the device.

The private key must be kept secret to ensure genuineness. It must have exactly 512 bits (^= 64 bytes). Note that it is transmitted unencrypted. Please note that you need the current private key as authorization. This implies that you will not be able to update the key, once you loose your current one.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>key_old</i>	Pointer to an array containing the old key in raw binary data.
in	<i>key_new</i>	Pointer to an array containing the new key in raw binary data.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

currently only available with XRPad2

3.2.1.51 HIS_RETURN Acquisition_SetSDCardForceFsck (HACQDESC hAcqDesc)

Set flag to force check of filesystem on SD card. After setting the flag the detector has to be rebooted to execute the fsck.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.52 HIS_RETURN Acquisition_SetSDCardTimeout (HACQDESC *hAcqDesc*, unsigned short *sdcard_timeout*)

Sets the timeout after which an unacknowledged image is stored to the SD card.

The timeout applies only to images that have to be acknowledged by the application. When the image is stored, the client is notified by an event as far as an event- callback function is specified by using Acquisition_SetEventCallback. An image may be stored before this timeout is reached under the following circumstances:

- There are more than two internal buffers already in use. The oldest image is then stored to SD, automatically
- The handle was closed using Acquisition_Close or Acquisition_CloseAll
- The connection to the detector was lost/timed out

This might happen

1. if the detector shuts down for any reason
2. if the client application shuts down for any reason

Note

Please note that, due to SD card performance, there is some delay between the timeout and the moment the file is fully written.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>sdcard_timeout</i>	SD card timeout value in seconds.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.53 HIS_RETURN Acquisition_SetSystemTime (HACQDESC *hAcqDesc*, char * *cDateTime*)

sets the system time of the detector.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>cDateTime</i>	Pointer to a datetime string

format: YYYY.MM.DD-hh:mm:ss (24h format), the timezone of the detector is always UTC

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.54 HIS_RETURN Acquisition_SetTailTimeforTriggerMode (HACQDESC hAcqDesc, unsigned short usTailTime, XIS_DetectorTriggerMode eTriggerMode)

set the tail time in front of the readout in msec

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>usTailTime</i>	Tail time to set in msec
in	<i>eTrigger-Mode</i>	TriggerMode for which the tail time shall be set

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

Default for Framewise is 0msec default for DDD, AED, Phototimed is 20msec. Only available for XRpad detectors

3.2.1.55 HIS_RETURN Acquisition_SetTemperatureThresholds (HACQDESC hAcqDesc, unsigned int threshold_warning, unsigned int threshold_critical)

Sets the warning level and the critical threshold for the surface temperature of the detector.

The surface temperature of the detector is determined by several temperature sensors inside the device. Use this function to set the threshold to set the warning level and the critical threshold.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>threshold_warning</i>	Warning temperature in milldegree celsius. If the surface temperature is above this temperature, the API will trigger a warning event.
in	<i>threshold_critical</i>	Critical temperature in milldegree celsius. If the surface temperature is above this temperature, the detector will shut down after a predefined timeout.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.56 HIS_RETURN Acquisition_SetTemperatureTimeout (HACQDESC *hAcqDesc*, unsigned short *timeout*)

This function sets the timeout period the detector waits before the thermal shutdown.

If the temperature of the detector is above the critical threshold, it waits for a predefined time before it shuts down. This shutdown timeout might be reset by using Acquisition_ResetTemperatureTimeout or is always resetted if an image is acquired.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>timeout</i>	Timeout value in seconds. A value less than 20 seconds is not permitted.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.57 HIS_RETURN Acquisition_VerifyGenuineness (HACQDESC *hAcqDesc*, char(*) *msg*[128], size_t * *msg_len*, unsigned char(*) *md*[20])

Acquisition_VerifyGenuineness Verifies if the device is genuine.

The genuineness is verified by utilizing the HMAC-SHA1 algorithm. The host name of the device serves as message and is delivered as output parameter *msg*. The key is by default an array of 64 bytes of zeros. The key might be changed by using Acquisition_SetPrivateKey.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
out	<i>msg</i>	Pointer to an int[128] array. Retrieves the message. This is equal to the host name of the device. This string is NULL-terminated. The terminating NULL-byte itself is not included in the calculation of the hash.
out	<i>msg_len</i>	The message length.
out	<i>md</i>	The message digest (hash) of <i>msg</i> as described above in raw binary data.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.58 HIS_RETURN Acquisition_wpe_ChangeNetworkConfig (const char * *ipAddress*, int *configIndex*, struct networkConfiguration * *config*)

Change the network configuration via API call.

Parameters

<i>ipAddress</i>	ID address of the detector
<i>configIndex</i>	Configuration index to change
<i>config</i>	networkConfiguration structure to store

Note

config ([networkConfiguration](#)) has some readonly members, please check reference for that!!

Returns

returns HIS_ALL_OK when successfull, An ErroCode is returned when unsuccessful

3.2.1.59 HIS_RETURN Acquisition_wpe_FillDefaultNetworkConfiguration (struct networkConfiguration * *config*)

this function will fill a structure with default network settings

Parameters

<i>config</i>	pointer to the "struct networkConfiguration" to use.
---------------	--

Returns

HIS_ALL_OK if ok or a error code.
The error code is a HIS_ERROR_... or the abs of an WPE_ERR_... code.

3.2.1.60 HIS_RETURN Acquisition_wpe_ForceIP (const char * *macAddress*, struct networkConfiguration * *config*, int *port*, int * *isAnswered*)

Send a device a "force IP" request as broadcast to be used a temporary IP-adress.

Parameters

<i>macAddress</i>	the devices MAC address using the format "00:0A:35:11:DE:D1"
<i>config</i>	pointer to the "struct networkConfiguration" to use.
<i>port</i>	The port number to use, if 0, the default port is used.
<i>isAnswered</i>	return 1 if the request was answered or 0 if not

Note

config ([networkConfiguration](#)) has some readonly members, please check reference for that!!

Returns

HIS_ALL_OK if ok or a negative error code

3.2.1.61 HIS_RETURN Acquisition_wpe_GetExamFlag (const char * *ipAddress*, unsigned long * *pExamFlag*)

This function retrieved the status of the exam flag which is set during a running acquisition.

Parameters

<i>ipAddress</i>	IP-Address of the device to control
<i>pExamFlag</i>	pointer to an unsigned long value to retrieve the exam flag status (1 - exam ongoing; any other - no exam ongoing)

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.2.1.62 HIS_RETURN Acquisition_wpe_SetMaxOnboardCorrValue (HACQDESC *hAcqDesc*, unsigned short *usMax*, unsigned short *usReplace*)

Set Max Value for onboard corrections.

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
in	<i>usMax</i>	Max Values for corrections of current value exceeds usMax its replaced by usReplace
in	<i>usReplace</i>	used to replace values > usMax

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.2.1.63 HIS_RETURN Acquisition_wpe_SetUniqueImageTag (HACQDESC *hAcqDesc*, const char * *imageTag*)

Sets an Image Tag via wpe200 library.

If available for this detector, an image tag is set for the next image. The tag is used for the Ack-Mechanism, available via `Acquisition_Set_OnboardOptions`.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>imageTag</i>	Image tag to be used as filename for autosave when not acknowledged and to acknowledge the image

Returns

An `ErroCode` is returned when unsuccessful

Warning

Always use a unique identifier as image tag. If the image tag is ambiguous, the image may not be explicit identified and may remain on the SD Card storage of the detector. The SD Card might run out of memory in that case.

3.3 Common Init functions

Functions

- HIS_RETURN [Acquisition_Close](#) (HACQDESC hAcqDesc)
Hardware and XISL are closed by this routine. The acquisition descriptor is no longer valid.
- HIS_RETURN [Acquisition_CloseAll](#) ()
This function closes / shuts down all connections to detectors via frame grabbers and IP connections currently allocated by the XISL. All acquisition descriptor structures returned by other functions are invalid after this function call.
- HIS_RETURN [Acquisition_EnableLogging](#) (BOOL onOff)
The purpose of this function is to toggle (enable/disable) the internal logging of the Xisl.
- HIS_RETURN [Acquisition_EnumSensors](#) (UINT *pdwNumSensors, BOOL b-EnableIRQ, BOOL bAlwaysOpen)
This function enumerates all currently connected sensors. All recognized sensors are initialized automatically. To get the HACQDESC of every sensor, use Acquisition_GetNextSensor. For a programming example see the initialization part of the XISL demonstration. Note: In case of Network/IP Sensors only sensors with standart-gateway equal to zero are initialized automatically.
- HIS_RETURN [Acquisition_GetCommChannel](#) (HACQDESC hAcqDesc, UINT *pdwChannelType, int *pnChannelNr)
This function returns the type of the communication device that is used to transfer data from the detector into the PC RAM.
- HIS_RETURN [Acquisition_GetConfiguration](#) (HACQDESC hAcqDesc, UINT *dwFrames, UINT *dwRows, UINT *dwColumns, UINT *dwDataType, UINT *dwSortFlags, BOOL *bIRQEnabled, DWORD *dwAcqType, DWORD *dw-SystemId, DWORD *dwSyncMode, DWORD *dwHwAccess)
This function retrieves all important acquisition parameters, that can be set by - Acquisition_Init or that are set by the self configuration mechanisms of the XISL.
- HIS_RETURN [Acquisition_GetHwHeaderInfo](#) (HACQDESC hAcqDesc, CHw-HeaderInfo *pInfo)
This function returns the contents of the camera's hardware header in a CHwHeader-Info structure.
- HIS_RETURN [Acquisition_GetHwHeaderInfoEx](#) (HACQDESC hAcqDesc, CHw-HeaderInfo *pInfo, CHwHeaderInfoEx *pInfoEx)
This function acquires the frame header of the connected detector. If dwHeaderID in the CHwHeaderInfo structure is 14 (1621detectors) pInfoEx will retrieve the extended header, otherwise the structure will be filled with 0xFFFF.
- HIS_RETURN [Acquisition_GetIntTimes](#) (HACQDESC hAcqDesc, double *dblInt-Time, int *nIntTimes)
This function retrieves the current integration times.
- HIS_RETURN [Acquisition_GetNextSensor](#) (ACQDESCPOS *Pos, HACQDESC *phAcqDesc)
You can use this function to iterate through all recognized sensors in the system. for a programming example see the initialization part of the XISL demonstration.

- HIS_RETURN [Acquisition_Init](#) ([HACQDESC](#) *phAcqDesc, DWORD dwBoardType, int nChannelNr, BOOL bEnableIRQ, UINT Rows, UINT Columns, UINT dwSortFlag, BOOL bSelfInit, BOOL bAlwaysOpen)

The Acquisition_Init function initializes the frame grabber board and the corresponding driver. It enables desired hardware interrupts, prepares acquisition threads, defines callback functions to react on acquisition status changes and tests for sufficient memory space for DMA (direct memory access).

- HIS_RETURN [Acquisition_SetCallbacksAndMessages](#) ([HACQDESC](#) hAcqDesc, HWND hWnd, UINT dwErrorMsg, UINT dwLoosingFramesMsg, void(CALLBACK *lpfnEndFrameCallback)([HACQDESC](#)), void(CALLBACK *lpfnEndAcqCallback)([HACQDESC](#)))

The Acquisition_SetCallbacksAndMessages function defines callback functions to react on acquisition status changes. For a programming example see the initialization part of the XISL demonstration.

- HIS_RETURN [Acquisition_SetLogLevel](#) ([XisLoggingLevels](#) xislLogLvl)

This function is to set the logging level of the Xisl logging.

- *****HIS_RETURN [Acquisition_SetLogOutput](#) (const char *filePath, BOOL consoleOnOff)

*****/

- HIS_RETURN [Acquisition_TogglePerformanceLogging](#) (BOOL onOff)

This function is used to toggle the performance logging. The performance logged will measure the time it takes for a particular function call to execute.

3.3.1 Function Documentation

3.3.1.1 HIS_RETURN Acquisition_Close (HACQDESC hAcqDesc)

Hardware and XISL are closed by this routine. The acquisition descriptor is no longer valid.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
-----------------	--

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.2 HIS_RETURN Acquisition_CloseAll ()

This function closes / shuts down all connections to detectors via frame grabbers and IP connections currently allocated by the XISL. All acquisition descriptor structures returned by other functions are invalid after this function call.

Returns

Value If the function is successful it returns zero, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.3 HIS_RETURN Acquisition_EnableLogging (BOOL *onOff*)

The purpose of this function is to toggle (enable/disable) the internal logging of the Xisl.

Parameters

<i>onOff</i>	Boolean value used to determine whether to turn internal logging on or off
--------------	--

Returns

If the function is successful it returns zero, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.4 HIS_RETURN Acquisition_EnumSensors (UINT * *pdwNumSensors*, BOOL *bEnableIRQ*, BOOL *bAlwaysOpen*)

This function enumerates all currently connected sensors. All recognized sensors are initialized automatically. To get the HACQDESC of every sensor, use Acquisition_GetNextSensor. For a programming example see the initialization part of the XISL demonstration. Note: In case of Network/IP Sensors only sensors with standart-gateway equal to zero are initialized automatically.

Parameters

<i>pdwNumSensors</i>	Address of a 4 byte integer that receives the number of recognized sensors.
<i>bEnableIRQ</i>	If you want to run the acquisition in polling mode set this parameter to zero. If you want to enable hardware interrupts set the parameter to one.
<i>bAlwaysOpen</i>	If this parameter is TRUE the XISL is capturing all communication port regardless if this port is already opened by other processes running on the system. The use of this option is only recommended in debug versions of your applications because it is not possible to free all system resources allocated by another process.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.5 HIS_RETURN Acquisition_GetCommChannel (HACQDESC *hAcqDesc*, UINT * *pdwChannelType*, int * *pnChannelNr*)

This function returns the type of the communication device that is used to transfer data from the detector into the PC RAM.

Parameters

<i>hAcqDesc</i>	Pointer to HACQDESC.
<i>pdwChannelType</i>	Address of a 4 byte integer that receives an id of the currently open communication device.
<i>pnChannelNr</i>	Address of a 4 byte integer that receives the number of communication channel. If the above mentioned communication device is a frame grabber this number is unique to identify the grabber if more than one grabber of one type is installed on the system. (see frame grabber installation description). If the communication device is an RS232 interface then this number contains the COM port number

Note

- 0 HIS_BOARD_TYPE_NONE no device (not valid)
- 1 HIS_BOARD_TYPE_ELTEC_XRD-FG or XRD-FGe Frame Grabber
- 8 HIS_BOARD_TYPE_ELTEC_XRD_FGX XRD-FGX frame grabber
- 16 HIS_BOARD_TYPE_ELTEC_XRD_FGE_Opto XRD-FGe Opto
- 32 HIS_BOARD_TYPE_ELTEC_GbIF GigabitEthernet
- 96 HIS_BOARD_TYPE_ELTEC_EMBEDDED Embedded Detector (e.g. X-Rpad)

Returns

If the initialization is successful HIS_ALL_OK is returned, otherwise the return value is greater. Call Acquisition_GetErrorCode to get extended information.

3.3.1.6 HIS_RETURN Acquisition_GetConfiguration (HACQDESC *hAcqDesc*, UINT * *dwFrames*, UINT * *dwRows*, UINT * *dwColumns*, UINT * *dwDataType*, UINT * *dwSortFlags*, BOOL * *bIRQEnabled*, DWORD * *dwAcqType*, DWORD * *dwSystemId*, DWORD * *dwSyncMode*, DWORD * *dwHwAccess*)

This function retrieves all important acquisition parameters, that can be set by - Acquisition_Init or that are set by the self configuration mechanisms of the XISL.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwFrames</i>	Number of frames of acquisition buffer.
<i>dwRows</i>	Number of rows of the sensor.

<i>dwColumns</i>	Number of columns of the sensor.
<i>dwDataType</i>	Type of data of acquisition buffer (should always be two = unsigned short).
<i>dwSortFlags</i>	Type of sorting. This value depends on camera and used sensor (see sorting schemes).
<i>bIRQ-Enabled</i>	Retrieves a flag that indicates if interrupts are enabled (see Acquisition_Init, hardware interrupts) or if the hardware is running in polling mode. In interrupt mode this parameter is equal to one and zero in the other case.
<i>dwAcqType</i>	Only for internal use.
<i>dwSystemId</i>	PROM identification number of the used camera. This number is only important if the camera operates objectionably and you need any support from PerkinElmer.
<i>dwSync-Mode</i>	This parameter receives the sync mode, see table.
<i>dwHw-Access</i>	This parameter receives the hardware access parameter, that is programmed into the camera. If you have to use this parameter (that depends from your contract with PerkinElmer) please contact PerkinElmer for the possible values.

Note

- HIS_SYNCMODE_FREE_RUNNING The sensor is operating in free running mode.
- HIS_SYNCMODE_EXTERNAL_TRIGGER The sensor is operating in triggered mode. Frames are only send if an external trigger signal is applied to the camera.
- HIS_SYNCMODE_INTERNAL_TIMER The synchronization signal can be generated by the internal timer of the frame grabber (see Acquisition_SetTimerSync)
- HIS_SYNCMODE_SOFT_TRIGGER The synchronization signal can be generated by software (see Acquisition_SetFrameSync).

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.7 HIS_RETURN Acquisition_GetHwHeaderInfo (HACQDESC hAcqDesc, CHwHeaderInfo * pInfo)

This function returns the contents of the camera's hardware header in a [CHwHeaderInfo](#) structure.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pInfo</i>	Pointer to a Structure of type CHwHeaderInfo that contains the contents of the camera's hardware header necessary for self configuration features.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.8 HIS_RETURN Acquisition_GetHwHeaderInfoEx (HACQDESC hAcqDesc, CHwHeaderInfo * pInfo, CHwHeaderInfoEx * pInfoEx)

This function acquires the frame header of the connected detector. If dwHeaderID in the [CHwHeaderInfo](#) structure is 14 (1621detectors) pInfoEx will retrieve the extended header, otherwise the structure will be filled with 0xFFFF.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pInfo</i>	Pointer to Structure of type CHwHeaderInfo to retrieve the detector's hardware header.
<i>pInfoEx</i>	Pointer to Structure of type CHwHeaderInfoEx to retrieve the detector's hardware header when available. Can be NULL.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

* //val 2009-08-13 wait for frame end before retrieving the header

3.3.1.9 HIS_RETURN Acquisition_GetIntTimes (HACQDESC hAcqDesc, double * dblIntTime, int * nIntTimes)

This function retrieves the current integration times.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dblIntTime</i>	Pointer to an array of 8 byte floating point numbers. This array must contain at least 8 entries.
<i>nIntTimes</i>	This parameter contains the number of maximum entries in the array of 8 byte floating point numbers pointed to by *dblIntTime. After return of the function this variable provides the number of available integration times.

```

double dblIntTimes[8];
int nIntTimes = 8;
if (Acquisition_GetIntTimes(hAcqDesc, dblIntTimes, &nIntTimes) != HIS_ALL_OK)
{
    //error handling
}
printf("Number of available integration times: %d\n", nIntTimes);
for (int i=0; i<nIntTimes; i++)
{
    printf("%d: %f\n", i, dblIntTimes[i]);
}

```

Returns

Value If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.3.1.10 HIS_RETURN Acquisition_GetNextSensor (ACQDESCPOS * Pos, HACQDESC * phAcqDesc)

You can use this function to iterate through all recognized sensors in the system. for a programming example see the initialization part of the XISL demonstration.

Parameters

<i>Pos</i>	Pointer to an unsigned 4 byte integer that receives informations that are needed for subsequent calls of this function. To receive the acquisition descriptor (HACQDESC) for the first recognized sensor set Pos to NULL. In case of the 64bit library ACQDESCPOS is a void pointer like defined in Acq.h of the driver SDK.
<i>phAcqDesc</i>	Handle of a structure that contains all needed parameters for acquisition (HACQDESC). If you call Acquisition_Init the first time set hAcqDesc to NULL, in subsequent calls use the former returned value.

Returns

If the initialization is successful zero is returned, otherwise the return value is greater. Call Acquisition_GetErrorCode to get extended information.

3.3.1.11 HIS_RETURN Acquisition_Init (HACQDESC * phAcqDesc, DWORD dwBoardType, int nChannelNr, BOOL bEnableIRQ, UINT Rows, UINT Columns, UINT dwSortFlag, BOOL bSelfInit, BOOL bAlwaysOpen)

The Acquisition_Init function initializes the frame grabber board and the corresponding driver. It enables desired hardware interrupts, prepares acquisition threads, defines callback functions to react on acquisition status changes and tests for sufficient memory space for DMA (direct memory access).

Parameters

<i>phAcqDesc</i>	Handle of a structure that contains all needed parameters for acquisition (HACQDESC). If you call Acquisition_Init the first time set hAcqDesc to NULL, in subsequent calls use the former returned value.
<i>dwBoard-Type</i>	This parameter defines on which communication device the sensor is located. Only one type of frame grabber can be used at the same time.
<i>nChannelNr</i>	This parameter defines the device number. Its possible values depend from dwBoardType and the number of the installed components. For instance if you installed 2 frame grabber boards and you want to acquire data from that one, on that the hardware board selector is set to three, set dwChannelNr equal to 3.
<i>bEnableIRQ</i>	If you want to run the acquisition in polling mode set this parameter to zero. If you want to enable hardware interrupts set the parameter to one.
<i>Rows</i>	Number of sensor columns and rows
<i>Columns</i>	Number of sensor columns and rows
<i>dwSortFlag</i>	Depending on the sensor different sorting schemes are needed because the data come in incorrect order from the detector. dwSortFlags can be one of the following values: The sorting is done automatically by XISL during acquisition. The sorting routines are written in machine code and are therefore very fast.
<i>bSelfInit</i>	If bSelfInit is set to true the function retrieves the detector parameters (Rows, Columns, SortFlags) automatically. If bSelfInit is set to false the configuration parameters supplied by Rows, Columns, dwSortFlags are used.
<i>bAlways-Open</i>	If this parameter is TRUE the XISL is capturing the requested communication port regardless if this port is already opened by other processes running on the system. The use of this option is only recommended in debug versions of your applications because it is not possible to free all system resources allocated by another process.

Note

It enables desired hardware interrupts, prepares acquisition threads, defines call-back functions to react on acquisition status changes and tests for sufficient memory space for DMA (direct memory access). Supported interfaces and channel types:

- 1 HIS_BOARD_TYPE_ELTEC The communication interface to the detector is an XRD-FG or XRD-FGe frame grabber.
- 8 HIS_BOARD_TYPE_ELTEC_XRD_FGX The communication interface to the detector is an XRD-FGX frame grabber.
- 16 HIS_BOARD_TYPE_ELTEC_XRD_FGE_Opto The communication interface to the detector is an XRD-FGe Opto frame grabber.
- 32 HIS_BOARD_TYPE_ELTEC_GbIF The detector communicates over GigabitEthernet with the host PC. (For GbIF and XRpad Detectors - Acquisition_GbIF_Init must be used)
- 96 HIS_BOARD_TYPE_ELTEC_EMBEDDED The detector communicates over GigabitEthernet or WLAN with the host PC. (For GbIF and XRpad -

Detectors Acquisition_GbIF_Init must be used) This value is an ored combination of HIS_BOARD_TYPE_ELTEC_GbIF and HIS_BOARD_TYPE_ELTEC_WPE showing that extended Functionality to the GbIF is available. The value will be set automatically after initialization and can be retrieved by - Acquisition_GetCommChannel(..) Sorting Flags:

- HIS_SORT_NOSORT (0x0) no sorting / XRpad series / 0822 / 1622 / 1642 / 1611
- HIS_SORT_QUAD (0x1) RID128
- HIS_SORT_COLUMN (0x2) RID256
- HIS_SORT_COLUMNQUAD (0x3) RID128-400
- HIS_SORT_QUAD_INVERSE (0x4) RID1024-100
- HID_SORT_QUAD_TILE (0x5) RID512-400 A0
- HIS_SORT_QUAD_TILE_INVERSE (0x6) XRD512-400 A1/A2 XRD 0840
- HIS_SORT_QUAD_TILE_INVERSE_SCRAMBLE (0x7) XRD 512-400 E
- HIS_SORT_OCT_TITLE_INVERSE (0x8) XRD 1640 A , XRD 1620 A , XRD 0820
- HIS_SORT_HEX_TILE_INVERSE (11) XRD 1620/21 AM/AN
- HIS_SORT_HEX_CS(12) XRD 1620/40 AN CS
- HIS_SORT_TOP_BOTTOM(15) XRD 4343RF

Returns

HIS_ALL_OK function successful otherwise an error code

3.3.1.12 HIS_RETURN Acquisition_SetCallbacksAndMessages (HACQDESC
hAcqDesc, *HWND hWnd*, *UINT dwErrorMsg*, *UINT dwLoosingFramesMsg*,
void(CALLBACK *lpfnEndFrameCallback)(HACQDESC) , void(CALLBACK
**lpfnEndAcqCallback)(HACQDESC))*

The Acquisition_SetCallbacksAndMessages function defines callback functions to react on acquisition status changes. For a programming example see the initialization part of the XISL demonstration.

Parameters

<i>hAcqDesc</i>	Handle of a structure that contains all needed parameters for acquisition (HACQDESC). If you call Acquisition_Init the first time set hAcqDesc to NULL, in subsequent calls use the former returned value.
<i>hWnd</i>	If the HSL recognizes an end of DMA transfer and it is ready with sorting, it checks if the application called Acquisition_SetReady after re-drawing. If the application did not call the function, an user defined message (dwLoosingFramesMsg) is posted to hWnd for further handling. If an error occurred during acquisition also a user defined message (dwErrorMsg) is posted to hWnd.

<i>dwErrorMsg</i>	Defines a user message that is posted to hWnd if an error occurs during acquisition.
<i>dwLoosing-FramesMsg</i>	Defines a user message that is posted to hWnd if Acquisition_SetReady wasn't called by the application at the end of sorting.
<i>lpfnEnd-Frame-Callback</i>	Defines a function pointer that is called after the XISL did the sorting. The prototype for the function is given by: In this routine you can do corrections, on-line image processing and redrawing of your data images. Be careful with sending messages from this callback to your application. lpfnEndFrameCallback and lpfnEndAcqCallback are called from a separate thread which is dissimilar to the applications main thread. That should cause problems if you send messages to your main thread via SendMessage. If this causes problems use PostMessage instead. If this parameter is set to NULL it is ignored
<i>lpfnEndAcq-Callback</i>	Defines a function pointer that is called after the XISL did the sorting. The prototype for the function is given by: In this routine you can perform any clean up at acquisition end. If this parameter is set to NULL, it is ignored.

Note

The prototype for the function is given by: void CALLBACK OnEndAcqCallback(H-ACQDESC hAcqDesc); In this routine you can perform any clean up at acquisition end. If this parameter is set to NULL, it is ignored.

Returns

HIS_ALL_OK if successful

3.3.1.13 HIS_RETURN Acquisition_SetLogLevel (XislLoggingLevels *xislLogLvl*)

This function is to set the logging level of the Xisl logging.

Parameters

<i>xislLogLvl</i>	Enum value of which logging level to be used (i.e. LEVEL_DEBUG to display Debug logging and below)
-------------------	--

Returns

If the function is successful it returns zero, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.14 * * * * * HIS_RETURN Acquisition_SetLogOutput (const char * *filePath*, BOOL *consoleOnOff*)

*****/

*****/

This function will create a log file based on the file name provided and will enable console logging if desired.

Parameters

<i>filePath</i>	Customizable file name for logging output file. If the filepath parameter is NULL a default logging file will be created.
<i>consoleOn-Off</i>	Parameter used to enable or disable logging output to the console

Returns

If the function is successful it returns zero, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.3.1.15 HIS_RETURN Acquisition_TogglePerformanceLogging (BOOL onOff)

This function is used to toggle the performance logging. The performance logged will measure the time it takes for a particular function call to execute.

Parameters

<i>onOff</i>	Boolean value used to determine whether to turn internal logging on or off
--------------	--

Returns

If the function is successful it returns zero, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.4 Init Functions for GbIF

Functions

- HIS_RETURN [Acquisition_GbIF_CheckNetworkSpeed](#) (HACQDESC hAcqDesc, WORD *wTiming, long *lPacketDelay, long lMaxNetworkLoadPercent)

This function determines which Timing and Packetdelay the Detector can be set for the current system and network configuration. Please note that this function is intended for one detector only. If more than one detector is connected to the network adapter (detectors in LAN), the parameter lMaxNetworkPercent must divided by the number of connected sensors. Since the network load might change during operation this function cannot guarantee optimal performance. Use Acquisition_GbIF_SetPacketDelay and Acquisition_SetCameraMode(..) to apply the determined settings. For X-Rpad detectors a fixed value is returned for LAN/WLAN transmission. (80% network load are recommended)

- HIS_RETURN [Acquisition_GbIF_DiscoverDetectors](#) ()

Discover all NICs for PKI devices via IP broadcast.

- HIS_RETURN [Acquisition_GbIF_DiscoveredDetectorByIndex](#) (long lIndex, GBIF_DEVICE_PARAM *pDevice)

Retrieves a GBIF_DEVICE_PARAM structure by index.

- HIS_RETURN [Acquisition_GbIF_DiscoveredDetectorCount](#) (long *pDeviceCount)

Retrieves the count of the discovered devices.

- HIS_RETURN [Acquisition_GbIF_ForceIP](#) (GBIF_STRING_DATATYPE *cMAC, GBIF_STRING_DATATYPE *cDefIP, GBIF_STRING_DATATYPE *cDefSubNetMask, GBIF_STRING_DATATYPE *cStdGateway)

To configure the device it can be helpful to force the device temporarily to connect with a certain IP _Adress, usually one out of the same subnet and with the same Std-Gateway like the network card of your computer system. With restart of the detector the device will loose the temporary IP and behave as configured (e.g. IP per DHCP or LLA). This function in not available for the XRpad.

- HIS_RETURN [Acquisition_GbIF_GetConnectionSettings](#) (GBIF_STRING_DATATYPE *ucMAC, unsigned long *ulBootOptions, GBIF_STRING_DATATYPE *ucDefIP, GBIF_STRING_DATATYPE *ucDefSubNetMask, GBIF_STRING_DATATYPE *ucStdGateway)

This function retrieves the connection parameters of a GbIF detector.

- HIS_RETURN [Acquisition_GbIF_GetDetectorProperties](#) (HACQDESC hAcqDesc, GBIF_Detector_Properties *pDetectorProperties)

This function fills the GBIF_Detector_Properties structure, which contains permanently stored information of the connected device.

- HIS_RETURN [Acquisition_GbIF_GetDevice](#) (GBIF_STRING_DATATYPE *ucAddress, DWORD dwAddressType, GBIF_DEVICE_PARAM *pDevice)

This function retrieves the device specified by the passed MAC address.

- HIS_RETURN [Acquisition_GbIF_GetDeviceCnt](#) (long *plNrOfboards)

This function retrieves the total number of sensors found in the network by a network broadcast. If more than one network adapter are installed, a broadcast will be performed on all of them.

- HIS_RETURN [Acquisition_GbIF_GetDeviceList](#) (GBIF_DEVICE_PARAM *pGbIF_DEVICE_PARAM, int nDeviceCnt)

This function retrieves a list of GbIF detector devices found by network broadcast. If multiple network adapters are used in the host system, all of them are checked whether GbIF detectors are connected.

- HIS_RETURN [Acquisition_GbIF_GetDeviceParams](#) (HACQDESC hAcqDesc, GBIF_DEVICE_PARAM *pDevice)

This function retrieves the device parameters of a board that has already been opened.

- HIS_RETURN [Acquisition_GbIF_GetFilterDrvState](#) (HACQDESC hAcqDesc)

Returns the Status of the GbIF filter Driver (if installed) otherwise an Error Code.

- HIS_RETURN [Acquisition_GbIF_GetPacketDelay](#) (HACQDESC hAcqDesc, long *lPacketdelay)

Retrieve the InterPacket Delay, which is set for the current data connection.

- HIS_RETURN [Acquisition_GbIF_GetVersion](#) (int *pMajor, int *pMinor, int *pRelease, char *pStrVersion, int iStrLength)

Retrieve version information about the used GbIF library.

- HIS_RETURN [Acquisition_GbIF_Init](#) (HACQDESC *phAcqDesc, int nChannelNr, BOOL bEnableIRQ, UINT uiRows, UINT uiColumns, BOOL bSelfInit, BOOL bAlwaysOpen, long lInitType, GBIF_STRING_DATATYPE *ucAddress)

The function Acquisition_GbIF_Init initializes the Ethernet connected detectors and the corresponding drivers. It prepares acquisition threads, defines callback functions to react on acquisition status changes.

- HIS_RETURN [Acquisition_GbIF_SetConnectionSettings](#) (GBIF_STRING_DATATYPE *cMAC, unsigned long ulBootOptions, GBIF_STRING_DATATYPE *cDefIP, GBIF_STRING_DATATYPE *cDefSubNetMask, GBIF_STRING_DATATYPE *cStdGateway)

This function provides the parameters to configure how the detector connects to the network adapter. This function is not available for the XRpad.

- HIS_RETURN [Acquisition_GbIF_SetDiscoveryTimeout](#) (long timeout)

Sets the discovery timeout in milliseconds.

- HIS_RETURN [Acquisition_GbIF_SetPacketDelay](#) (HACQDESC hAcqDesc, long lPacketdelay)

The Inter-Packet Delay can be set flexibly to balance out the workload of the IP connection between detector and network adapter. It is recommended to be configured depending on the network load. The value can be retrieved by calling Acquisition_GbIF_CheckNetworkSpeed.

- HIS_RETURN [Acquisition_GbIF_SetPortRange](#) (long lStartPort, long lNrOfports)

enables frame acknowledgement

- HIS_RETURN [Acquisition_wpe_GetVersion](#) (int *major, int *minor, int *release, int *build)

This function returns the version of the used wpe200 library if available and loadable.

3.4.1 Function Documentation

3.4.1.1 HIS_RETURN Acquisition_GbIF_CheckNetworkSpeed (HACQDESC hAcqDesc, WORD * wTiming, long * lPacketDelay, long lMaxNetworkLoadPercent)

This function determines which Timing and Packetdelay the Detector can be set for the current system and network configuration. Please note that this function is intended for one detector only. If more than one detector is connected to the network adapter (detectors in LAN), the parameter lMaxNetworkPercent must divided by the number of connected sensors. Since the network load might change during operation this function cannot guarantee optimal performance. Use Acquisition_GbIF_SetPacketDelay and - Acquisition_SetCameraMode(..) to apply the determined settings. For XRpad detectors a fixed value is returned for LAN/WLAN transmission. (80% network load are recommended)

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>wTiming</i>	Pointer to suggested Free Running Timing for actual System - Performance
<i>lPacketDelay</i>	Pointer to calculated max InterPacketDelay for suggested timing
<i>lMax-Network-LoadPercent</i>	Percentage of Network load for which the Packet Delay shall be cheked. To allow a stable data transmissions with some space for packet resend select~80 percent

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.2 HIS_RETURN Acquisition_GbIF_DiscoverDetectors ()

Discover all NICs for PKI devices via IP broadcast.

Note

Use Acquisition_GbIF_DiscoveredDevel_GetGBifConnectionStatusiceCount and - Acquisition_GbIF_DiscoveredDeviceByIndex in subsequent calls to retrieve information about the discovered devices.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.4.1.3 HIS_RETURN Acquisition_GbIF_DiscoveredDetectorByIndex (long lIndex, GBIF_DEVICE_PARAM * pDevice)

Retrieves a [GBIF_DEVICE_PARAM](#) structure by index.

Note

Use Acquisition_GbIF_DiscoverDevices in a preceding call to discover the network.

Parameters

<i>Index</i>	Index of the device to retrieve the parameters. Valid indices are in the range $0 \leq \text{index} \leq (\text{count} - 1)$.
<i>pDevice</i>	Pointer to a GBIF_DEVICE_PARAM structure to retrieve the parameters of the device.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.4.1.4 HIS_RETURN Acquisition_GbIF_DiscoveredDetectorCount (long * *pDeviceCount*)

Retrieves the count of the discovered devices.

Note

Use Acquisition_GbIF_DiscoverDevices in a preceding call to discover the network.

Parameters

<i>pDevice-Count</i>	Pointer to a variable that retrieves the device count.
----------------------	--

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.4.1.5 HIS_RETURN Acquisition_GbIF_ForceIP (GBIF_STRING_DATATYPE * *cMAC*, GBIF_STRING_DATATYPE * *cDefIP*, GBIF_STRING_DATATYPE * *cDefSubNetMask*, GBIF_STRING_DATATYPE * *cStdGateway*)

To configure the device it can be helpful to force the device temporarily to connect with a certain IP _Adress, usually one out of the same subnet and with the same StdGateway like the network card of your computer system. With restart of the detector the device will loose the temporary IP and behave as configured (e.g. IP per DHCP or LLA). This function is not available for the XRpad.

Parameters

<i>cMAC</i>	MAC address of the device to retrieve a temporary IP.
<i>cDefIP</i>	Temporary IP.

<i>cDefSubNetMask</i>	Temporary Subnet Mask.
<i>cStdGateway</i>	Temporary Gateway.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.6 HIS_RETURN Acquisition_GbIF_GetConnectionSettings (
GBIF_STRING_DATATYPE * ucMAC, unsigned long * ulBootOptions,
GBIF_STRING_DATATYPE * ucDefIP, GBIF_STRING_DATATYPE * ucDefSubNetMask,
GBIF_STRING_DATATYPE * ucStdGateway)

This function retrieves the connection parameters of a GbIF detector.

Parameters

<i>ucMAC</i>	MAC address of the device having the requested settings.
<i>ulBootOptions</i>	OR-able flag indicating the type of connection bitwisely
<i>ucDefIP</i>	Retrieves the IP address the device is connected with.
<i>ucDefSubNetMask</i>	Retrieves the Subnet the device is in.
<i>ucStdGateway</i>	Retrieves the Standard Gateway of the connection to the device.

Note

- HIS_GbIF_IP_STATIC - Use Static IP address stored in the sensor
- HIS_GbIF_IP_LLA - Sensor will propose a Local Link Address
- HIS_GbIF_IP_DHCP - Sensor will receive IP address by DHCP server

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.7 HIS_RETURN Acquisition_GbIF_GetDetectorProperties (HACQDESC
hAcqDesc, GBIF_Detector_Properties * pDetectorProperties)

This function fills the [GBIF_Detector_Properties](#) structure, which contains permanently stored information of the connected device.

Parameters

<i>hAcqDesc</i>	Acquisition descriptor structure.
<i>pDetector-Properties</i>	Pointer to DetectorProperties structure. The memory for the object must be allocated before calling this function.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.8 HIS_RETURN Acquisition_GbIF_GetDevice (GBIF_STRING_DATATYPE * *ucAddress*, DWORD *dwAddressType*, GBIF_DEVICE_PARAM * *pDevice*)

This function retrieves the device specified by the passed MAC address.

Parameters

<i>ucAddress</i>	Address (array of 16 characters) to open the specified board. It can represent the MAC address, IP address or device name of the sensor.
<i>dwAddress-Type</i>	To identify of which type the parameter <i>ucAddress</i> is, <i>lInitType</i> can have the following values All values are defined within the file <i>acq.h</i> . IP-Address, MAC-Address and Detector Name can be retrieved by - <i>Acquisition_GbIF_GetDeviceList</i> .
<i>pDevice</i>	Struct which retrieves attributes and configuration of the device defined by <i>ucAddress</i> .

Note

- HIS_GbIF_IP The sensor is selected by the IP-Address passed in *cAddress*
- HIS_GbIF_MAC The sensor is selected by the MAC-Adress passed in *c-Address*
- HIS_GbIF_NAME The sensor is selected by the Detector-Name passed in *cAddress*

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.9 HIS_RETURN Acquisition_GbIF_GetDeviceCnt (long * *pINrOfboards*)

This function retrieves the total number of sensors found in the network by a network broadcast. If more than one network adapter are installed, a broadcast will be performed on all of them.

Parameters

<i>pINr-Ofboards</i>	Retrieves the number of sensors found in the network broadcasting.
----------------------	--

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.10 HIS_RETURN Acquisition_GbIF_GetDeviceList (GBIF_DEVICE_PARAM * pGbIF_DEVICE_PARAM, int nDeviceCnt)

This function retrieves a list of GbIF detector devices found by network broadcast. If multiple network adapters are used in the host system, all of them are checked whether GbIF detectors are connected.

Parameters

<i>pGbIF_DEVICE_PARAM</i>	Returns a list of GBIF_DEVICE_PARAM elements, with nDeviceCnt entries. For that there has to be memory allocated of the size nDeviceCnt*sizeof(GBIF_DEVICE_PARAM) before calling the function
<i>nDeviceCnt</i>	is the number of devices, which are found within the network. This parameter has to be known before calling Acquisition_GbIF_GetDeviceList and has to be passed to the function. It can be retrieved by calling Acquisition_GbIF_GetDeviceCnt.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.11 HIS_RETURN Acquisition_GbIF_GetDeviceParams (HACQDESC hAcqDesc, GBIF_DEVICE_PARAM * pDevice)

This function retrieves the device parameters of a board that has already been opened.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pDevice</i>	Pointer to structure which retrieves attributes and configuration of the device defined by hAcqDesc.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.12 HIS_RETURN Acquisition_GbIF_GetFilterDrvState (HACQDESC hAcqDesc)

Returns the Status of the GbIF filter Driver (if installed) otherwise an Error Code.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
-----------------	---

Returns

If the function is successful it returns the status of the Filter driver

- 1 for active
- -1 for disabled / not installed
- if the function is not included in gbif.dll or the HACQDESC is not valid an ErroCode is returned

Note

The Filter driver is deprecated and can not be used with current library versions

3.4.1.13 HIS_RETURN Acquisition_GbIF_GetPacketDelay (HACQDESC *hAcqDesc*, long * *IPacketdelay*)

Retrieve the InterPacket Delay, which is set for the current data connection.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>IPacketdelay</i>	Retrieves the currently set value for Inter-Packet Delay

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.14 HIS_RETURN Acquisition_GbIF_GetVersion (int * *pMajor*, int * *pMinor*, int * *pRelease*, char * *pStrVersion*, int *iStrLength*)

Retrieve version information about the used GbIF library.

Parameters

<i>pMajor</i>	Pointer to a variable to retrieve the major version number.
<i>pMinor</i>	Pointer to a variable to retrieve the minor version number.
<i>pRelease</i>	Pointer to a variable to retrieve the release number.
<i>pStrVersion</i>	Pointer to an array of characters to retrieve a unique hash id string.
<i>iStrLength</i>	Length of the pStrVersion array.

Note

Since the hash could become 40 characters long (even if it is very unlikely that this will ever happen), we consider a length of 64 characters for the pStrVersion array.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.4.1.15 HIS_RETURN Acquisition_GbIF_Init (HACQDESC * *phAcqDesc*, int *nChannelNr*, BOOL *bEnableIRQ*, UINT *uiRows*, UINT *uiColumns*, BOOL *bSelfInit*, BOOL *bAlwaysOpen*, long *lInitType*, GBIF_STRING_DATATYPE * *ucAddress*)

The function Acquisition_GbIF_Init initializes the Ethernet connected detectors and the corresponding drivers. It prepares acquisition threads, defines callback functions to react on acquisition status changes.

Parameters

<i>phAcqDesc</i>	Handle of a structure that contains all needed parameters for acquisition (HACQDESC). If you call Acquisition_GbIF_Init the first time set hAcqDesc to NULL, in subsequent calls use the former returned value.
<i>nChannelNr</i>	This parameter defines a number to refer to the initialized device. For each GbIF sensor to be initialized an individual channel number has to be assigned. If you try to access multiple sensors using the same channel number, only the first one will be successfully initialized.
<i>bEnableIRQ</i>	To run the acquisition in polling mode set this parameter to zero. To enable hardware interrupts set the parameter to one.
<i>uiRows</i>	Number of rows of the sensor.
<i>uiColumns</i>	Number of columns of the sensor.
<i>bSelfInit</i>	If bSelfInit is set to true the function retrieves the detector parameters (Rows, Columns, SortFlags) automatically. If bSelfInit is set to false the configuration parameters supplied by Rows, Columns, dwSortFlags are used.
<i>bAlwaysOpen</i>	If this parameter is TRUE the XISL is capturing the requested communication port regardless if this port is already opened by other processes running on the system. The use of this option is only recommended in debug versions of your applications because it is not possible to free all system resources allocated by another process.
<i>lInitType</i>	To identify of which type the parameter cAddress is, lInitType can have the following values: <ul style="list-style-type: none"> HIS_GbIF_FIRST_CAM 0 HIS_GbIF_IP 1 HIS_GbIF_MAC 2 HIS_GbIF_NAME 3 All values are defined within the file acq.h. I-P-Address, MAC-Address and Detector Name can be retrieved by Acquisition_GbIF_GetDeviceList
<i>ucAddress</i>	Address (array of 16 characters) to open the specified board. It can represent the MAC address, IP address or device name of the sensor.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.4.1.16 HIS_RETURN Acquisition_GbIF_SetConnectionSettings (
GBIF_STRING_DATATYPE * cMAC, unsigned long ulBootOptions,
GBIF_STRING_DATATYPE * cDefIP, GBIF_STRING_DATATYPE * cDefSubNetMask,
GBIF_STRING_DATATYPE * cStdGateway)

This function provides the parameters to configure how the detector connects to the network adapter. This function is not available for the XRpad.

Parameters

<i>cMAC</i>	MAC address of the device to be configured.
<i>ulBootOptions</i>	<ul style="list-style-type: none"> OR-able flag to set the type of connection bitwisely:
<i>cDefIP</i>	In case ulBootOptions is equal to HIS_GbIF_IP_STATIC, cDefIP can be used to set a new value as static IP. For that, cDefIP has to contain an IP address when calling the function.
<i>cDefSubNetMask</i>	In case ulBootOptions is equal to HIS_GbIF_IP_STATIC, cDefSubNetMask can be used to set a new value as static IP. For that, cDefSubNetMask has to contain an IP address when calling the function.
<i>cStdGateway</i>	In case ulBootOptions is equal to HIS_GbIF_IP_STATIC, cStdGateway can be used to set a new value as static IP. For that, cStdGateway has to contain an IP address when calling the function. When cStdGateway is zero the detector will be able to be initialized by Acquisition_EnumSensors(..)

Note

- HIS_GbIF_IP_STATIC - Use Static IP address stored in the sensor
- HIS_GbIF_IP_LLA - Sensor will propose a Local Link Address
- HIS_GbIF_IP_DHCP - Sensor will receive IP address by DHCP server

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.17 HIS_RETURN Acquisition_GbIF_SetDiscoveryTimeout (long timeout)

Sets the discovery timeout in milliseconds.

The discovery timeout is the time the library waits for responses of devices in the network, after sending out a broadcast. You should increase this value, if you encounter problems by discovering devices. This applies especially for noisy WLAN environments.

The default value is 750 milliseconds.

Parameters

<i>in</i>	<i>timeout</i>	Discovery timeout in milliseconds. Must be a positive value.
-----------	----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.4.1.18 HIS_RETURN Acquisition_GbIF_SetPacketDelay (HACQDESC *hAcqDesc*, long *IPacketdelay*)

The Inter-Packet Delay can be set flexibly to balance out the workload of the IP connection between detector and network adapter. It is recommended to be configured depending on the network load. The value can be retrieved by calling Acquisition_GbIF_CheckNetworkSpeed.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>IPacketdelay</i>	Value for Inter-Packet Delay, which is to be set in the unit TICKS

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.4.1.19 HIS_RETURN Acquisition_GbIF_SetPortRange (long *IStartPort*, long *INrOfports*)

enables frame acknowledgement

Parameters

<i>hAcqDesc</i>	
<i>Enable</i>	

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

acknowledges last frame

Parameters

<i>hAcqDesc</i>	
<i>reserved</i>	

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

set port range to be used with Ethernet detectors to IStartPort, INrOfports

Parameters

<i>IStartPort</i>	Starting Port (should be 1024)
<i>INrOfports</i>	Number of Devices max (should be 16)

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.4.1.20 HIS_RETURN Acquisition_wpe_GetVersion (int * *major*, int * *minor*, int * *release*, int * *build*)

This function returns the version of the used wpe200 library if available and loadable.

Parameters

<i>major</i>	major version
<i>minor</i>	minor version
<i>release</i>	release version
<i>build</i>	build version

Returns

HIS_ALL_OK function successful otherwise an error code

3.5 Acquire and correction functions

Functions

- HIS_RETURN [Acquisition_Abort](#) (HACQDESC hAcqDesc)
This routine aborts a currently running acquisition.
- HIS_RETURN [Acquisition_AbortCurrentFrame](#) (HACQDESC hAcqDesc)
This routine aborts the transfer of the current frame from the detector to the frame grabber. The detector will be reset then and a new frame transfer will be started immediately. The detector should run in the same mode as the image correction files have been obtained. Therefore it is not recommended to use this function during a measurement.
- HIS_RETURN [Acquisition_Acquire_GainImage](#) (HACQDESC hAcqDesc, WORD *pOffsetData, DWORD *pdwGainData, UINT nRows, UINT nColumns, UINT nFrames)
This function acquires nFrames which are all offset corrected by data stored in pOffsetData. After that the gain data are added in a 32 bit buffer and after acquisition the data values are divided by nFrames (averaging). After averaging the data are further processed for subsequent gain correction of image data. The last acquired data at frame end time are available via pGainData. At the end of the acquisition time the gain data are also accessible via pGainData. The offset data are necessary to derive a valid gain image. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.
- HIS_RETURN [Acquisition_Acquire_GainImage_Ex](#) (HACQDESC hAcqDesc, WORD *pOffsetData, DWORD *pGainData, UINT nRows, UINT nCols, UINT nFrames, UINT dwOpt)
This function acquires nFrames which are all offset corrected by data stored in pOffsetData. After that the gain data are added in a 32 bit buffer and after acquisition the data values are divided by nFrames (averaging). After averaging the data are further processed for subsequent gain correction of image data. The last acquired data at frame end time are available via pGainData. At the end of the acquisition time the gain data are also accessible via pGainData. The offset data are necessary to derive a valid gain image. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.
- HIS_RETURN [Acquisition_Acquire_GainImage_Ex_ROI](#) (HACQDESC hAcqDesc, WORD *pOffsetData, DWORD *pdwGainData, UINT nRows, UINT nColumns, UINT nFrames, UINT dwOpt, UINT uiULX, UINT uiULY, UINT uiBRX, UINT uiBRY, UINT uiMode)
This function is similar to the Acquisition_Acquire_GainImage(..)-function. The function provides the possibility to use a well-defined region of interest for the median determination. The median is used to calculate the gain of single pixels. (see Mathematical description of corrections) The function should be used to acquire the gain image when not the whole panel is illuminated.
- HIS_RETURN [Acquisition_Acquire_GainImage_Ex_ROI_PreloadCorr](#) (HACQDESC hAcqDesc, DWORD *pdwGainData, UINT nRows, UINT nColumns, UINT nFrames, UINT dwOpt, UINT uiULX, UINT uiULY, UINT uiBRX, UINT uiBRY, UINT uiMode)

The function provides the same functionality as Acquisition_Acquire_GainImage_EX_ROI(...) except loading the image correction data. Please use Acquisition_SetCorrData_Ex(..) to set the correction data before.

- HIS_RETURN Acquisition_Acquire_GainImage_PreloadCorr (HACQDESC hAcqDesc, DWORD *pGainData, UINT nRows, UINT nCols, UINT nFrames)

The function provides the same functionality as Acquisition_Acquire_Gain_Image(...) except loading the image correction data. Use Acquisition_SetCorrData_Ex before to set Offset Correction data only..*

- HIS_RETURN Acquisition_Acquire_Image (HACQDESC hAcqDesc, UINT dwFrames, UINT dwSkipFrms, UINT dwOpt, unsigned short *pwOffsetData, DWORD *pdwGainData, DWORD *pdwPxICorrList)

This function acquires dwFrames frames and performs offset, gain and pixel corrections automatically. The routine returns immediately. If you want to be informed about frame end or acquisition end, then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

- HIS_RETURN Acquisition_Acquire_Image_Ex (HACQDESC hAcqDesc, UINT dwFrames, UINT dwSkipFrms, UINT dwOpt, unsigned short *pwOffsetData, UINT dwGainFrames, unsigned short *pwGainData, unsigned short *pwGainAvgData, DWORD *pdwGainData, DWORD *pdwPxICorrList)

This function acquires dwFrames frames and performs offset, gain and pixel corrections automatically. The routine returns immediately. If you want to be informed about frame end or acquisition end, then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

- HIS_RETURN Acquisition_Acquire_Image_PreloadCorr (HACQDESC hAcqDesc, UINT dwFrames, UINT dwSkipFrms, UINT dwOpt)

The function provides the same functionality as Acquisition_Acquire_Image(...) except loading the image correction data. Use Acquisition_SetCorrData_Ex before to set the correction data.

- HIS_RETURN Acquisition_Acquire_OffsetImage (HACQDESC hAcqDesc, WORD *pwOffsetData, UINT nRows, UINT nColumns, UINT nFrames)

This function acquires nFrames, adds them in a 32 bit buffer and after acquisition the data values are divided by nFrames (averaging). The last acquired data at frame end time are available via pOffsetData. At the end of the acquisition time the averaged data are also accessible via pOffsetData. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

- HIS_RETURN Acquisition_Acquire_OffsetImage_Ex (HACQDESC hAcqDesc, WORD *pwOffsetData, UINT nRows, UINT nColumns, UINT nFrames, UINT dwOpt)

This function acquires nFrames, adds them in a 32 bit buffer and after acquisition the data values are divided by nFrames (averaging). The last acquired data at frame end time are available via pOffsetData. At the end of the acquisition time the averaged data are also accessible via pOffsetData. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

- HIS_RETURN Acquisition_Acquire_OffsetImage_PreloadCorr (HACQDESC hAcqDesc, WORD *pwOffsetData, UINT nRows, UINT nColumns, UINT nFrames, UINT dwOpt)

The function provides the same functionality as `Acquisition_Acquire_Offset_Image(...)` except loading the image correction data. Use `Acquisition_SetCorrData_Ex` before to all correction data to NULL.

- HIS_RETURN [Acquisition_CreateGainMap](#) (WORD *pGainData, WORD *pGainAVG, int nCount, int nFrame)

This function creates a list of median values for gain correction using an bright offset corr image to be used with the function `Acquisition_Acquire_Image_Ex`, `Acquisition_SetCorrData_Ex` or `Acquisition_DoOffsetGainCorrection_Ex`. One value for each image in the pGainData-sequence.

- HIS_RETURN [Acquisition_CreateGainMap32](#) (unsigned long *pGainData, unsigned long *pGainAVG, int nCount, int nFrame)

This function creates a list of median values to be used with the function `Acquisition_Acquire_Image_Ex`. One value for each image in the pGainData-sequence.

- HIS_RETURN [Acquisition_CreateOnboardPixelMaskFrom16BitPixelMask](#) (unsigned short *uspPixelMaskSrc, DWORD dwRows, DWORD dwColumns, unsigned char *bpOnboardPixelMask)

This function creates an on-board style (i.e. 1 bit/pixel) pixel mask from a standard (16 bit/pixel) pixel mask image to be used on XRpad2 detectors.

- HIS_RETURN [Acquisition_CreatePixelMap](#) (WORD *pwData, int nDataRows, int nDataColumns, int *pCorrList, int *pnCorrListSize)

This function specifies the size of or creates a pixel correction list (depending on pwData) that one can use in `Acquisition_Acquire_Image` or `Acquisition_DoPixelCorrection`.

- HIS_RETURN [Acquisition_DefineDestBuffers](#) (HACQDESC hAcqDesc, unsigned short *pProcessedData, UINT nFrames, UINT nRows, UINT nColumns)

This function defines the pointers of the destination buffer for `Acquisition_Acquire_Image`. The data are written into this buffer after sorting. The buffer must have a proper size depending on acquisition mode. To acquire one image the buffer must have the size $\text{sensor rows} * \text{sensor columns} * 2$. To acquire a sequence the buffer must have the size $\text{sensor rows} * \text{sensor columns} * 2 * \text{frames}$. In the case of continuous acquisition the buffer must have the size $\text{sensor rows} * \text{sensor columns} * 2 * \text{frames}$ of the ring buffer.

- HIS_RETURN [Acquisition_DoOffsetCorrection](#) (WORD *pSource, WORD *pDest, WORD *pOffsetData, int nCount)

This function performs an offset correction for the data defined in pSource. A suitable place to call this function is the end of frame callback function or the end of acquisition callback.

- HIS_RETURN [Acquisition_DoOffsetGainCorrection](#) (WORD *pSource, WORD *pDest, WORD *pOffsetData, DWORD *pGainData, int nCount)

This function performs an offset and a gain correction for the data defined in pSource at once. A suitable place to call this function is the end of frame callback function or the end acq callback.

- HIS_RETURN [Acquisition_DoOffsetGainCorrection32](#) (unsigned long *pSource, unsigned long *pDest, unsigned long *pOffsetData, unsigned long *pGainData, int nCount)

Similar to `Acquisition_DoOffsetGainCorrection`.

- HIS_RETURN [Acquisition_DoOffsetGainCorrection_Ex](#) (WORD *pSource, WORD *pDest, WORD *pOffsetData, WORD *pGainData, WORD *pGainAVG, int nCount, int nFrame)

This function performs an offset and a multi gain correction using offset corrected bright images as gain for the data defined in pSource at once. A suitable place to call this function is the end of frame callback function or the end of acq callback.

- HIS_RETURN [Acquisition_DoPixelCorrection](#) (WORD *pData, int *pCorrList)

This function performs a defect pixel correction on the data defined by pData using the pCorrList. A suitable place to call this function is the end of frame callback or after the end acq callback.

- HIS_RETURN [Acquisition_GetCorrData](#) (HACQDESC hAcqDesc, unsigned short **ppwOffsetData, DWORD **ppdwGainData, DWORD **ppdwPxLCorrList)

This function returns the addresses of the applied correction buffers.

- HIS_RETURN [Acquisition_GetCorrData_Ex](#) (HACQDESC hAcqDesc, unsigned short **ppwOffsetData, unsigned short **ppwGainData, unsigned short **ppwGainAvgData, UINT **nGainFrames, DWORD **ppdwGainData, DWORD **ppdwPxLCorrList)

This function retrieves the current correction data during a running acquisition.

- HIS_RETURN [Acquisition_GetLatestFrameHeader](#) (HACQDESC hAcqDesc, C-HwHeaderInfo *pInfo, CHwHeaderInfoEx *pInfoEx)

Use this function to retrieve the last acquired frame header of the connected detector. If dwHeaderID in the CHwHeaderInfo structure is 14 pInfoEx will retrieve the extended header, otherwise the structure will be filled with 0xFFFF.

- HIS_RETURN [Acquisition_SetCorrData](#) (HACQDESC hAcqDesc, unsigned short *pwOffsetData, DWORD *pdwGainData, DWORD *pdwPxLCorrList)

This function switches the correction buffers during a running acquisition. It is also possible to switch off corrections by setting the pOffsetData, pGainData and pPixelCorrList to NULL.

- HIS_RETURN [Acquisition_SetCorrData_Ex](#) (HACQDESC hAcqDesc, unsigned short *pwOffsetData, unsigned short *pwGainData, unsigned short *pwGainAvgData, UINT nGainFrames, DWORD *pdwGainData, DWORD *pdwPxLCorrList)

This function switches the correction buffers during a running acquisition. It is also possible to switch off corrections by setting the pOffsetData, pGainData pwGainData, pwGainAvgData, nGainFrames, and pPixelCorrList to NULL.

3.5.1 Function Documentation

3.5.1.1 HIS_RETURN Acquisition_Abort (HACQDESC hAcqDesc)

This routine aborts a currently running acquisition.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
-----------------	--

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.2 HIS_RETURN Acquisition_AbortCurrentFrame (HACQDESC *hAcqDesc*)

This routine aborts the transfer of the current frame from the detector to the frame grabber. The detector will be reset then and a new frame transfer will be started immediately. The detector should run in the same mode as the image correction files have been obtained. Therefore it is not recommended to use this function during a measurement.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
-----------------	--

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.3 HIS_RETURN Acquisition_Acquire_GainImage (HACQDESC *hAcqDesc*, WORD * *pwOffsetData*, DWORD * *pdwGainData*, UINT *nRows*, UINT *nColumns*, UINT *nFrames*)

This function acquires *nFrames* which are all offset corrected by data stored in *pOffsetData*. After that the gain data are added in a 32 bit buffer and after acquisition the data values are divided by *nFrames* (averaging). After averaging the data are further processed for subsequent gain correction of image data. The last acquired data at frame end time are available via *pGainData*. At the end of the acquisition time the gain data are also accessible via *pGainData*. The offset data are necessary to derive a valid gain image. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pwOffsetData</i>	Pointer that contains offset data. (see Acquisition_Acquire_OffsetImage). It is recommended to acquire the Offset shortly before calling Acquisition_Acquire_GainImage
<i>pdwGainData</i>	Pointer to buffer that receives the gain data.
<i>nRows</i>	Number of rows of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nColumns</i>	Number of columns of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nFrames</i>	Number of frames to acquire.

Note

For the optical frame grabbers its recommended to use the Gainsequence correction using offset corrected bright images as a gain (even a single point gain correction is used) since the onboard correction works with offset corrected bright images. This will increase the speed to upload the images to the onboard correction buffer

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.4 HIS_RETURN Acquisition_Acquire_GainImage_Ex (HACQDESC *hAcqDesc*, WORD * *pOffsetData*, DWORD * *pGainData*, UINT *nRows*, UINT *nCols*, UINT *nFrames*, UINT *dwOpt*)

This function acquires *nFrames* which are all offset corrected by data stored in *pOffsetData*. After that the gain data are added in a 32 bit buffer and after acquisition the data values are divided by *nFrames* (averaging). After averaging the data are further processed for subsequent gain correction of image data. The last acquired data at frame end time are available via *pGainData*. At the end of the acquisition time the gain data are also accessible via *pGainData*. The offset data are necessary to derive a valid gain image. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pOffsetData</i>	Pointer that contains offset data. (see Acquisition_Acquire_OffsetImage). It is recommended to acquire the Offset shortly before calling Acquisition_Acquire_GainImage.
<i>pGainData</i>	Pointer to buffer that receives the gain data.
<i>nRows</i>	Number of rows of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nCols</i>	Number of columns of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nFrames</i>	Number of frames to acquire.
<i>dwOpt</i>	Options

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

**3.5.1.5 HIS_RETURN Acquisition_Acquire_GainImage_Ex_ROI (HACQDESC
hAcqDesc, WORD * *pwOffsetData*, DWORD * *pdwGainData*, UINT *nRows*, UINT
nColumns, UINT *nFrames*, UINT *dwOpt*, UINT *uiULX*, UINT *uiULY*, UINT *uiBRX*, UINT
uiBRY, UINT *uiMode*)**

This function is similar to the Acquisition_Acquire_GainImage(..)-function. The function provides the possibility to use a well-defined region of interest for the median determination. The median is used to calculate the gain of single pixels. (see Mathematical description of corrections) The function should be used to acquire the gain image when not the whole panel is illuminated.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>pwOffset-Data</i>	Pointer that contains offset data. (see Acquisition_Acquire_Offset-Image). It is recommended to acquire the Offset shortly before calling Acquisition_Acquire_GainImage
<i>pdwGain-Data</i>	Pointer to buffer that receives the gain data
<i>nRows</i>	Number of rows and columns of the offset data buffer. If the values are not suitable to the current connected sensor the function return with an error.
<i>nColumns</i>	Number of rows and columns of the offset data buffer. If the values are not suitable to the current connected sensor the function return with an error.
<i>nFrames</i>	Number of frames to acquire.
<i>dwOpt</i>	must be 0
<i>uiULX</i>	UpperLeftX define a rect which is used to calculate the Median for the pixel-gain calculation
<i>uiULY</i>	UpperLeftY define a rect which is used to calculate the Median for the pixel-gain calculation
<i>uiBRX</i>	BottomRightX define a rect which is used to calculate the Median for the pixel-gain calculation
<i>uiBRY</i>	BottomRightY define a rect which is used to calculate the Median for the pixel-gain calculation
<i>uiMode</i>	see table

Note

- 0 - normal Gain whole image used for median determination.
- 1 - Median for pixel-gain calculation from ROI (defined by uiULX. . .)
- 2 - Median for pixel-gain calculation from ROI each pixel-gain outside ROI will be set to 1

- 3 - Median for pixel-gain calculation from ROI each pixel-gain outside ROI will be set to 0

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

**/ val 20051012*

3.5.1.6 HIS_RETURN Acquisition_Acquire_GainImage_Ex_ROI_PreloadCorr (HACQDESC *hAcqDesc*, DWORD * *pdwGainData*, UINT *nRows*, UINT *nColumns*, UINT *nFrames*, UINT *dwOpt*, UINT *uiULX*, UINT *uiULY*, UINT *uiBRX*, UINT *uiBRY*, UINT *uiMode*)

The function provides the same functionality as Acquisition_Acquire_GainImage_EX - ROI(...) except loading the image correction data. Please use Acquisition_SetCorrData_Ex(..) to set the correction data before.

Note

Please find the parameter description above

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.7 HIS_RETURN Acquisition_Acquire_GainImage_PreloadCorr (HACQDESC *hAcqDesc*, DWORD * *pGainData*, UINT *nRows*, UINT *nCols*, UINT *nFrames*)

The function provides the same functionality as Acquisition_Acquire_Gain_Image(...) except loading the image correction data. Use Acquisition_SetCorrData_Ex before to set Offset Correction data only.*.

Note

Please find the parameter description at Acquisition_Acquire_GainImage(..)

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.8 HIS_RETURN Acquisition_Acquire_Image (HACQDESC *hAcqDesc*, UINT *dwFrames*, UINT *dwSkipFrms*, UINT *dwOpt*, unsigned short * *pwOffsetData*, DWORD * *pdwGainData*, DWORD * *pdwPxICorrList*)

This function acquires *dwFrames* frames and performs offset, gain and pixel corrections automatically. The routine returns immediately. If you want to be informed about frame end or acquisition end, then define in *Acquisition_Init* the suitable Callback functions and post from there a corresponding message to your application.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwFrames</i>	Number of frames to acquire is one of the sequence options is set for <i>dwOpt</i> . If the continuous option is set this value gives the number of frames in a ring buffer that is used for continuous data acquisition.
<i>dwSkipFrms</i>	Number of frames to skip before a frames is copied into the acquisition buffer.
<i>dwOpt</i>	Options for sequence acquisition.
<i>pwOffset-Data</i>	Pointer that contains offset data. (see <i>Acquisition_Acquire_OffsetImage</i>). The Offset must be actual. It is recommended to acquire them shortly before calling <i>Acquisition_Acquire</i> . If you don't want to perform an offset correction set this parameter to NULL.
<i>pdwGain-Data</i>	Pointer that contains gain data. (see <i>Acquisition_Acquire_GainImage</i>). If you don't want to perform a gain correction set this parameter to NULL.
<i>pdwPxICorr-List</i>	Pointer to a buffer that contains pixel correction data. <i>pdwPixelData</i> points to a linear array of data. Its size is given through ((number of wrong pixels) * 10 + 1) * sizeof(int). The first entry in a group of nine contains the offset of the pixel from the base pointer of the data array. The other eight entries equal the offset of the correction pixels from the base pointer. If you want to use less than eight pixels for correction, then set the remaining entries to -1. The value of the pixel is replaced by the mean value of the correction pixels. If you don't want to perform a pixel correction set this parameter to NULL. An easier way to create a pixel map is the use of the XISL function <i>Acquisition_CreatePixelMap</i> .

Note

- HIS_SEQ_TWO_BUFFERS 0x1 Storage of the sequence into two buffers. Secure image acquisition by separated data transfer and later performed image correction.
- HIS_SEQ_ONE_BUFFER 0x2 Storage of the sequence into one buffer. Direct acquisition and linked correction into one buffer.
- HIS_SEQ_AVERAGE 0x4 All acquired single images are directly added into one buffer and after acquisition divided by the number of frames, including linked correction files.
- HIS_SEQ_DEST_ONE_FRAME 0x8 Sequence of frames using the same image buffer
- HIS_SEQ_COLLATE 0x10 Skip frames after acquiring frames in a ring buffer

- HIS_SEQ_CONTINUOUS 0x100 Continuous acquisition Frames are continuously acquired into a ring buffer of dwFrames

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.9 HIS_RETURN Acquisition_Acquire_Image_Ex (HACQDESC hAcqDesc, UINT dwFrames, UINT dwSkipFrms, UINT dwOpt, unsigned short * pwOffsetData, UINT dwGainFrames, unsigned short * pwGainData, unsigned short * pwGainAvgData, DWORD * pdwGainData, DWORD * pdwPxlCorrList)

This function acquires dwFrames frames and performs offset, gain and pixel corrections automatically. The routine returns immediately. If you want to be informed about frame end or acquisition end, then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwFrames</i>	Number of frames to acquire is one of the sequence options is set for dwOpt. If the continuous option is set this value gives the number of frames in a ring buffer that is used for continuous data acquisition.
<i>dwSkipFrms</i>	Number of frames to skip before a frames is copied into the acquisition buffer.
<i>dwOpt</i>	Options for sequence acquisition.
<i>pwOffsetData</i>	Pointer that contains offset data. (see Acquisition_Acquire_Offset-Image). The Offset must be actual. It is recommended to acquire them shortly before calling Acquisition_Acquire. If you don't want to perform an offset correction set this parameter to NULL.
<i>dwGainFrames</i>	Number of frames in pwGainData
<i>pwGainData</i>	pointer to the median-list created by Acquisition_CreateGainMap
<i>pwGainAvgData</i>	List of Medians representing the frames of the pwGainData sequence
<i>pdwGainData</i>	Pointer that contains gain data. (see Acquisition_Acquire_GainImage). If you don't want to perform a gain correction set this parameter to NULL
<i>pdwPxlCorrList</i>	Pointer to a buffer that contains pixel correction data.pdwPixelData points to a linear array of data. Its size is given through ((number of wrong pixels) * 10 + 1) * sizeof(int). The first entry in a group of nine contains the offset of the pixel from the base pointer of the data array. The other eight entries equal the offset of the correction pixels from the base pointer. If you want to use less than eight pixels for correction, then set the remaining entries to -1. The value of the pixel is replaced by the mean value of the correction pixels. If you don't want to perform a pixel correction set this parameter to NULL. An easier way to create a pixel map is the use of the XISL function Acquisition_CreatePixelMap.

Note

- **HIS_SEQ_TWO_BUFFERS** 0x1 Storage of the sequence into two buffers. Secure image acquisition by separated data transfer and later performed image correction.
- **HIS_SEQ_ONE_BUFFER** 0x2 Storage of the sequence into one buffer. Direct acquisition and linked correction into one buffer.
- **HIS_SEQ_AVERAGE** 0x4 All acquired single images are directly added into one buffer and after acquisition divided by the number of frames, including linked correction files.
- **HIS_SEQ_DEST_ONE_FRAME** 0x8 Sequence of frames using the same image buffer
- **HIS_SEQ_COLLATE** 0x10 Skip frames after acquiring frames in a ring buffer
- **HIS_SEQ_CONTINUOUS** 0x100 Continuous acquisition Frames are continuously acquired into a ring buffer of dwFrames

Returns

If the function is successful it returns **HIS_ALL_OK**, otherwise an error code. To get extended information call **Acquisition_GetErrorCode**.

3.5.1.10 **HIS_RETURN Acquisition_Acquire_Image_PreloadCorr (HACQDESC hAcqDesc, UINT dwFrames, UINT dwSkipFrms, UINT dwOpt)**

The function provides the same functionality as **Acquisition_Acquire_Image(...)** except loading the image correction data. Use **Acquisition_SetCorrData_Ex** before to set the correction data.

Note

Please find the parameter description at **Acquisition_Acquire_Image(..)**

Returns

If the function is successful it returns **HIS_ALL_OK**, otherwise an error code. To get extended information call **Acquisition_GetErrorCode**.

3.5.1.11 **HIS_RETURN Acquisition_Acquire_OffsetImage (HACQDESC hAcqDesc, WORD * pOffsetData, UINT nRows, UINT nColumns, UINT nFrames)**

This function acquires nFrames, adds them in a 32 bit buffer and after acquisition the data values are divided by nFrames (averaging). The last acquired data at frame end time are available via pOffsetData. At the end of the acquisition time the averaged data are also accessible via pOffsetData. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in **Acquisition_Init** the suitable Callback functions and post from there a corresponding message to your application.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pwOffset-Data</i>	Pointer to a acquisition buffer for offset data.
<i>nRows</i>	Number of rows of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nColumns</i>	Number of columns of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nFrames</i>	Number of frames to acquire.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.12 HIS_RETURN Acquisition_Acquire_OffsetImage_Ex (HACQDESC *hAcqDesc*, WORD * *pwOffsetData*, UINT *nRows*, UINT *nColumns*, UINT *nFrames*, UINT *dwOpt*)

This function acquires *nFrames*, adds them in a 32 bit buffer and after acquisition the data values are divided by *nFrames* (averaging). The last acquired data at frame end time are available via *pOffsetData*. At the end of the acquisition time the averaged data are also accessible via *pOffsetData*. The routine returns immediately. If you want to be informed about frame end or acquisition end then define in Acquisition_Init the suitable Callback functions and post from there a corresponding message to your application.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pwOffset-Data</i>	Pointer to a acquisition buffer for offset data.
<i>nRows</i>	Number of rows of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nColumns</i>	Number of columns of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nFrames</i>	Number of frames to acquire.
<i>dwOpt</i>	

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.13 HIS_RETURN Acquisition_Acquire_OffsetImage_PreloadCorr (HACQDESC *hAcqDesc*, WORD * *pwOffsetData*, UINT *nRows*, UINT *nColumns*, UINT *nFrames*, UINT *dwOpt*)

The function provides the same functionality as Acquisition_Acquire_Offset_Image(...) except loading the image correction data. Use Acquisition_SetCorrData_Ex before to all correction data to NULL.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pwOffsetData</i>	Pointer to a acquisition buffer for offset data.
<i>nRows</i>	Number of rows of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nColumns</i>	Number of columns of the offset data buffer. If the value is not suitable to the current connected sensor the function return with an error.
<i>nFrames</i>	Number of frames to acquire.
<i>dwOpt</i>	

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.5.1.14 HIS_RETURN Acquisition_CreateGainMap (WORD * *pGainData*, WORD * *pGainAVG*, int *nCount*, int *nFrame*)

This function creates a list of median values for gain correction using an bright offset corr image to be used with the function Acquisition_Acquire_Image_Ex, Acquisition_SetCorrData_Ex or Acquisition_DoOffsetGainCorrection_Ex. One value for each image in the pGainData-sequence.

Parameters

<i>pGainData</i>	pointer to a sequence of offset corrected data. The images in the array must be sorted by median.
<i>pGainAVG</i>	pointer to a nFrame sized array of type word
<i>nCount</i>	Number of pixels to per image.
<i>nFrame</i>	Number of Frames in pGainData

Returns

If the function is successful it returns true, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

**3.5.1.15 HIS_RETURN Acquisition_CreateGainMap32 (unsigned long * *pGainData*,
unsigned long * *pGainAVG*, int *nCount*, int *nFrame*)**

This function creates a list of median values to be used with the function Acquisition_Acquire_Image_Ex . One value for each image in the pGainData-sequence.

Parameters

<i>pGainData</i>	Pointer to a sequence of offset corrected data. The images in the array must be sorted by median.
<i>pGainAVG</i>	Pointer to a nFrame sized array of type word
<i>nCount</i>	Number of pixels to per image.
<i>nFrame</i>	Number of Frames in pGainData

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.16 HIS_RETURN Acquisition_CreateOnboardPixelMaskFrom16BitPixelMask (unsigned short * *uspPixelMaskSrc*, DWORD *dwRows*, DWORD *dwColumns*, unsigned char * *bpOnboardPixelMask*)

This function creates an on-board style (i.e. 1 bit/pixel) pixel mask from a standard (16 bit/pixel) pixel mask image to be used on XRpad2 detectors.

Parameters

<i>uspPixelMaskSrc</i>	Input (16 bit/pixel) pixel mask image buffer
<i>dwRows</i>	Number of rows in the image
<i>dwColumns</i>	Number of columns in the image
<i>bpOnboardPixelMask</i>	Output, on-board style (1 bit/pixel) pixel mask image buffer

Returns

HIS_ALL_OK

3.5.1.17 HIS_RETURN Acquisition_CreatePixelMap (WORD * *pwData*, int *nDataRows*, int *nDataColumns*, int * *pCorrList*, int * *pnCorrListSize*)

This function specifies the size of or creates a pixel correction list (depending on *pwData*) that one can use in *Acquisition_Acquire_Image* or *Acquisition_DoPixelCorrection*.

Parameters

<i>pwData</i>	Pointer to a data source buffer. Defective pixels are marked by setting their value to 1 (0xFFFF). All other pixel values are recognized as good ones.
<i>nDataRows</i>	Number of rows of the data source.
<i>nDataColumns</i>	Number of columns of the data source.
<i>pCorrList</i>	Pointer to an array (pixel map buffer) that receives the pixel correction data. If <i>pCorrList</i> is set to NULL then only the required size of the pixel map buffer is returned in <i>pnCorrListSize</i> .
<i>pnCorrListSize</i>	Pointer to an integer that receives the required size of the pixel map buffer if <i>pCorrList</i> is set to NULL otherwise it contains the size of the pixel buffer at function call time.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call *Acquisition_GetErrorCode*.

3.5.1.18 HIS_RETURN Acquisition_DefineDestBuffers (HACQDESC *hAcqDesc*, unsigned short * *pProcessedData*, UINT *nFrames*, UINT *nRows*, UINT *nColumns*)

This function defines the pointers of the destination buffer for *Acquisition_Acquire_Image*. The data are written into this buffer after sorting. The buffer must have a proper size depending on acquisition mode. To acquire one image the buffer must have the size *sensor rows * sensor columns * 2*. To acquire a sequence the buffer must have the size *sensor rows * sensor columns * 2 * frames*. In the case of continuous acquisition the buffer must have the size *sensor rows * sensor columns * 2 * frames* of the ring buffer.

Parameters

<i>hAcqDesc</i>	Pointer to HACQDESC.
<i>pProcessedData</i>	Pointer to the destination buffer.
<i>nFrames</i>	Number of frames of the destination buffer. It must be greater than zero.

<i>nRows</i>	Number of rows of the destination buffer. If this number is not suitable to the sensor the function return with an error code. If you need extended information call Acquisition_GetErrorCode.
<i>nColumns</i>	Number of columns of the destination buffer. If this number is not suitable to the sensor the function return with an error code. If you need extended information call Acquisition_GetErrorCode.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.19 HIS_RETURN Acquisition_DoOffsetCorrection (WORD * *pSource*, WORD * *pDest*, WORD * *pOffsetData*, int *nCount*)

This function performs an offset correction for the data defined in *pSource*. A suitable place to call this function is the end of frame callback function or the end of acquisition callback.

Parameters

<i>pSource</i>	Pointer to data source buffer.
<i>pDest</i>	Pointer to data destination buffer. This parameter can be equal to <i>pSource</i> .
<i>pOffsetData</i>	Pointer that contains offset data. (see Acquisition_Acquire_Offset-Image). These data should be lately acquired, so that it is up-to-date. It is recommended to acquire them shortly before calling Acquisition_Do-OffsetCorrection.
<i>nCount</i>	Number of pixels to correct.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.5.1.20 HIS_RETURN Acquisition_DoOffsetGainCorrection (WORD * *pSource*, WORD * *pDest*, WORD * *pOffsetData*, DWORD * *pGainData*, int *nCount*)

This function performs an offset and a gain correction for the data defined in *pSource* at once. A suitable place to call this function is the end of frame callback function or the end acq callback.

Parameters

<i>pSource</i>	Pointer to data source buffer.
<i>pDest</i>	Pointer to data destination buffer. This parameter can be equal to <i>pSource</i> .

<i>pOffsetData</i>	Pointer that contains offset data. (see Acquisition_Acquire_OffsetImage). These data should be lately acquired, so that it is up-to-date. It is recommended to acquire them shortly before calling Acquisition_DoOffsetCorrection.
<i>pGainData</i>	Pointer that contains gain data. (see Acquisition_Acquire_GainImage).
<i>nCount</i>	Number of data entries to correct.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.5.1.21 HIS_RETURN Acquisition_DoOffsetGainCorrection32 (unsigned long * pSource, unsigned long * pDest, unsigned long * pOffsetData, unsigned long * pGainData, int nCount)

Similar to Acquisition_DoOffsetGainCorrection.

3.5.1.22 HIS_RETURN Acquisition_DoOffsetGainCorrection_Ex (WORD * pSource, WORD * pDest, WORD * pOffsetData, WORD * pGainData, WORD * pGainAVG, int nCount, int nFrame)

This function performs an offset and a multi gain correction using offset corrected bright images as gain for the data defined in pSource at once. A suitable place to call this function is the end of frame callback function or the end of acq callback.

Parameters

<i>pSource</i>	Pointer to data source buffer.
<i>pDest</i>	Pointer to data destination buffer. This parameter can be equal to pSource.
<i>pOffsetData</i>	Pointer that contains offset data. (see Acquisition_Acquire_OffsetImage). These data should be lately acquired, so that it is up-to-date. It is recommended to acquire them shortly before calling Acquisition_DoOffsetCorrection.
<i>pGainData</i>	Pointer to a sequence of offset corrected data. The images in the array must be sorted by median.
<i>pGainAVG</i>	Pointer to the median-list created by Acquisition_CreateGainMap
<i>nCount</i>	Number of pixels to correct.
<i>nFrame</i>	Number of frames in pGainData

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.5.1.23 HIS_RETURN Acquisition_DoPixelCorrection (WORD * *pData*, int * *pCorrList*)

This function performs a defect pixel correction on the data defined by *pData* using the *pCorrList*. A suitable place to call this function is the end of frame callback or after the end acq callback.

Parameters

<i>pData</i>	Pointer to data.
<i>pCorrList</i>	Pointer that contains correction data. <i>pCorrList</i> points to a linear array of data. Its size is given through ((number of pixels) * 9 + 1). The first entry in a group of nine contains the offset of the pixel from the base pointer of the data array. The other eight entries equal the offset of the correction pixels from the base pointer. If you want to use less than eight pixels for correction, then set the remaining entries to -1. The value of the pixel is replaced by the mean value of the correction pixels. The end of the pixel correction list is indicated by a value of 1 as the last entry in the pixel map.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.24 HIS_RETURN Acquisition_GetCorrData (HACQDESC *hAcqDesc*, unsigned short ** *ppwOffsetData*, DWORD ** *ppdwGainData*, DWORD ** *ppdwPxlCorrList*)

This function returns the addresses of the applied correction buffers.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>ppwOffsetData</i>	Point to offset data.
<i>ppdwGainData</i>	Pointer to gain data.
<i>ppdwPxlCorrList</i>	Pointer to pixel correction list.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.25 HIS_RETURN Acquisition_GetCorrData_Ex (HACQDESC hAcqDesc, unsigned short ** ppwOffsetData, unsigned short ** ppwGainData, unsigned short ** ppwGainAvgData, UINT ** nGainFrames, DWORD ** ppdwGainData, DWORD ** ppdwPxlCorrList)

This function retrieves the current correction data during a running acquisition.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>ppwOffsetData</i>	Point to offset data. (see Acquisition_Acquire_OffsetImage).
<i>ppwGainData</i>	Pointer to a sequence of offset corrected data. The images in the array must be sorted by median
<i>ppwGainAvgData</i>	Pointer to the median-list created by Acquisition_CreateGainMap
<i>nGainFrames</i>	UINT *pointer retrieving the number of frames in pwGainData
<i>ppdwGainData</i>	Pointer that contains gain data (see Acquisition_Acquire_GainImage).
<i>ppdwPxlCorrList</i>	Pointer to a pixel correction list (see Acquisition_DoPixelCorrection)

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.26 HIS_RETURN Acquisition_GetLatestFrameHeader (HACQDESC hAcqDesc, CHwHeaderInfo * pInfo, CHwHeaderInfoEx * pInfoEx)

Use this function to retrieve the last acquired frame header of the connected detector. If dwHeaderID in the [CHwHeaderInfo](#) structure is 14 pInfoEx will retrieve the extended header, otherwise the structure will be filled with 0xFFFF.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pInfo</i>	Pointer to Structure of type CHwHeaderInfo to retrieve the detector's hardware header
<i>pInfoEx</i>	Pointer to Structure of type CHwHeaderInfoEx to retrieve the detector's hardware header when available. Can be NULL

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.27 HIS_RETURN Acquisition_SetCorrData (HACQDESC hAcqDesc, unsigned short * pwOffsetData, DWORD * pdwGainData, DWORD * pdwPxlCorrList)

This function switches the correction buffers during a running acquisition. It is also possible to switch off corrections by setting the pOffsetData, pGainData and pPixelCorrList to NULL.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>pwOffsetData</i>	Point to offset data (see Acquisition_Acquire_OffsetImage).
<i>pdwGainData</i>	Pointer that contains gain data (see Acquisition_Acquire_GainImage).
<i>pdwPxlCorrList</i>	Pointer to a pixel correction list (see Acquisition_DoPixelCorrection).

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.5.1.28 HIS_RETURN Acquisition_SetCorrData_Ex (HACQDESC hAcqDesc, unsigned short * pwOffsetData, unsigned short * pwGainData, unsigned short * pwGainAvgData, UINT nGainFrames, DWORD * pdwGainData, DWORD * pdwPxlCorrList)

This function switches the correction buffers during a running acquisition. It is also possible to switch off corrections by setting the pOffsetData, pGainData pwGainData, pwGainAvgData, nGainFrames, and pPixelCorrList to NULL.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>pwOffsetData</i>	Pointer to offset data (see Acquisition_Acquire_OffsetImage).
<i>pwGainData</i>	Pointer to a sequence of offset corrected data. The images in the array must be sorted by median.
<i>pwGainAvgData</i>	Pointer to the median-list created by Acquisition_CreateGainMap
<i>nGainFrames</i>	number of frames in pwGainData
<i>pdwGainData</i>	Pointer that contains gain data (see Acquisition_Acquire_GainImage).

<i>pdwPxICorr- List</i>	Pointer to a pixel correction list (see Acquisition_DoPixelCorrection)
-----------------------------	--

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6 Setting up the Device

Functions

- HIS_RETURN [Acquisition_GetCameraBinningMode](#) (HACQDESC hAcqDesc, -WORD *wMode)
Use this function to retrieve the detectors binning mode.
- HIS_RETURN [Acquisition_GetCameraFOVMode](#) (HACQDESC hAcqDesc, WORD *wMode)
Use this function to retrieve the detectors FOV mode.
- HIS_RETURN [Acquisition_GetCameraROI](#) (HACQDESC hAcqDesc, unsigned short *usActivateGrp)
This function returns the currently activated Region of Interest in Sectional Readout Mode.
- HIS_RETURN [Acquisition_GetCameraTriggerMode](#) (HACQDESC hAcqDesc, WORD *wMode)
Retreives TriggerMode for Detectors with Header > 14 and Detectortype > 2.
- HIS_RETURN [Acquisition_GetIpAdress](#) (HACQDESC hAcqDesc, const char **ipAddress)
This function retrieves the current ip address of the detector (if any).
- HIS_RETURN [Acquisition_ResetFrameCnt](#) (HACQDESC hAcqDesc)
This function is used to set internal frame counter to zero. Note: If this function is used during active acquisition the detector Readout time may be affected. (Only XISL version > 3-2-0-9, HeaderId 14 and Cameratype 1 and 2)
- HIS_RETURN [Acquisition_SetCameraBinningMode](#) (HACQDESC hAcqDesc, -WORD wMode)
Use this function to set the detectors binning mode. Not all detectors support image binning. Please refer to the detector manual for a table of supported binning modes.
- HIS_RETURN [Acquisition_SetCameraFOVMode](#) (HACQDESC hAcqDesc, WORD wMode)
This function selects the field of view for XRD 4343RF detectors.
- HIS_RETURN [Acquisition_SetCameraGain](#) (HACQDESC hAcqDesc, WORD wMode)
This function can be used to select the gain of the detector.
- HIS_RETURN [Acquisition_SetCameraMode](#) (HACQDESC hAcqDesc, UINT dwMode)
This function sets the acquisition timing mode of the detector. Currently eight fixed frame times of the detector are provided.
- HIS_RETURN [Acquisition_SetCameraROI](#) (HACQDESC hAcqDesc, unsigned short usActivateGrp)
This function selects a Defined Region of Interest for Readout.
- HIS_RETURN [Acquisition_SetCameraTriggerMode](#) (HACQDESC hAcqDesc, WORD wMode)
The function sets the internal trigger scheme for Detectors.
- HIS_RETURN [Acquisition_SetDACOffset](#) (HACQDESC hAcqDesc, WORD wDACOffsetValue)

Use this function to sets the DAC for offset floor level within the detector for XRD 4343 detectors.

- HIS_RETURN [Acquisition_SetDACOffsetBinningFPS](#) (HACQDESC hAcqDesc, - WORD wBinningMode, double dblFps, WORD *pwValueToFPGA)

XRD 4343 only. Use this function to sets the DAC offset value within the detector by passing the expected frames/sec and the used binning mode. This function may only be used for evaluation purpose. Please contact the application team to get information about how to set the DACOffset.

- HIS_RETURN [Acquisition_SetFrameSyncMode](#) (HACQDESC hAcqDesc, DWORD dwMode)

This function sets the synchronization mode of the detector.

- HIS_RETURN [Acquisition_SetFrameSyncTimeMode](#) (HACQDESC hAcqDesc, unsigned int uiMode, unsigned int dwDelayTime)

This function sets the synchronization mode of the detector to triggered mode and sets the delay time for "DDD/AED " triggered mode with defined integration time.

- HIS_RETURN [Acquisition_SetTimerSync](#) (HACQDESC hAcqDesc, DWORD *dwCycleTime)

This function configures the CycleTime for internal triggered mode.

- HIS_RETURN [Acquisition_SetTriggerOutSignalOptions](#) (HACQDESC hAcqDesc, unsigned short usTiggerOutSignalMode, unsigned short usEP_SeqLength, unsigned short usEP_FirstBrightFrm, unsigned short usEP_LastBrightFrm, unsigned short usEP_Delay1, unsigned short usEP_Delay2, unsigned short usDDD_Delay, int iTriggerOnRisingEdgeEnable, int iSaveAsDefault)

This function defines behavior of the 'TrigOut' - signal of the detector trigger connector and defines the exposure delay.

3.6.1 Function Documentation

3.6.1.1 HIS_RETURN Acquisition_GetCameraBinningMode (HACQDESC hAcqDesc, WORD * wMode)

Use this function to retrieve the detectors binning mode.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	Pointer to value to retrieve the actual binning mode.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.2 HIS_RETURN Acquisition_GetCameraFOVMode (HACQDESC hAcqDesc, WORD * wMode)

Use this function to retrieve the detectors FOV mode.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	Pointer to value to retrieve the actual Field Of View mode.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.3 HIS_RETURN Acquisition_GetCameraROI (HACQDESC *hAcqDesc*, unsigned short * *usActivateGrp*)

This function returns the currently activated Region of Interest in Sectional Readout Mode.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>usActivate-Grp</i>	pointer to unsigend short value representing a bitwise Selection of - Groups to Readout and Transmit

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.4 HIS_RETURN Acquisition_GetCameraTriggerMode (HACQDESC *hAcqDesc*, WORD * *wMode*)

Retreives TriggerMode for Detectors with Header >14 and Detectortype >2.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	Command 42: <ul style="list-style-type: none"> • 0 - DDD • 1 - DDD without clearance scan • 2 - Start Stop • [3] - Trigger Frames (Triggern 1.Art)

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.5 HIS_RETURN Acquisition_GetIpAddress (HACQDESC hAcqDesc, const char ** ipAddress)

This function retrieves the current ip address of the detector (if any).

Parameters

in	<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
out	<i>ipAddress</i>	Pointer to character array. This will be point to a statically allocated chunk of memory, that may be invalidated after subsequent or concurrent calls.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.6.1.6 HIS_RETURN Acquisition_ResetFrameCnt (HACQDESC hAcqDesc)

This function is used to set internal frame counter to zero. Note: If this function is used during active acquisition the detector Readout time may be affected. (Only XISL version > 3-2-0-9, HeaderId 14 and Cameratype 1 and 2)

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
-----------------	---

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.7 HIS_RETURN Acquisition_SetCameraBinningMode (HACQDESC hAcqDesc, WORD wMode)

Use this function to set the detectors binning mode. Not all detectors support image binning. Please refer to the detector manual for a table of supported binning modes.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	Binning Mode to be set <ul style="list-style-type: none"> • value 1 : no binning (default) • value 2 : 2x2 binning • value 3 : 4x4 binning (3x3 R&F) • value 4 : 1x2 binning • value 5 : 1x4 binning

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.6.1.8 HIS_RETURN Acquisition_SetCameraFOVMode (HACQDESC *hAcqDesc*, WORD *wMode*)

This function selects the field of view for XRD 4343RF detectors.

2013-04-22 Val GetTriggerStatus GblF

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	Binning Mode to be set to value, see table

Note

- value 1 : no FOV (default) 2880 x 2880 binning1 1440 x 1440 binning2 960 x 960 binning3
- value 2 : 1920 x 1920 binning1 960 x 960 binning2 640 x 640 binning3
- value 3 : 1440 x 1440 binning1 720 x 736 binning2 480 x 480 binning3
- value 4 : 1440 x 2880 binning1 720 x 1440 binning2 480 x 960 binning3
- value 5 : 480 x 2880 binning1 240 x 1440 binning2 160 x 960 binning3
- Refer to manual for a FOV mode table.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.9 HIS_RETURN Acquisition_SetCameraGain (HACQDESC *hAcqDesc*, WORD *wMode*)

This function can be used to select the gain of the detector.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	Gain factor to set. For the XRpad and XRD 4343RF please refer to the detector reference manual. For the AM-Type the values of all capacities are added. All bitwise combinations are valid. For example: 3 => 1.3p-F. For the xN/xO/xP-Type the Value in the table is set when the detector provides the functionality (refer to detector specification).

Note

Detector xN/xO/xP

- 0 0.25pF
- 1 0.5pF
- 2 1 pF
- 3 2 pF
- 4 4 pF
- 5 8 pF

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.10 HIS_RETURN Acquisition_SetCameraMode (HACQDESC *hAcqDesc*, UINT *dwMode*)

This function sets the acquisition timing mode of the detector. Currently eight fixed frame times of the detector are provided.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwMode</i>	Must be a number between 0 and 7. The corresponding frame time depends on the used detector setting.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.6.1.11 HIS_RETURN Acquisition_SetCameraROI (HACQDESC *hAcqDesc*, unsigned short *usActivateGrp*)

This function selects a Defined Region of Interest for Readout.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>usActivate-Grp</i>	Bitwise Selection of Groups to Readout and Transmit

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.12 HIS_RETURN Acquisition_SetCameraTriggerMode (HACQDESC *hAcqDesc*, WORD *wMode*)

The function sets the internal trigger scheme for Detectors.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wMode</i>	see table

Note

Command 42:

- 0 - WM / DDD
- 1 - WM / DDD without clearance scan
- 2 - Start Stop
- [3] - Trigger Frames (Triggern 1.Art)
- 4 - AutoTrigger Frames
- 5 - Trigger on Row Tag (Framewise with Filter on Trigger input) / RnF frame-Wise flow controlled
- 6 - single shot with post Offset by single trigger (second trigger aborts readout of bright image) - XRPAD2 only
- 7 - dual energy with post offset - XRPAD2 only The selected mode will be used once the detector is set to triggered mode using Acquisition_SetFrame-SyncMode(..)

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.6.1.13 HIS_RETURN Acquisition_SetDACOffset (HACQDESC *hAcqDesc*, WORD *wDACOffsetValue*)

Use this function to sets the DAC for offset floor level within the detector for XRD 4343 detectors.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wDACOffsetValue</i>	wADCOffsetValue to be set. This value has to be between 1 and 4095

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.14 HIS_RETURN Acquisition_SetDACOffsetBinningFPS (HACQDESC *hAcqDesc*, WORD *wBinningMode*, double *dblFps*, WORD * *pwValueToFPGA*)

XRD 4343 only. Use this function to sets the DAC offset value within the detector by passing the expected frames/sec and the used binning mode This function may only be used for evaluation purpose. Please contact the application team to get information about how to set the DACOffset.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>wBinningMode</i>	current binning mode
<i>dblFps</i>	number of frames per sec
<i>pwValueToFPGA</i>	pointer to a variable to retrieve the value send to the fpga

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.15 HIS_RETURN Acquisition_SetFrameSyncMode (HACQDESC *hAcqDesc*, DWORD *dwMode*)

This function sets the synchronization mode of the detector.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwMode</i>	This parameter can values see table

Note**dwMode**

- HIS_SYNCMODE_FREE_RUNNING - no sync, the detector will run without trigger / power save mode in XRpad detectors
- HIS_SYNCMODE_EXTERNAL_TRIGGER - the detector waits for an external signal for readout
- HIS_SYNCMODE_INTERNAL_TIMER - the detector will be triggered by an internal generator (interval to be set with Acquisition_SetTimerSync)

- HIS_SYNCMODE_SOFT_TRIGGER - the detector will be triggered by an internal signal generated via software (Acquisition_SetFrameSync)
- HIS_SYNCMODE_AUTO_TRIGGER - the detector will read out an image when xr-ray exposure is detected
- HIS_SYNCMODE_EXTERNAL_TRIGGER_FG - the external trigger port of the framegrabber will be used as trigger input

HIS_SYNCMODE_EXTERNAL_TRIGGER_FG is required by RnF and 1611 detectors to route the external triggersignal attached to the frame grabber to the detector

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.16 HIS_RETURN Acquisition_SetFrameSyncTimeMode (HACQDESC hAcqDesc, unsigned int uiMode, unsigned int dwDelayTime)

This function sets the synchronization mode of the detector to triggered mode and sets the delay time for "DDD/AED " triggered mode with defined integration time.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>uiMode</i>	Timing Mode must be 0 .. 7
<i>dwDelayTime</i>	Additional delay time in milliseconds (can be 0-ms up to (5000 ms-int-Time)

```
// set special triggermode to 'Data Delivered on demand without clearance scan'
Acquisition_SetCameraTriggerMode (hAcqDesc,1);
```

```
// set Framegrabber to Soft_Trigger-Mode
Acquisition_SetFrameSyncMode(hAcqDesc,HIS_SYNCMODE_SOFT_TRIGGER);
```

```
//set detector to timing 0 delay 1sec
Acquisition_SetFrameSyncTimeMode(hAcqDesc,0,1000);
```

```
hevEndAcq = CreateEvent(NULL, FALSE, FALSE, NULL);
if ((nRet=Acquisition_Acquire_Image(hAcqDesc, 1, 0, HIS_SEQ_ONE_BUFFER, NULL,
    NULL, NULL))!=HIS_ALL_OK)
{
    //error handling
}
```

```
// start the readout
Acquisition_SetFrameSync(hAcqDesc);
```


Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.17 HIS_RETURN Acquisition_SetTimerSync (HACQDESC hAcqDesc, DWORD * dwCycleTime)

This function configures the CycleTime for internal triggered mode.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwCycleTime</i>	Pointer to a 32 bit integer that provides the required cycle time in microsec. After returning of the function this parameter contains the realized cycle time. Before calling this function you have to set the frame grabber to a suitable synchronization mode by a call of Acquisition_SetFrameSyncMode. Some frame grabbers can realize synchronization cycles only in discreet steps. Thats why the function returns the realized cycle time in dwCycleTime.

The parameter dwCycleTime should typically be smaller or equal to 5000000 microsec.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.6.1.18 HIS_RETURN Acquisition_SetTriggerOutSignalOptions (HACQDESC hAcqDesc, unsigned short usTiggerOutSignalMode, unsigned short usEP_SeqLength, unsigned short usEP_FirstBrightFrm, unsigned short usEP_LastBrightFrm, unsigned short usEP_Delay1, unsigned short usEP_Delay2, unsigned short usDDD_Delay, int iTriggerOnRisingEdgeEnable, int iSaveAsDefault)

This function defines behavior of the 'TrigOut' - signal of the detector trigger connector and defines the exposure delay.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>usTiggerOutSignalMode</i>	see table
<i>usEP_SeqLength</i>	Defines the Sequence Length for EP mode
<i>usEP_FirstBrightFrm</i>	Defines the First Frame in the Sequence where the TrigOut Signal is activated
<i>usEP_LastBrightFrm</i>	Defines the Last Frame in the Sequence where the TrigOut Signal is activated

<i>usEP_Delay1</i>	Defines the Delay1 from Begin of Frame in the EP_Frames until the TrigOut Signal gets activated
<i>usEP_Delay2</i>	Defines the Delay2 from Begin of Frame in the EP_Frames until the TrigOut Signal gets deactivated
<i>usDDD_Delay</i>	Additional Delay in DDD/AED/Phototimed Mode
<i>iTriggerOnRisingEdgeEnable</i>	<ul style="list-style-type: none"> • 0 -Trigger falling Edge

- 1 -trigger on rising edge

Parameters

<i>iSaveAsDefault</i>	1 - Save usTiggerOutSignalMode and iTriggerOnRisingEdgeEnable in EEPROM to be availabel after the next power cycle
-----------------------	--

Note

trigger out signal modes: 0 - FRM_EN_PWM 1 - FRM_EN_PWM_INV 2 - EP 3 - EP_INV 4 - DDD_Pulse 5 - DDD_Pulse_INV 6 - GND 7 - VCC

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.7 Callback function

Functions

- HIS_RETURN [Acquisition_GetAcqData](#) (HACQDESC hAcqDesc, DWORD *dwAcqData)
This routine returns a pointer to value that can be set by Acquisition_SetAcqData. - These two functions are useful to avoid global variables to put through parameters to the end of frame and end of acquisition callback functions setted by Acquisition_Init.
- HIS_RETURN [Acquisition_GetActFrame](#) (HACQDESC hAcqDesc, DWORD *dwActAcqFrame, DWORD *dwActSecBuffFrame)
This function retrieves the current acquisition frames.
- HIS_RETURN [Acquisition_GetReady](#) (HACQDESC hAcqDesc)
This function checks whether the passed Acquisition Descriptor is valid.
- HIS_RETURN [Acquisition_GetWinHandle](#) (HACQDESC hAcqDesc, HWND *hWnd)
This function retrieves the current acquisition window handle if defined before.
- HIS_RETURN [Acquisition_SetAcqData](#) (HACQDESC hAcqDesc, DWORD dwAcqData)
This routine sets a value that can be received with Acquisition_GetAcqData. These two functions are useful to avoid global variables to put through parameters to the end of frame and end of acquisition callback functions setted by Acquisition_Init. You can use this function to pass parameters to your callback routines. Thereby it is possible to avoid global variables. The recommended way of use is to define a data structure/value/class and to pass the address of an instance of it by Acquisition_SetAcqData. This has the advantage that a variety of parameters, which are defined within the structure/value/class, are accessible within the callback functions.
- HIS_RETURN [Acquisition_SetReady](#) (HACQDESC hAcqDesc, BOOL bFlag)
This function must be called when the application finished redrawing of the new acquired data. A good place to call this function is the end of frame callback function defined in Acquisition_Init or the message handler to for the redraw message after redrawing.

3.7.1 Function Documentation

3.7.1.1 HIS_RETURN Acquisition_GetAcqData (HACQDESC hAcqDesc, DWORD * dwAcqData)

This routine returns a pointer to value that can be set by Acquisition_SetAcqData. - These two functions are useful to avoid global variables to put through parameters to the end of frame and end of acquisition callback functions setted by Acquisition_Init.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwAcqData</i>	Data to be received. To put through more than one parameters define a structure with the required parameters and cast the pointer to dwData.

Note

For the 64bit version of the XiSL.dll the return type of this parameter has changed to void*

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.7.1.2 HIS_RETURN Acquisition_GetActFrame (HACQDESC hAcqDesc, DWORD * dwActAcqFrame, DWORD * dwActSecBuffFrame)

This function retrieves the current acquisition frames.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwActAcq-Frame</i>	Actual frame of acquisition buffer. The acquisition buffer is allocated internally by the frame grabber driver and is not accessible externally. However, the index can be used to verify the consistency of image sequences; it runs repeatedly from 1 to 8.
<i>dwActSec-BuffFrame</i>	Actual frame for second buffer that is needed to acquire sequences of averaged images. It is the frame count of the acquisition buffer defined by Acquisition_DefineDestBuffers.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.7.1.3 HIS_RETURN Acquisition_GetReady (HACQDESC hAcqDesc)

This function checks whether the passed Acquisition Descriptor is valid.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
-----------------	--

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.7.1.4 HIS_RETURN Acquisition_GetWinHandle (HACQDESC hAcqDesc, HWND * hWnd)

This function retrieves the current acquisition window handle if defined before.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>hWnd</i>	Retrieves the window handle defined in Acquisition_Init.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.7.1.5 HIS_RETURN Acquisition_SetAcqData (HACQDESC *hAcqDesc*, DWORD *dwAcqData*)

This routine sets a value that can be received with Acquisition_GetAcqData. These two functions are useful to avoid global variables to put through parameters to the end of frame and end of acquisition callback functions setted by Acquisition_Init.

You can use this function to pass parameters to your callback routines. Thereby it is possible to avoid global variables. The recommended way of use is to define a data structure/value/class and to pass the address of an instance of it by Acquisition_SetAcqData. This has the advantage that a variety of parameters, which are defined within the structure/value/class, are accessible within the callback functions.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwAcqData</i>	32bit: Data to be accessible within callback functions. dwData can represent a single value as well as an address. In that case a structure/-value/class has to be defined with the required parameters and cast dwData to the pointer. The value/ address can be retrieved in the End-FrameCallback/ EndAcqCallback by using Acquisition_GetAcqData(..). 64bit: Data to be accessible within callback functions: Pass a pointer to a user defined structure/value/class which can be retrieved in the End-FrameCallback/ EndAcqCallback by using Acquisition_GetAcqData(..)

Note

For the 64bit version of the XiSL.dll the return type of this parameter has changed to void*

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.7.1.6 HIS_RETURN Acquisition_SetReady (HACQDESC *hAcqDesc*, BOOL *bFlag*)

This function must be called when the application finished redrawing of the new acquired data. A good place to call this function is the end of frame callback function defined in Acquisition_Init or the message handler to for the redraw message after redrawing.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>bFlag</i>	Boolean value. Set to zero to signal XISL set redrawing isn't ready, otherwise set to one.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.8 Common error handling

Functions

- HIS_RETURN [Acquisition_GetErrorCode](#) (HACQDESC hAcqDesc, DWORD *dwHISerr, DWORD *dwBoardErr)
The function returns extended information if an error occurred during a XISL function call.
- HIS_RETURN [Acquisition_wpe_GetErrorCode](#) ()
This function retrieves the last error code from the wpe200.dll.
- HIS_RETURN [Acquisition_wpe_GetErrorCodeEx](#) (char *pBuffer, long len)
This function retrieves the last error string from the wpe200.dll.

3.8.1 Function Documentation

3.8.1.1 HIS_RETURN Acquisition_GetErrorCode (HACQDESC hAcqDesc, DWORD * dwHISerr, DWORD * dwBoardErr)

The function returns extended information if an error occurred during a XISL function call.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>dwHISerr</i>	Retrieves an error code regarding the XISL itself.
<i>dwBoardErr</i>	Retrieves an error code regarding the acquisition board. Please consult the corresponding documentation of your data acquisition board.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.8.1.2 HIS_RETURN Acquisition_wpe_GetErrorCode (void)

This function retrieves the last error code from the wpe200.dll.

Returns

Latest error code

3.8.1.3 HIS_RETURN Acquisition_wpe_GetErrorCodeEx (char * pBuffer, long len)

This function retrieves the last error string from the wpe200.dll.

Parameters

<i>pBuffer</i>	Pointer to a character array to receive the latest error message. (char[256] recommended).
<i>len</i>	length of the passed character array

Returns

HIS_ALL_OK when function call was successfull otherwise an error code

3.9 functions provided by Acq.h

Functions

- HIS_RETURN [Acq_wpe_GetSystemInformation](#) (const char *ipAddress, char *buffer, int bufferLen)
This function is used to retrieve system information from the detector.
- HIS_RETURN [Acq_WPE_Init](#) ()
Acq_WPE_Init called when dll is loaded, initialize wpe, global objects, etc...
- HIS_RETURN [Acquisition_ActivateServiceMode](#) (HACQDESC hAcqDesc, BOOL bActivate)
This function activates the Service Mode, which means Service Data is written into acquired images. This makes it easy to provide images for support.
- HIS_RETURN [Acquisition_CreateXISFileInMemory](#) (void *pMemoryFileBuffer, void *pDataBuffer, UINT dwRows, UINT dwColumns, UINT dwFrames, BOOL uiOnboardFileHeader, [XIS_FileType](#) filetype)
This function creates an XIS image file (including header) in the given memory buffer.
- HIS_RETURN [Acquisition_DeleteFile](#) ([XislFileHandle](#) fileHandle)
Acquisition_DeleteFile Deletes a file from a location.
- HIS_RETURN [Acquisition_DoOffsetCorrection32](#) (unsigned long *pSource, unsigned long *pDest, unsigned long *pOffsetData, int nCount)
Similar to Acquisition_DoOffsetCorrection.
- HIS_RETURN [Acquisition_DoOffsetGainCorrection_Ex32](#) (unsigned long *pSource, unsigned long *pDest, unsigned long *pOffsetData, unsigned long *pGainData, unsigned long *pGainAVG, int nCount, int nFrame)
Similar to Acquisition_DoOffsetGainCorrection_Ex.
- HIS_RETURN [Acquisition_GetConnectionStatus](#) (HACQDESC hAcqDesc)
retrieves the status of the current connection (if supported).
- HIS_RETURN [Acquisition_GetDACOffsetFloorValueFromFlash](#) (HACQDESC hAcqDesc, unsigned int uiMode, WORD *pwValue)
This function reads the desired offset value from the detector memory for the selected mode. It must be set using Acquisition_SetDACOffset or Acquisition_SetDACOffsetFloorValueByMode.
- HIS_RETURN [Acquisition_GetDetectorProperties](#) (HACQDESC hAcqDesc, [GBIF_Detector_Properties](#) *pDetectorProperties)
This function retrieved the [GBIF_Detector_Properties](#) structure, which contains permanently stored information of the connected device.
- HIS_RETURN [Acquisition_GetFileInfo](#) ([XislFileHandle](#) fileHandle, [XislFileInfo](#) *pFileInfo)
Acquisition_GetFileInfo Retrieves basic information about a file.
- HIS_RETURN [Acquisition_GetRotationAngle](#) (HACQDESC hAcqDesc, long *lRotAngle)
Get onboard Rotation Setting for FG-E Opto.
- HIS_RETURN [Acquisition_GetTriggerOutStatus](#) (HACQDESC hAcqDesc, int *iTriggerStatus)
This function retrieves the triggerout status from detector.

- HIS_RETURN [Acquisition_GetVersion](#) (int *major, int *minor, int *release, int *build)
reads the version of the xisl dll
- HIS_RETURN [Acquisition_GetXISFileBufferSize](#) (size_t *pFileSize, UINT dwRows, UINT dwColumns, UINT dwFrames, BOOL uiOnboardFileHeader, [XIS_FileType](#) filetype)
This function returns the requiered size(in bytes) of an XIS image file with header information.
- HIS_RETURN [Acquisition_IsAcquiringData](#) (HACQDESC hAcqDesc)
This function tests if the XISL is currently acquiring data.
- HIS_RETURN [Acquisition_LoadFile](#) ([XislFileHandle](#) fileHandle, unsigned char **buffer)
Acquisition_LoadFile Loads a file from several locations.
- HIS_RETURN [Acquisition_LoadXISFileToMemory](#) (const char *filename, void *pMemoryFileBuffer, size_t bufferSize)
This function loads an XIS Image file into the given data buffer.
- HIS_RETURN [Acquisition_SaveFile](#) (const char *filename, void *pImageBuffer, UINT dwRows, UINT dwColumns, UINT dwFrames, BOOL uiOnboardFileHeader, [XIS_FileType](#) filetype)
This function saves the given data buffer to file with an appropriately generated HIS-type header.
- HIS_RETURN [Acquisition_SaveRawData](#) (const char *filename, const unsigned char *buffer, size_t bufferSize)
This function saves the given buffer to file. This will either save a raw data file or an XIS image file if Acquisition_CreateXISFileInMemory is used.
- HIS_RETURN [Acquisition_SetConsoleLogging](#) (BOOL enableConsole)
This function will be used to enable or disable logging to the console window.
- HIS_RETURN [Acquisition_SetDACOffsetFloorValueByMode](#) (HACQDESC hAcqDesc, unsigned int uiMode)
This function changes the detector DAC Offset setting for the selected mode.
- HIS_RETURN [Acquisition_SetDACOffsetFloorValueInFlash](#) (HACQDESC hAcqDesc, unsigned int uiMode, WORD wValue)
This function writes the desired offset value to the detector memory for the selected mode. It must be set using Acquisition_SetDACOffset or Acquisition_SetDACOffsetFloorValueByMode.
- HIS_RETURN [Acquisition_SetDACOffsetFloorValueInFlashInternal](#) (HACQDESC hAcqDesc, unsigned int uiMode, WORD wValue)
This function writes the desired offset value to the detector memory for the selected mode. It must be set using Acquisition_SetDACOffset or Acquisition_SetDACOffsetFloorValueByMode.
- HIS_RETURN [Acquisition_SetFileLogging](#) (const char *filename, BOOL enableLogging)
This function will create a log file based on the file name provided and will enable console logging if desired.
- HIS_RETURN [Acquisition_SetFPGACameraMode](#) (HACQDESC hAcqDesc, [FPGAType](#) FPGACommand, BOOL bInverse)

This function sends an FPGA command to the detector. This functions is for internal use only.

- HIS_RETURN [Acquisition_SetRotationAngle](#) (HACQDESC hAcqDesc, long lRotAngle)

This function activates the onboard image rotation angle on optical frame-grabbers for images up to 2048x2048 pixels.

- HIS_RETURN [Acquisition_wpe_GetVersionNEW](#) (int *major, int *minor, int *release, int *build)

Retrieve version information about the used wpe library.

3.9.1 Function Documentation

3.9.1.1 HIS_RETURN Acq_wpe_GetSystemInformation (const char * ipAddress, char * buffer, int bufferLen)

This function is used to retrieve system information from the detector.

Parameters

<i>ipAddress</i>	IP-Address of the device to control
<i>buffer</i>	buffer to hold the return text
<i>bufferLen</i>	length of the buffer

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.9.1.2 HIS_RETURN Acq_WPE_Init ()

Acq_WPE_Init called when dll is loaded, initialize wpe, global objects, etc...

Returns

HIS_RETURN

3.9.1.3 HIS_RETURN Acquisition_ActivateServiceMode (HACQDESC hAcqDesc, BOOL bActivate)

This function activates the Service Mode, which means Service Data is written into acquired images. This makes it easy to provide images for support.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>bActivate</i>	Show Service Data: bActivate = TRUE. otherwise: bActivate = FALSE.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.9.1.4 HIS_RETURN Acquisition_CreateXISFileInMemory (void * *pMemoryFileBuffer*, void * *pDataBuffer*, UINT *dwRows*, UINT *dwColumns*, UINT *dwFrames*, BOOL *uiOnboardFileHeader*, XIS_FileType *filetype*)

This function creates an XIS image file (including header) in the given memory buffer.

Parameters

<i>pMemory-FileBuffer</i>	Output buffer that will contain the XIS image file (with header)
<i>pDataBuffer</i>	Input image data buffer (without header information)
<i>dwRows</i>	Number of rows in the image
<i>dwColumns</i>	Number of columns in the image
<i>dwFrames</i>	Number of frames in the image
<i>uiOnboard-FileHeader</i>	A boolean value indicating whether the file has the on-board style header (TRUE) or not (FALSE)
<i>filetype</i>	The data-type of the pixels in the image.

Returns

HIS_ALL_OK

3.9.1.5 HIS_RETURN Acquisition_DeleteFile (XisFileHandle *fileHandle*)

Acquisition_DeleteFile Deletes a file from a location.

Parameters

<i>in</i>	<i>fileHandle</i>	A valid file handle.
-----------	-------------------	----------------------

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

3.9.1.6 HIS_RETURN Acquisition_DoOffsetCorrection32 (unsigned long * *pSource*, unsigned long * *pDest*, unsigned long * *pOffsetData*, int *nCount*)

Similar to Acquisition_DoOffsetCorrection.

3.9.1.7 HIS_RETURN Acquisition_DoOffsetGainCorrection_Ex32 (unsigned long * *pSource*, unsigned long * *pDest*, unsigned long * *pOffsetData*, unsigned long * *pGainData*, unsigned long * *pGainAVG*, int *nCount*, int *nFrame*)

Similar to Acquisition_DoOffsetGainCorrection_Ex.

3.9.1.8 HIS_RETURN Acquisition_GetConnectionStatus (HACQDESC *hAcqDesc*)

retrievs the status of the current connection (if supported).

Parameters

in	<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor.
----	-----------------	---

Returns

Returns HIS_ALL_OK on success or an appropriate error code including HIS_ERROR_NO_FPGA_ACK otherwise.

Note

Currently only valid for network connected detectors

3.9.1.9 HIS_RETURN Acquisition_GetDACOffsetFloorValueFromFlash (HACQDESC *hAcqDesc*, unsigned int *uiMode*, WORD * *pwValue*)

This function reads the desired offset value from the detector memory for the selected mode. It must be set using Acquisition_SetDACOffset or Acquisition_SetDACOffsetFloorValueByMode.

Parameters

<i>hAcqDesc</i>	Acquisition descriptor structure.
<i>uiMode</i>	Mode to read (0-63).
<i>pwValue</i>	pwValue Pointer to a value to retrieve the offset value

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

Note

currently only available for 4343RF detectors

3.9.1.10 HIS_RETURN Acquisition_GetDetectorProperties (HACQDESC *hAcqDesc*, GBIF_Detector_Properties * *pDetectorProperties*)

This function retrieved the [GBIF_Detector_Properties](#) structure, which contains permanently stored information of the connected device.

Parameters

<i>hAcqDesc</i>	Acquisition descriptor structure.
<i>pDetectorProperties</i>	Pointer to instance of GBIF_Detector_Properties

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.9.1.11 HIS_RETURN Acquisition_GetFileInfo (XislFileHandle *fileHandle*, XislFileInfo * *fileInfo*)

Acquisition_GetFileInfo Retrieves basic information about a file.

Parameters

in	<i>fileHandle</i>	A valid file handle.
out	<i>fileInfo</i>	If the function succeeds, this parameter retrieves the handle of the missed image file.

Returns

Returns always HIS_ALL_OK. This function cannot fail.

3.9.1.12 HIS_RETURN Acquisition_GetRotationAngle (HACQDESC *hAcqDesc*, long * *IRotAngle*)

Get onboard Rotation Setting for FG-E Opto.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>IRotAngle</i>	Retrieves rotation angle.

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.9.1.13 HIS_RETURN Acquisition_GetTriggerOutStatus (HACQDESC *hAcqDesc*, int * *iTriggerStatus*)

This function retrieves the triggerout status from detector.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>iTriggerStatus</i>	see table

Note

- -1 Error, probably timeout
- 0 trigger in was low
- 1 trigger in was high only available for XRpad detectors

Returns

HIS_ALL_OK trigger status could be retrieved otherwise an error code

3.9.1.14 HIS_RETURN Acquisition_GetVersion (int * *major*, int * *minor*, int * *release*, int * *build*)

reads the version of the xisl dll

Parameters

<i>major</i>	Pointer to a variable to retrieve the major version number.
<i>minor</i>	Pointer to a variable to retrieve the minor version number.
<i>release</i>	Pointer to a variable to retrieve the release number.
<i>build</i>	Pointer to a variable to retrieve the build number.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.9.1.15 HIS_RETURN Acquisition_GetXISFileBufferSize (size_t * *pFileSize*, UINT *dwRows*, UINT *dwColumns*, UINT *dwFrames*, BOOL *uiOnboardFileHeader*, XIS_FileType *filetype*)

This function returns the required size(in bytes) of an XIS image file with header information.

Parameters

<i>pFileSize</i>	Requested file size in bytes
<i>dwRows</i>	Number of rows in the image
<i>dwColumns</i>	Number of columns in the image
<i>dwFrames</i>	Number of frames in the image
<i>uiOnboard-FileHeader</i>	A boolean value indicating whether the file has the on-board style header (TRUE) or not (FALSE)
<i>filetype</i>	The data-type of the pixels in the image.

Returns

HIS_ALL_OK

3.9.1.16 HIS_RETURN Acquisition_IsAcquiringData (HACQDESC hAcqDesc)

This function tests if the XISL is currently acquiring data.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
-----------------	--

Returns

If the function is successful it returns 1 if an acquisition is running - 0 if no acquisition is running otherwise an error code.

3.9.1.17 HIS_RETURN Acquisition_LoadFile (XislFileHandle fileHandle, unsigned char ** buffer)

Acquisition_LoadFile Loads a file from several locations.

Parameters

in	<i>fileHandle</i>	A valid file handle.
out	<i>buffer</i>	Pointer to buffer to retrieve the location of the data.

Returns

Returns HIS_ALL_OK on success or an appropriate error code otherwise.

Note

This function is blocking. Depending on size and location of a file, the calling thread might block for an unpredictable period.

3.9.1.18 HIS_RETURN Acquisition_LoadXISFileToMemory (const char * *filename*, void * *pMemoryFileBuffer*, size_t *bufferSize*)

This function loads an XIS Image file into the given data buffer.

Parameters

<i>filename</i>	Path to the input file
<i>pMemory-FileBuffer</i>	Output image data buffer. This should be appropriately allocated to fit all the image and header data.
<i>bufferSize</i>	The size of the output buffer (in bytes) this should be equal to the size of the image and header.

Returns

HIS_ALL_OK

3.9.1.19 HIS_RETURN Acquisition_SaveFile (const char * *filename*, void * *pImageBuffer*, UINT *dwRows*, UINT *dwColumns*, UINT *dwFrames*, BOOL *uiOnboardFileHeader*, XIS_FileType *filetype*)

This function saves the given data buffer to file with an appropriately generated HIS-type header.

Parameters

<i>filename</i>	Path to the output file
<i>pImage-Buffer</i>	Input image data buffer (without header information)
<i>dwRows</i>	Number of rows in the image
<i>dwColumns</i>	Number of columns in the image
<i>dwFrames</i>	Number of frames in the image
<i>uiOnboard-FileHeader</i>	A boolean value indicating whether the file is to have the on-board style header (TRUE) or not (FALSE)
<i>filetype</i>	The data-type of the pixels in the image.

Returns

HIS_ALL_OK

3.9.1.20 HIS_RETURN Acquisition_SaveRawData (const char * *filename*, const unsigned char * *buffer*, size_t *bufferSize*)

This function saves the given buffer to file. This will either save a raw data file or an XIS image file if Acquisition_CreateXISFileInMemory is used.

Parameters

<i>filename</i>	Path to the output file
<i>buffer</i>	Input image data buffer (with or without header information)
<i>bufferSize</i>	The size of the data buffer (including any headers used) in bytes

Returns

HIS_ALL_OK

3.9.1.21 HIS_RETURN Acquisition_SetConsoleLogging (BOOL *enableConsole*)

This function will be used to enable or disable logging to the console window.

Parameters

<i>enable-Console</i>	Parameter used to enable or disable logging output to the console
-----------------------	---

Returns

If the function is successful it returns zero, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.9.1.22 HIS_RETURN Acquisition_SetDACOffsetFloorValueByMode (HACQDESC *hAcqDesc*, unsigned int *uiMode*)

This function changes the detector DAC Offset setting for the selected mode.

Parameters

<i>hAcqDesc</i>	Acquisition descriptor structure.
<i>uiMode</i>	Mode that should be changed (0-63).

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

Note

currently only available for 4343RF detectors; the value for the selected mode is retrieved from detector flash

3.9.1.23 HIS_RETURN Acquisition_SetDACOffsetFloorValueInFlash (HACQDESC *hAcqDesc*, unsigned int *uiMode*, WORD *wValue*)

This function writes the desired offset value to the detector memory for the selected mode. It must be set using Acquisition_SetDACOffset or Acquisition_SetDACOffset-FloorValueByMode.

Parameters

<i>hAcqDesc</i>	Acquisition descriptor structure.
<i>uiMode</i>	Mode that should be changed (0-63).
<i>wValue</i>	(0-4095) wValue value to set for the mode

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

Note

currently only available for 4343RF detectors; only modes ≥ 32 are changeable

3.9.1.24 HIS_RETURN Acquisition_SetDACOffsetFloorValueInFlashInternal (HACQDESC *hAcqDesc*, unsigned int *uiMode*, WORD *wValue*)

This function writes the desired offset value to the detector memory for the selected mode. It must be set using Acquisition_SetDACOffset or Acquisition_SetDACOffset-FloorValueByMode.

Parameters

<i>hAcqDesc</i>	Acquisition descriptor structure.
<i>uiMode</i>	Mode that should be changed (0-63).
<i>wValue</i>	wValue value to set for the mode

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

Note

currently only available for 4343RF detectors

3.9.1.25 HIS_RETURN Acquisition_SetFileLogging (const char * *filename*, BOOL *enableLogging*)

This function will create a log file based on the file name provided and will enable console logging if desired.

Parameters

<i>filename</i>	Customizable file name for logging output file. If the filepath parameter is NULL a default logging file will be created.
<i>enable-Logging</i>	Parameter used to enable or disable logging output to the console

Returns

If the function is successful it returns zero, otherwise an error code.

3.9.1.26 HIS_RETURN Acquisition_SetFPGACameraMode (HACQDESC *hAcqDesc*, FPGAType *FPGACommand*, BOOL *blInverse*)

This function sends an FPGA command to the detector. This functions is for internal use only.

Parameters

<i>hAcqDesc</i>	Handle of a valid Acquisition Descriptor
<i>FPGA-Command</i>	command to send to the detector
<i>blInverse</i>	do a bitwise invers on the data field

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code.

3.9.1.27 HIS_RETURN Acquisition_SetRotationAngle (HACQDESC *hAcqDesc*, long *lRotAngle*)

This function activates the onboard image rotation angle on optical frame-grabbers for images up to 2048x2048 pixels.

Parameters

<i>hAcqDesc</i>	Pointer to acquisition descriptor structure
<i>lRotAngle</i>	Rotation angle: must be 0, 90 or -90

Returns

If the function is successful it returns HIS_ALL_OK, otherwise an error code. To get extended information call Acquisition_GetErrorCode.

3.9.1.28 HIS_RETURN Acquisition_wpe_GetVersionNEW (int * *major*, int * *minor*, int * *release*, int * *build*)

Retrieve version information about the used wpe library.

Parameters

<i>major</i>	Pointer to a variable to retrieve the major version number.
<i>minor</i>	Pointer to a variable to retrieve the minor version number.
<i>release</i>	Pointer to a variable to retrieve the release number.
<i>build</i>	Pointer to a variable to retrieve the build number.

Returns

Returns HIS_ALL_OK on success or an appropriate error code on failure.

3.10 enumerations provided by Acq.h

Classes

- struct [XRpad_BatteryStatus](#)
- struct [XRpad_ShockEvent](#)
- struct [XRpad_ShockSensorReport](#)
- struct [XRpad_TempSensor](#)
- struct [XRpad_TempSensorReport](#)
- struct [XRpad_VersionInfo](#)

Defines

- #define [HIS_ALL_OK](#) 0
- #define [HIS_ERROR_ABORT](#) 46
- #define [HIS_ERROR_ABORTCURRFRAME](#) 21
- #define [HIS_ERROR_ACKNOWLEDGE_IMAGE](#) 133
- #define [HIS_ERROR_ACQ](#) 43
- #define [HIS_ERROR_ACQ_ALREADY_RUNNING](#) 5
- #define [HIS_ERROR_ACQABORT](#) 12
- #define [HIS_ERROR_ACQUISITION](#) 13
- #define [HIS_ERROR_ALREADY_EXISTS](#) 74
- #define [HIS_ERROR_AVERAGED_LOST](#) 49
- #define [HIS_ERROR_BAD_SORTING_PARAM](#) 50
- #define [HIS_ERROR_BOARDINIT](#) 2
- #define [HIS_ERROR_BUFFERSPACE_NOT_SUFF](#) 29
- #define [HIS_ERROR_CORRBUFFER_INCOMPATIBLE](#) 4
- #define [HIS_ERROR_CREATE_MEMORYMAPPING](#) 35
- #define [HIS_ERROR_CREATE_MUTEX](#) 42
- #define [HIS_ERROR_CURL](#) 65
- #define [HIS_ERROR_DESC_NOT_LOCAL](#) 44
- #define [HIS_ERROR_DOES_NOT_EXIST](#) 75
- #define [HIS_ERROR_EMI_NOT_SET](#) 115
- #define [HIS_ERROR_ENABLE_INTERRUPTS](#) 108
- #define [HIS_ERROR_ENABLE_ONBOARD_GAINOFFSET](#) 68
- #define [HIS_ERROR_ENABLE_ONBOARD_MEAN](#) 67
- #define [HIS_ERROR_ENABLE_ONBOARD_OFFSET](#) 66
- #define [HIS_ERROR_ENABLE_ONBOARD_PREVIEW](#) 69
- #define [HIS_ERROR_FRAME_INV](#) 31
- #define [HIS_ERROR_FUNC_NOTIMPL](#) 40
- #define [HIS_ERROR_GET_CHARGE_MODE](#) 139
- #define [HIS_ERROR_GET_NUM_BOARDS](#) 33
- #define [HIS_ERROR_GET_ONBOARD_OFFSET](#) 64
- #define [HIS_ERROR_GETHWHEADERINFO](#) 22
- #define [HIS_ERROR_GETOSVERSION](#) 16
- #define [HIS_ERROR_HEADER_TIMEOUT](#) 56

- `#define HIS_ERROR_HW_ALREADY_OPEN_BY_ANOTHER_PROCESS` 34
- `#define HIS_ERROR_HW_BOARD_CHANNEL_ALREADY_USED` 146
- `#define HIS_ERROR_HWHEADER_INV` 23
- `#define HIS_ERROR_ILLEGAL_INDEX` 60
- `#define HIS_ERROR_INVALID_FILENAME` 77
- `#define HIS_ERROR_INVALID_FUNC_CALL` 20
- `#define HIS_ERROR_INVALID_HANDLE` 73
- `#define HIS_ERROR_INVALID_PARAM` 45
- `#define HIS_ERROR_INVALIDACQDESC` 7
- `#define HIS_ERROR_INVALIDBUFFERNR` 72
- `#define HIS_ERROR_LOAD_COORECTIONIMAGETOBUFFER` 71
- `#define HIS_ERROR_LOADDRIVER` 39
- `#define HIS_ERROR_MEMORY` 1
- `#define HIS_ERROR_MEMORY_MAPPING` 41
- `#define HIS_ERROR_MISSING_VERSION_INFORMATION` 143
- `#define HIS_ERROR_NO_BOARD_IN_SUBNET` 52
- `#define HIS_ERROR_NO_FPGA_ACK` 57
- `#define HIS_ERROR_NOCAMERA` 3
- `#define HIS_ERROR_NODESC_AVAILABLE` 28
- `#define HIS_ERROR_NOT_DISCOVERED` 62
- `#define HIS_ERROR_NOT_INITIALIZED` 61
- `#define HIS_ERROR_NR_OF_BOARDS_CHANGED` 58
- `#define HIS_ERROR_ONBOARDVGFFAILED` 63
- `#define HIS_ERROR_OPEN_FILE` 76
- `#define HIS_ERROR_READ_DATA` 26
- `#define HIS_ERROR_RETRIEVE_ENHANCED_HEADER` 107
- `#define HIS_ERROR_SET_CHARGE_MODE` 118
- `#define HIS_ERROR_SET_IDLE_TIMEOUT` 117
- `#define HIS_ERROR_SET_IMAGE_TAG` 104
- `#define HIS_ERROR_SET_IMAGE_TAG_LENGTH` 106
- `#define HIS_ERROR_SET_ONBOARD_BINNING` 70
- `#define HIS_ERROR_SET_PROC_SCRIPT` 105
- `#define HIS_ERROR_SETBAUDRATE` 27
- `#define HIS_ERROR_SETCAMERAMODE` 30
- `#define HIS_ERROR_SETDISCOVERYTIMEOUT` 78
- `#define HIS_ERROR_SETEXAMFLAG` 59
- `#define HIS_ERROR_SETFRMSYNC` 17
- `#define HIS_ERROR_SETFRMSYNCMODE` 18
- `#define HIS_ERROR_SETLINETRIG_MODE` 24
- `#define HIS_ERROR_SETTIMERSYNC` 19
- `#define HIS_ERROR_SLOW_SYSTEM` 32
- `#define HIS_ERROR_TIMEOUT` 6
- `#define HIS_ERROR_UNABLE_TO_ACCESS_DETECTOR_FLASH` 55
- `#define HIS_ERROR_UNABLE_TO_CLOSE_BOARD` 54
- `#define HIS_ERROR_UNABLE_TO_OPEN_BOARD` 53
- `#define HIS_ERROR_UNKNOWN_IP_MAC_NAME` 51

- #define [HIS_ERROR_VXD_REGISTER_DMA_ADDRESS](#) 36
- #define [HIS_ERROR_VXD_REGISTER_IRQ](#) 14
- #define [HIS_ERROR_VXD_REGISTER_STAT_ADDR](#) 37
- #define [HIS_ERROR_VXD_REGISTER_STATADR](#) 15
- #define [HIS_ERROR_VXD_UNMASK_IRQ](#) 38
- #define [HIS_ERROR_VXDGETDMAADR](#) 11
- #define [HIS_ERROR_VXDNOTFOUND](#) 8
- #define [HIS_ERROR_VXDNOTOPEN](#) 9
- #define [HIS_ERROR_VXDUNKNOWNERROR](#) 10
- #define [HIS_ERROR_WRITE_DATA](#) 25
- #define [HIS_ERROR_WRONG_CAMERA_MODE](#) 48
- #define [HIS_ERROR_WRONGBOARDSELECT](#) 47
- #define [HIS_ERROR_XRPD_CONNECT](#) 134
- #define [HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EVENT](#) 112
- #define [HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EVENT_CRIT](#) 119
- #define [HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EVENT_WARN](#) 120
- #define [HIS_ERROR_XRPD_FACTORY_RESET_SHOCK_EVENT](#) 121
- #define [HIS_ERROR_XRPD_GET_AUTOPOWERONLOCATIONS](#) 137
- #define [HIS_ERROR_XRPD_GET_CURRENT_VOLTAGE](#) 147
- #define [HIS_ERROR_XRPD_GET_SDCARD_INFO](#) 113
- #define [HIS_ERROR_XRPD_GET_SDCARD_TIMEOUT](#) 142
- #define [HIS_ERROR_XRPD_GET_TEMPERATURE_THRESHOLDS](#) 129
- #define [HIS_ERROR_XRPD_NO_EVENT_INTERFACE](#) 111
- #define [HIS_ERROR_XRPD_NO_EVENTCALLBACK_DEFINED](#) 130
- #define [HIS_ERROR_XRPD_NO_LOCATION](#) 116
- #define [HIS_ERROR_XRPD_NO_NETWORK](#) 122
- #define [HIS_ERROR_XRPD_NOT_CONNECTED](#) 144
- #define [HIS_ERROR_XRPD_REQUEST_POWERSTATE](#) 136
- #define [HIS_ERROR_XRPD_RESEND_ALL_MSG](#) 132
- #define [HIS_ERROR_XRPD_RESET_SHOCK](#) 135
- #define [HIS_ERROR_XRPD_RESET_TEMPERATURE_TIMEOUT](#) 127
- #define [HIS_ERROR_XRPD_SDCARDPERFORMANCE](#) 145
- #define [HIS_ERROR_XRPD_SESSION_ERROR](#) 109
- #define [HIS_ERROR_XRPD_SET_AUTOPOWERONLOCATIONS](#) 138
- #define [HIS_ERROR_XRPD_SET_CPUFREQ_GOVERNOR](#) 148
- #define [HIS_ERROR_XRPD_SET_DATE_TIME](#) 131
- #define [HIS_ERROR_XRPD_SET_EVENT](#) 110
- #define [HIS_ERROR_XRPD_SET_FORCE_FSCK](#) 140
- #define [HIS_ERROR_XRPD_SET_NETWORK](#) 123
- #define [HIS_ERROR_XRPD_SET_PRIVATE_KEY](#) 125
- #define [HIS_ERROR_XRPD_SET_SDCARD_TIMEOUT](#) 141
- #define [HIS_ERROR_XRPD_SET_TEMP_FAKE_MODE](#) 114
- #define [HIS_ERROR_XRPD_SET_TEMPERATURE_THRESHOLDS](#) 128
- #define [HIS_ERROR_XRPD_SET_TEMPERATURE_TIMEOUT](#) 126
- #define [HIS_ERROR_XRPD_VERIFY_GENUINENESS](#) 124

Typedefs

- typedef UINT [ACQDESCPOS](#)
- typedef HANDLE [HACQDESC](#)
- typedef enum [OnboardBinningMode](#) [OnboardBinningMode](#)
- typedef enum [ProcScriptOperation](#) [ProcScriptOperation](#)
- typedef enum [XIS_Acquisition_Event](#) [XIS_Acquisition_Event](#)
- typedef enum [XIS_Battery_Event](#) [XIS_Battery_Event](#)
- typedef enum [XIS_Detector_Event](#) [XIS_Detector_Event](#)
- typedef enum [XIS_Detector_TRIGOUT_SignalMode](#) [XIS_Detector_TRIGOUT_SignalMode](#)
- typedef enum [XIS_DetectorTriggerMode](#) [XIS_DetectorTriggerMode](#)
- typedef enum [XIS_Event](#) [XIS_Event](#)
- typedef enum [XIS_FileType](#) [XIS_FileType](#)
- typedef enum [XIS_Library_Event](#) [XIS_Library_Event](#)
- typedef enum [XIS_Sensor_Event](#) [XIS_Sensor_Event](#)
- typedef enum [XislFileEntryType](#) [XislFileEntryType](#)
- typedef void * [XislFileHandle](#)
- typedef enum [XislFileStorageLocation](#) [XislFileStorageLocation](#)
- typedef void * [XislFtpSession](#)
- typedef enum [XRpad_BatteryHealth](#) [XRpad_BatteryHealth](#)
- typedef enum [XRpad_BatteryPresence](#) [XRpad_BatteryPresence](#)
- typedef struct [XRpad_BatteryStatus](#) [XRpad_BatteryStatus](#)
- typedef enum [XRpad_ChargeMode](#) [XRpad_ChargeMode](#)
- typedef enum [XRpad_DataInterfaceControlEnum](#) [XRpad_DataInterfaceControlEnum](#)
- typedef struct [XRpad_ShockEvent](#) [XRpad_ShockEvent](#)
- typedef struct [XRpad_ShockSensorReport](#) [XRpad_ShockSensorReport](#)
- typedef struct [XRpad_TempSensor](#) [XRpad_TempSensor](#)
- typedef struct [XRpad_TempSensorReport](#) [XRpad_TempSensorReport](#)
- typedef struct [XRpad_VersionInfo](#) [XRpad_VersionInfo](#)

Enumerations

- enum [OnboardBinningMode](#) { [ONBOARDBINNING2x1](#) = 0, [ONBOARDBINNING2x2](#) = 1, [ONBOARDBINNING4x1](#) = 2, [ONBOARDBINNING4x4](#) = 3, [ONBOARDBINNING3x3](#) = 4, [ONBOARDBINNING9to4](#) = 5 }
- enum [ProcScriptOperation](#) { [PREBINNING](#), [PREMEAN](#), [PRESTOREBUFFER](#), [OFFSET](#), [GAIN](#), [MEAN](#), [PREVIEW](#), [BINNING](#), [STOREBUFFER](#), [STORESD](#), [SEND](#) }
- enum [XIS_Acquisition_Event](#) { [XAE_TRIGOUT](#) = 0x00000002, [XAE_READOUT](#) = 0x00000004 }
- enum [XIS_Battery_Event](#) { [XBE_BATTERY_REPORT](#) = 0x00000001, [XBE_BATTERY_WARNING](#) = 0x00000002 }
- enum [XIS_Detector_Event](#) { [XDE_BUFFERS_IN_USE](#) = 0x00000001, [XDE_STORED_IMAGE](#) = 0x00000002, [XDE_DROPPED_IMAGE](#) = 0x00000003 }

- enum `XIS_Detector_TRIGOUT_SignalMode` { `TRIGOUT_SIGNAL_FRM_EN_PWM`, `TRIGOUT_SIGNAL_FRM_EN_PWM_INV`, `TRIGOUT_SIGNAL_EP`, `TRIGOUT_SIGNAL_EP_INV`, `TRIGOUT_SIGNAL_DDD_Pulse`, `TRIGOUT_SIGNAL_DDD_Pulse_INV`, `TRIGOUT_SIGNAL_GND`, `TRIGOUT_SIGNAL_VCC` }
- enum `XIS_DetectorTriggerMode` { `TRIGGERMODE_DDD`, `TRIGGERMODE_DDD_WO_CLEARANCE`, `TRIGGERMODE_STARTSTOP`, `TRIGGERMODE_FRAMEWISE`, `TRIGGERMODE_AED`, `TRIGGERMODE_ROWTAG`, `TRIGGERMODE_DDD_POST_OFFSET`, `TRIGGERMODE_DDD_DUAL_POST_OFFSET` }
- enum `XIS_Event` { `XE_ACQUISITION_EVENT` = 0x00000001, `XE_SENSOR_EVENT` = 0x00000002, `XE_SDCARD_EVENT` = 0x00000004, `XE_BATTERY_EVENT` = 0x00000005, `XE_LOCATION_EVENT` = 0x00000006, `XE_NETWORK_EVENT` = 0x00000007, `XE_DETECTOR_EVENT` = 0x00000008, `XE_LIBRARY_EVENT` = 0x00000009, `XE_SDCARD_FSCK_EVENT` = 0x0000000A }
- enum `XIS_FileType` { `PKI_RESERVED` = 1, `PKI_DOUBLE` = 2, `PKI_SHORT` = 4, `PKI_SIGNED` = 8, `PKI_ERRORMAPONBOARD` = 16, `PKI_LONG` = 32, `PKI_SIGNEDSHORT` = `PKI_SHORT` | `PKI_SIGNED`, `PKI_SIGNEDLONG` = `PKI_LONG` | `PKI_SIGNED`, `PKI_FAULTMASK` = `PKI_LONG` | `PKI_RESERVED` }
- enum `XIS_Library_Event` { `XLE_HIS_ERROR_PACKET_LOSS` = 0x00000001 }
- enum `XIS_Sensor_Event` { `XSE_HALL` = 0x00000001, `XSE_SHOCK` = 0x00000010, `XSE_TEMPERATURE` = 0x00000020, `XSE_TEMPERATURE_BACK_TO_NORMAL` = 0x00000021, `XSE_THERMAL_SHUTDOWN` = 0x00000022 }
- enum `XislFileEntryType` { `XFT_File` = 1, `XFT_Directory` = 2, `XFT_Link` = 4, `XFT_Other` = 0x80000000, `XFT_Any` = 0xFFFFFFFF }
- enum `XislFileStorageLocation` { `XFSL_Local` = 0, `XFSL_FTP` = 1 }
- enum `XislLoggingLevels` { `LEVEL_TRACE` = 0, `LEVEL_DEBUG`, `LEVEL_INFO`, `LEVEL_WARN`, `LEVEL_ERROR`, `LEVEL_FATAL`, `LEVEL_ALL`, `LEVEL_NONE` }
- enum `XRpad_BatteryHealth` { `XRpad_BATTERY_OK` = 0x0000, `XRpad_COMMUNICATION_ERROR` = 0x0001, `XRpad_TERMINATE_DISCHARGE_ALARM` = 0x0002, `XRpad_UNDERVOLTAGE_ALARM` = 0x0004, `XRpad_OVERVOLTAGE_ALARM` = 0x0008, `XRpad_OVERTEMPERATURE_ALARM` = 0x0010, `XRpad_BATTERY_UNKNOWN_ERROR` = 0x0080 }
- enum `XRpad_BatteryPresence` { `XRpad_NO_BATTERY` = 0, `XRpad_BATTERY_INSERTED` = 1, `XRpad_DUMMY_INSERTED` = 2 }
- enum `XRpad_ChargeMode` { `XRpad_NOT_CHARGING` = 0, `XRpad_CHARGING_SLOW` = 1, `XRpad_CHARGING_NORMAL` = 2, `XRpad_CHARGING_FAST` = 3, `XRpad_FULLY_CHARGED` = 4, `XRpad_DISCHARGING` = 5 }
- enum `XRpad_DataInterfaceControlEnum` { `XRpad_DATA_VIA_LAN` = 0, `XRpad_DATA_VIA_WLAN` = 1 }
- enum `XRpad_SystemControlEnum` { `XRpad_SYSTEM_CONTROL_REBOOT` = 0, `XRpad_SYSTEM_CONTROL_RESTART_NETWORK` = 1, `XRpad_SYSTEM_CONTROL_SHUTDOWN` = 2, `XRpad_SYSTEM_CONTROL_SET_DEEP_SLEEP` = 3, `XRpad_SYSTEM_CONTROL_SET_IDLE` = 4 }

3.10.1 Define Documentation

3.10.1.1 **#define HIS_ALL_OK 0**

No error

3.10.1.2 **#define HIS_ERROR_ABORT 46**

Error during abort acquisition function.

3.10.1.3 **#define HIS_ERROR_ABORTCURRFRAME 21**

Aborting current frame failed.

3.10.1.4 **#define HIS_ERROR_ACKNOWLEDGE_IMAGE 133**

Error acknowledging the image.

3.10.1.5 **#define HIS_ERROR_ACQ 43**

Error starting the acquisition.

3.10.1.6 **#define HIS_ERROR_ACQ_ALREADY_RUNNING 5**

Acquisition is already running.

3.10.1.7 **#define HIS_ERROR_ACQABORT 12**

An unexpected acquisition abort occurred.

3.10.1.8 **#define HIS_ERROR_ACQUISITION 13**

error occurred during data acquisition.

3.10.1.9 **#define HIS_ERROR_ALREADY_EXISTS 74**

Error Invalid filename, file already exists.

3.10.1.10 **#define HIS_ERROR_AVERAGED_LOST 49**

The number of images for frame grabber onboard averaging must be 2 to the power of n.

3.10.1.11 #define HIS_ERROR_BAD_SORTING_PARAM 50

Parameter for (onboard) sorting not valid.

3.10.1.12 #define HIS_ERROR_BOARDINIT 2

Unable to initialize board.

3.10.1.13 #define HIS_ERROR_BUFFERSPACE_NOT_SUFF 29

Buffer space not sufficient.

3.10.1.14 #define HIS_ERROR_CORRBUFFER_INCOMPATIBLE 4

Your correction files do not have a proper size.

3.10.1.15 #define HIS_ERROR_CREATE_MEMORYMAPPING 35

Error creating memory mapped file.

3.10.1.16 #define HIS_ERROR_CREATE_MUTEX 42

Could not create Mutex.

3.10.1.17 #define HIS_ERROR_CURL 65

Error CURL.

3.10.1.18 #define HIS_ERROR_DESC_NOT_LOCAL 44

Acquisition descriptor is not local.

3.10.1.19 #define HIS_ERROR_DOES_NOT_EXIST 75

Error Invalid filename type does not exist.

3.10.1.20 #define HIS_ERROR_EMI_NOT_SET 115

Error the requested EMI readout mode was not reported by the detector.

3.10.1.21 `#define HIS_ERROR_ENABLE_INTERRUPTS 108`

Error enabling XRPD interrupts.

3.10.1.22 `#define HIS_ERROR_ENABLE_ONBOARD_GAINOFFSET 68`

Error setting onboard gain corr mode.

3.10.1.23 `#define HIS_ERROR_ENABLE_ONBOARD_MEAN 67`

Error setting onboard mean corr mode.

3.10.1.24 `#define HIS_ERROR_ENABLE_ONBOARD_OFFSET 66`

Error setting onboard offset corr mode.

3.10.1.25 `#define HIS_ERROR_ENABLE_ONBOARD_PREVIEW 69`

Error setting onboard preview mode.

3.10.1.26 `#define HIS_ERROR_FRAME_INV 31`

Frame invalid.

3.10.1.27 `#define HIS_ERROR_FUNC_NOTIMPL 40`

Function is not implemented.

3.10.1.28 `#define HIS_ERROR_GET_CHARGE_MODE 139`

Error retrieving the requested charge mode from the detector.

3.10.1.29 `#define HIS_ERROR_GET_NUM_BOARDS 33`

Error during getting number of boards.

3.10.1.30 `#define HIS_ERROR_GET_ONBOARD_OFFSET 64`

Error getting onboard offset.

3.10.1.31 #define HIS_ERROR_GETHWHEADERINFO 22

Getting hardware header failed.

3.10.1.32 #define HIS_ERROR_GETOSVERSION 16

Getting version of operating system failed.

3.10.1.33 #define HIS_ERROR_HEADER_TIMEOUT 56

No frame header received from Detector.

3.10.1.34 #define HIS_ERROR_HW_ALREADY_OPEN_BY_ANOTHER_PROCESS 34

Communication channel already opened by another process.

3.10.1.35 #define HIS_ERROR_HW_BOARD_CHANNEL_ALREADY_USED 146

Requested channel is already openend.

3.10.1.36 #define HIS_ERROR_HWHEADER_INV 23

Hardware header is invalid.

3.10.1.37 #define HIS_ERROR_ILLEGAL_INDEX 60

Error Function called with an illegal index number.

3.10.1.38 #define HIS_ERROR_INVALID_FILENAME 77

Error Invalid filename for image tag or log file.

3.10.1.39 #define HIS_ERROR_INVALID_FUNC_CALL 20

Invalid function call.

3.10.1.40 #define HIS_ERROR_INVALID_HANDLE 73

Error Invalid SHOCKID.

3.10.1.41 **#define HIS_ERROR_INVALID_PARAM 45**

Invalid Parameter.

3.10.1.42 **#define HIS_ERROR_INVALIDACQDESC 7**

Acquisition descriptor invalid.

3.10.1.43 **#define HIS_ERROR_INVALIDBUFFERNR 72**

Error Invalid pointer/buffer passed as parameter.

3.10.1.44 **#define HIS_ERROR_LOAD_COORECTIONIMAGETOBUFFER 71**

Error Loading image from SD to onboard buffer.

3.10.1.45 **#define HIS_ERROR_LOADDRIVER 39**

Unable to load driver.

3.10.1.46 **#define HIS_ERROR_MEMORY 1**

Memory couldn't be allocated.

3.10.1.47 **#define HIS_ERROR_MEMORY_MAPPING 41**

Unable to create memory mapping.

3.10.1.48 **#define HIS_ERROR_MISSING_VERSION_INFORMATION 143**

Error not connected to on detector XRPD process.

3.10.1.49 **#define HIS_ERROR_NO_BOARD_IN_SUBNET 52**

Detector could not be found in the Subnet.

3.10.1.50 **#define HIS_ERROR_NO_FPGA_ACK 57**

Command not acknowledged.

3.10.1.51 #define HIS_ERROR_NOCAMERA 3

Got a time out. May be no detector present.

3.10.1.52 #define HIS_ERROR_NODESC_AVAILABLE 28

No acquisition descriptor available.

3.10.1.53 #define HIS_ERROR_NOT_DISCOVERED 62

Error No detectors discovered yet.

3.10.1.54 #define HIS_ERROR_NOT_INITIALIZED 61

Error Function or function environment not correctly initialised.

3.10.1.55 #define HIS_ERROR_NR_OF_BOARDS_CHANGED 58

Number of boards within network changed during broadcast.

3.10.1.56 #define HIS_ERROR_ONBOARD AVG FAILED 63

Error onbaord averaging failed.

3.10.1.57 #define HIS_ERROR_OPEN_FILE 76

Error Invalid filename for image tag or log file.

3.10.1.58 #define HIS_ERROR_READ_DATA 26

Reading data failed.

3.10.1.59 #define HIS_ERROR_RETRIEVE_ENHANCED_HEADER 107

Error retrieving the enhanced header.

3.10.1.60 #define HIS_ERROR_SET_CHARGE_MODE 118

Error setting the software requested charge mode.

3.10.1.61 **#define HIS_ERROR_SET_IDLE_TIMEOUT 117**

Error setting the on detector idle timeout.

3.10.1.62 **#define HIS_ERROR_SET_IMAGE_TAG 104**

Error setting the onboard image tag.

3.10.1.63 **#define HIS_ERROR_SET_IMAGE_TAG_LENGTH 106**

Error Image tag length exceeded 128char (including path: autosave/).

3.10.1.64 **#define HIS_ERROR_SET_ONBOARD_BINNING 70**

Error setting onboard binning mode.

3.10.1.65 **#define HIS_ERROR_SET_PROC_SCRIPT 105**

Error setting the onboard process script.

3.10.1.66 **#define HIS_ERROR_SETBAUDRATE 27**

Setting baud rate failed.

3.10.1.67 **#define HIS_ERROR_SETCAMERAMODE 30**

Setting detector mode failed.

3.10.1.68 **#define HIS_ERROR_SETDISCOVERYTIMEOUT 78**

Error setting gbif discovery timeout.

3.10.1.69 **#define HIS_ERROR_SETEXAMFLAG 59**

Unable to set the exam flag.

3.10.1.70 **#define HIS_ERROR_SETFRMSYNC 17**

Can not set frame sync.

3.10.1.71 #define HIS_ERROR_SETFRMSYNCMODE 18

Can not set frame sync mode.

3.10.1.72 #define HIS_ERROR_SETLINETRIG_MODE 24

Setting line trigger mode failed.

3.10.1.73 #define HIS_ERROR_SETTIMERSYNC 19

Can not set timer sync.

3.10.1.74 #define HIS_ERROR_SLOW_SYSTEM 32

System too slow.

3.10.1.75 #define HIS_ERROR_TIMEOUT 6

Got a time out from hardware.

3.10.1.76 #define HIS_ERROR_UNABLE_TO_ACCESS_DETECTOR_FLASH 55

Unable to access the flash memory of Detector.

3.10.1.77 #define HIS_ERROR_UNABLE_TO_CLOSE_BOARD 54

Unable to close connection to Network Detector.

3.10.1.78 #define HIS_ERROR_UNABLE_TO_OPEN_BOARD 53

Unable to open connection to Network Detector.

3.10.1.79 #define HIS_ERROR_UNKNOWN_IP_MAC_NAME 51

Connection to Network Detector cannot be opened due to invalid IP address / MAC / Detector name.

3.10.1.80 #define HIS_ERROR_VXD_REGISTER_DMA_ADDRESS 36

Error registering DMA address.

3.10.1.81 **#define HIS_ERROR_VXD_REGISTER_IRQ 14**

Unable to register interrupt.

3.10.1.82 **#define HIS_ERROR_VXD_REGISTER_STAT_ADDR 37**

Error registering static address.

3.10.1.83 **#define HIS_ERROR_VXD_REGISTER_STATADR 15**

Register status address failed.

3.10.1.84 **#define HIS_ERROR_VXD_UNMASK_IRQ 38**

Unable to unmask interrupt.

3.10.1.85 **#define HIS_ERROR_VXDGETDMAADR 11**

VxD Error: GetDmaAddr failed.

3.10.1.86 **#define HIS_ERROR_VXDNOTFOUND 8**

Unable to find VxD.

3.10.1.87 **#define HIS_ERROR_VXDNOTOPEN 9**

Unable to open VxD.

3.10.1.88 **#define HIS_ERROR_VXDUNKNOWNERROR 10**

Unknown error during VxD loading.

3.10.1.89 **#define HIS_ERROR_WRITE_DATA 25**

Writing data failed.

3.10.1.90 **#define HIS_ERROR_WRONG_CAMERA_MODE 48**

Change of Detector Mode during Acquisition.

3.10.1.91 #define HIS_ERROR_WRONGBOARDSELECT 47

The wrong board is selected.

3.10.1.92 #define HIS_ERROR_XRPD_CONNECT 134

Error connecting to the on detector XRPD process.

3.10.1.93 #define HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EVENT 112

Error creating fake shock events.

3.10.1.94 #define HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EVENT_CRIT 119

Error creating critical level fake shock events.

3.10.1.95 #define HIS_ERROR_XRPD_CREATE_FAKE_SHOCK_EVENT_WARN 120

Error creating warning level fake shock events.

3.10.1.96 #define HIS_ERROR_XRPD_FACTORY_RESET_SHOCK_EVENT 121

Error resetting the shock events to factory values.

3.10.1.97 #define HIS_ERROR_XRPD_GET_AUTOPOWERONLOCATIONS 137

Error retrieving the auto power on locations from the detector.

3.10.1.98 #define HIS_ERROR_XRPD_GET_CURRENT_VOLTAGE 147

Error retrieving the Voltage or Current from detector.

3.10.1.99 #define HIS_ERROR_XRPD_GET_SDCARD_INFO 113

Error retrieving the sd card info.

3.10.1.100 #define HIS_ERROR_XRPD_GET_SDCARD_TIMEOUT 142

Error getting the sd card timeout from the detector.

3.10.1.101 **#define HIS_ERROR_XRPD_GET_TEMPERATURE_THRESHOLDS 129**

Error getting the temperature thresholds from the detector.

3.10.1.102 **#define HIS_ERROR_XRPD_NO_EVENT_INTERFACE 111**

Error No interface to communicate event messages active.

3.10.1.103 **#define HIS_ERROR_XRPD_NO_EVENTCALLBACK_DEFINED 130**

Error no eventcallback defined for irq messages.

3.10.1.104 **#define HIS_ERROR_XRPD_NO_LOCATION 116**

Error retrieving the location info from the detector.

3.10.1.105 **#define HIS_ERROR_XRPD_NO_NETWORK 122**

Error getting the on detector LAN network speed.

3.10.1.106 **#define HIS_ERROR_XRPD_NOT_CONNECTED 144**

Error not connected to on detector XRPD process.

3.10.1.107 **#define HIS_ERROR_XRPD_REQUEST_POWERSTATE 136**

Error setting the power state on the detector.

3.10.1.108 **#define HIS_ERROR_XRPD_RESEND_ALL_MSG 132**

Error triggering the resend of all current messages by the XRPD.

3.10.1.109 **#define HIS_ERROR_XRPD_RESET_SHOCK 135**

Error resetting the shock event.

3.10.1.110 **#define HIS_ERROR_XRPD_RESET_TEMPERATURE_TIMEOUT 127**

Error resetting the temperature timeout counter.

3.10.1.111 `#define HIS_ERROR_XRPD_SDCARDPERFORMANCE` 145

Error retrieving the SD card performance.

3.10.1.112 `#define HIS_ERROR_XRPD_SESSION_ERROR` 109

Error XRPD session Error.

3.10.1.113 `#define HIS_ERROR_XRPD_SET_AUTOPOWERONLOCATIONS` 138

Error setting the auto power on locations on the detector.

3.10.1.114 `#define HIS_ERROR_XRPD_SET_CPUFREQ_GOVERNOR` 148

Unable to set on detector CPU governor.

3.10.1.115 `#define HIS_ERROR_XRPD_SET_DATE_TIME` 131

Error setting on detectors date and time.

3.10.1.116 `#define HIS_ERROR_XRPD_SET_EVENT` 110

Error No interface to communicate event messages active.

3.10.1.117 `#define HIS_ERROR_XRPD_SET_FORCE_FSCK` 140

Error requesting an fscheck on next boot.

3.10.1.118 `#define HIS_ERROR_XRPD_SET_NETWORK` 123

Error setting the on detector LAN network speed.

3.10.1.119 `#define HIS_ERROR_XRPD_SET_PRIVATE_KEY` 125

Error setting the private key for genuiness.

3.10.1.120 `#define HIS_ERROR_XRPD_SET_SDCARD_TIMEOUT` 141

Error setting the sd card timeout on the detector.

3.10.1.121 **#define HIS_ERROR_XRPD_SET_TEMP_FAKE_MODE 114**

Error activating the fake temperatur mode on the detector.

3.10.1.122 **#define HIS_ERROR_XRPD_SET_TEMPERATURE_THRESHOLDS 128**

Error setting the temperature thresholds on the detector.

3.10.1.123 **#define HIS_ERROR_XRPD_SET_TEMPERATURE_TIMEOUT 126**

Error setting the temperature timeout.

3.10.1.124 **#define HIS_ERROR_XRPD_VERIFY_GENUINENESS 124**

Error verifying the private key for genuiness.

3.10.2 Typedef Documentation

3.10.2.1 **typedef UINT ACQDESCPOS**

3.10.2.2 **typedef HANDLE HACQDESC**

AcquisitionDesc defines a data structure that is used by all functions of XISL. It contains all required parameters for the acquisition. Access to the data fields is only possible via the XISL API functions. HACQDESC defines a HANLDE to the acquisition descriptor.

See also

AcquisitionDesc

3.10.2.3 **typedef enum OnboardBinningMode OnboardBinningMode**

3.10.2.4 **typedef enum ProcScriptOperation ProcScriptOperation**

3.10.2.5 **typedef enum XIS_Acquisition_Event XIS_Acquisition_Event**

3.10.2.6 **typedef enum XIS_Battery_Event XIS_Battery_Event**

3.10.2.7 **typedef enum XIS_Detector_Event XIS_Detector_Event**

3.10.2.8 **typedef enum XIS_Detector_TRIGOUT_SignalMode
XIS_Detector_TRIGOUT_SignalMode**

3.10.2.9 **typedef enum XIS_DetectorTriggerMode XIS_DetectorTriggerMode**

- 3.10.2.10 `typedef enum XIS_Event XIS_Event`
- 3.10.2.11 `typedef enum XIS_FileType XIS_FileType`
- 3.10.2.12 `typedef enum XIS_Library_Event XIS_Library_Event`
- 3.10.2.13 `typedef enum XIS_Sensor_Event XIS_Sensor_Event`
- 3.10.2.14 `typedef enum XislFileEntryType XislFileEntryType`
- 3.10.2.15 `typedef void* XislFileHandle`
- 3.10.2.16 `typedef enum XislFileStorageLocation XislFileStorageLocation`
- 3.10.2.17 `typedef void* XislFtpSession`
- 3.10.2.18 `typedef enum XRpad_BatteryHealth XRpad_BatteryHealth`
- 3.10.2.19 `typedef enum XRpad_BatteryPresence XRpad_BatteryPresence`
- 3.10.2.20 `typedef struct XRpad_BatteryStatus XRpad_BatteryStatus`
- 3.10.2.21 `typedef enum XRpad_ChargeMode XRpad_ChargeMode`
- 3.10.2.22 `typedef enum XRpad_DataInterfaceControlEnum
XRpad_DataInterfaceControlEnum`

Possible image transfer channels for XRpad

- 3.10.2.23 `typedef struct XRpad_ShockEvent XRpad_ShockEvent`
- 3.10.2.24 `typedef struct XRpad_ShockSensorReport XRpad_ShockSensorReport`
- 3.10.2.25 `typedef struct XRpad_TempSensor XRpad_TempSensor`
- 3.10.2.26 `typedef struct XRpad_TempSensorReport XRpad_TempSensorReport`
- 3.10.2.27 `typedef struct XRpad_VersionInfo XRpad_VersionInfo`

3.10.3 Enumeration Type Documentation

- 3.10.3.1 `enum OnboardBinningMode`

Enumerator:

ONBOARDBINNING2x1
ONBOARDBINNING2x2
ONBOARDBINNING4x1

ONBOARDBINNING4x4

ONBOARDBINNING3x3

ONBOARDBINNING9to4

3.10.3.2 enum ProcScriptOperation

Enumerator:

PREBINNING

PREMEAN

PRESTOREBUFFER

OFFSET

GAIN

MEAN

PREVIEW

BINNING

STOREBUFFER

STORESD

SEND

3.10.3.3 enum XIS_Acquisition_Event

Enumerator:

XAE_TRIGOUT

XAE_READOUT

3.10.3.4 enum XIS_Battery_Event

Enumerator:

XBE_BATTERY_REPORT

XBE_BATTERY_WARNING

3.10.3.5 enum XIS_Detector_Event

Enumerator:

XDE_BUFFERS_IN_USE

XDE_STORED_IMAGE

XDE_DROPPED_IMAGE

3.10.3.6 enum XIS_Detector_TRIGOUT_SignalMode

Enumerator:

TRIGOUT_SIGNAL_FRM_EN_PWM
TRIGOUT_SIGNAL_FRM_EN_PWM_INV
TRIGOUT_SIGNAL_EP
TRIGOUT_SIGNAL_EP_INV
TRIGOUT_SIGNAL_DDD_Pulse
TRIGOUT_SIGNAL_DDD_Pulse_INV
TRIGOUT_SIGNAL_GND
TRIGOUT_SIGNAL_VCC

3.10.3.7 enum XIS_DetectorTriggerMode

Enumerator:

TRIGGERMODE_DDD
TRIGGERMODE_DDD_WO_CLEARANCE
TRIGGERMODE_STARTSTOP
TRIGGERMODE_FRAMEWISE
TRIGGERMODE_AED
TRIGGERMODE_ROWTAG
TRIGGERMODE_DDD_POST_OFFSET
TRIGGERMODE_DDD_DUAL_POST_OFFSET

3.10.3.8 enum XIS_Event

Enumerator:

XE_ACQUISITION_EVENT
XE_SENSOR_EVENT
XE_SDCARD_EVENT
XE_BATTERY_EVENT
XE_LOCATION_EVENT
XE_NETWORK_EVENT
XE_DETECTOR_EVENT
XE_LIBRARY_EVENT
XE_SDCARD_FSCK_EVENT

3.10.3.9 enum XIS_FileType

Enumerator:

PKI_RESERVED
PKI_DOUBLE
PKI_SHORT
PKI_SIGNED
PKI_ERRORMAPONBOARD
PKI_LONG
PKI_SIGNEDSHORT
PKI_SIGNEDLONG
PKI_FAULTMASK

3.10.3.10 enum XIS_Library_Event

Enumerator:

XLE_HIS_ERROR_PACKET_LOSS Frame is lost. not all network packages could be received.

3.10.3.11 enum XIS_Sensor_Event

Enumerator:

XSE_HALL
XSE_SHOCK
XSE_TEMPERATURE
XSE_TEMPERATURE_BACK_TO_NORMAL
XSE_THERMAL_SHUTDOWN

3.10.3.12 enum XisIFileEntryType

Enumerator:

XFT_File
XFT_Directory
XFT_Link
XFT_Other
XFT_Any

3.10.3.13 enum XislFileStorageLocation

Enumerator:

XFSL_Local

XFSL_FTP

3.10.3.14 enum XislLoggingLevels

Enumerator:

LEVEL_TRACE

LEVEL_DEBUG

LEVEL_INFO

LEVEL_WARN

LEVEL_ERROR

LEVEL_FATAL

LEVEL_ALL

LEVEL_NONE

3.10.3.15 enum XRpad_BatteryHealth

Enumerator:

XRpad_BATTERY_OK

XRpad_COMMUNICATION_ERROR

XRpad_TERMINATE_DISCHARGE_ALARM

XRpad_UNDERVOLTAGE_ALARM

XRpad_OVERVOLTAGE_ALARM

XRpad_OVERTEMPERATURE_ALARM

XRpad_BATTERY_UNKNOWN_ERROR

3.10.3.16 enum XRpad_BatteryPresence

Enumerator:

XRpad_NO_BATTERY

XRpad_BATTERY_INSERTED

XRpad_DUMMY_INSERTED

3.10.3.17 enum XRpad_ChargeMode

Enumerator:

XRpad_NOT_CHARGING
XRpad_CHARGING_SLOW
XRpad_CHARGING_NORMAL
XRpad_CHARGING_FAST
XRpad_FULLY_CHARGED
XRpad_DISCHARGING

3.10.3.18 enum XRpad_DataInterfaceControlEnum

Possible image transfer channels for XRpad

Enumerator:

XRpad_DATA_VIA_LAN
XRpad_DATA_VIA_WLAN

3.10.3.19 enum XRpad_SystemControlEnum

Possible system control actions for XRpad

Enumerator:

XRpad_SYSTEM_CONTROL_REBOOT
XRpad_SYSTEM_CONTROL_RESTART_NETWORK
XRpad_SYSTEM_CONTROL_SHUTDOWN
XRpad_SYSTEM_CONTROL_SET_DEEP_SLEEP
XRpad_SYSTEM_CONTROL_SET_IDLE

Chapter 4

Class Documentation

4.1 CHwHeaderInfo Struct Reference

Public Attributes

- int [bAddRow](#)
- BOOL [bAddRow](#)
- int [bPwrSave](#)
- BOOL [bPwrSave](#)
- int [bSyncMode](#)
- BOOL [bSyncMode](#)
- unsigned long [dwAccess](#)
- DWORD [dwAccess](#)
- unsigned long [dwAcqMode](#)
- DWORD [dwAcqMode](#)
- unsigned long [dwBias](#)
- DWORD [dwBias](#)
- unsigned long [dwDataSorting](#)
- DWORD [dwDataSorting](#)
- unsigned long [dwDataType](#)
- DWORD [dwDataType](#)
- unsigned long [dwFrmFillRowIntervalls](#)
- DWORD [dwFrmFillRowIntervalls](#)
- unsigned long [dwFrmNrRows](#)
- DWORD [dwFrmNrRows](#)
- unsigned long [dwFrmRowType](#)
- DWORD [dwFrmRowType](#)
- unsigned long [dwGain](#)
- DWORD [dwGain](#)
- unsigned long [dwHeaderID](#)
- DWORD [dwHeaderID](#)
- unsigned long [dwLeakRows](#)

- DWORD [dwLeakRows](#)
- unsigned long [dwNrColumns](#)
- DWORD [dwNrColumns](#)
- unsigned long [dwNrOfFillingRows](#)
- DWORD [dwNrOfFillingRows](#)
- unsigned long [dwNrRows](#)
- DWORD [dwNrRows](#)
- unsigned long [dwOffset](#)
- DWORD [dwOffset](#)
- unsigned long [dwPROMID](#)
- DWORD [dwPROMID](#)
- unsigned long [dwTiming](#)
- DWORD [dwTiming](#)
- unsigned long [dwZoomBRColumn](#)
- DWORD [dwZoomBRColumn](#)
- unsigned long [dwZoomBRRow](#)
- DWORD [dwZoomBRRow](#)
- unsigned long [dwZoomULColumn](#)
- DWORD [dwZoomULColumn](#)
- unsigned long [dwZoomULRow](#)
- DWORD [dwZoomULRow](#)

4.1.1 Member Data Documentation

4.1.1.1 int CHwHeaderInfo::bAddRow

4.1.1.2 BOOL CHwHeaderInfo::bAddRow

4.1.1.3 int CHwHeaderInfo::bPwrSave

4.1.1.4 BOOL CHwHeaderInfo::bPwrSave

4.1.1.5 int CHwHeaderInfo::bSyncMode

4.1.1.6 BOOL CHwHeaderInfo::bSyncMode

4.1.1.7 unsigned long CHwHeaderInfo::dwAccess

4.1.1.8 DWORD CHwHeaderInfo::dwAccess

4.1.1.9 unsigned long CHwHeaderInfo::dwAcqMode

4.1.1.10 DWORD CHwHeaderInfo::dwAcqMode

4.1.1.11 unsigned long CHwHeaderInfo::dwBias

- 4.1.1.12 **DWORD CHwHeaderInfo::dwBias**
- 4.1.1.13 **unsigned long CHwHeaderInfo::dwDataSorting**
- 4.1.1.14 **DWORD CHwHeaderInfo::dwDataSorting**
- 4.1.1.15 **unsigned long CHwHeaderInfo::dwDataType**
- 4.1.1.16 **DWORD CHwHeaderInfo::dwDataType**
- 4.1.1.17 **unsigned long CHwHeaderInfo::dwFrmFillRowIntervalls**
- 4.1.1.18 **DWORD CHwHeaderInfo::dwFrmFillRowIntervalls**
- 4.1.1.19 **unsigned long CHwHeaderInfo::dwFrmNrRows**
- 4.1.1.20 **DWORD CHwHeaderInfo::dwFrmNrRows**
- 4.1.1.21 **unsigned long CHwHeaderInfo::dwFrmRowType**
- 4.1.1.22 **DWORD CHwHeaderInfo::dwFrmRowType**
- 4.1.1.23 **unsigned long CHwHeaderInfo::dwGain**
- 4.1.1.24 **DWORD CHwHeaderInfo::dwGain**
- 4.1.1.25 **unsigned long CHwHeaderInfo::dwHeaderID**
- 4.1.1.26 **DWORD CHwHeaderInfo::dwHeaderID**
- 4.1.1.27 **unsigned long CHwHeaderInfo::dwLeakRows**
- 4.1.1.28 **DWORD CHwHeaderInfo::dwLeakRows**
- 4.1.1.29 **unsigned long CHwHeaderInfo::dwNrColumns**
- 4.1.1.30 **DWORD CHwHeaderInfo::dwNrColumns**
- 4.1.1.31 **unsigned long CHwHeaderInfo::dwNrOfFillingRows**
- 4.1.1.32 **DWORD CHwHeaderInfo::dwNrOfFillingRows**
- 4.1.1.33 **unsigned long CHwHeaderInfo::dwNrRows**
- 4.1.1.34 **DWORD CHwHeaderInfo::dwNrRows**
- 4.1.1.35 **unsigned long CHwHeaderInfo::dwOffset**

- 4.1.1.36 `DWORD CHwHeaderInfo::dwOffset`
- 4.1.1.37 `unsigned long CHwHeaderInfo::dwPROMID`
- 4.1.1.38 `DWORD CHwHeaderInfo::dwPROMID`
- 4.1.1.39 `unsigned long CHwHeaderInfo::dwTiming`
- 4.1.1.40 `DWORD CHwHeaderInfo::dwTiming`
- 4.1.1.41 `unsigned long CHwHeaderInfo::dwZoomBRColumn`
- 4.1.1.42 `DWORD CHwHeaderInfo::dwZoomBRColumn`
- 4.1.1.43 `unsigned long CHwHeaderInfo::dwZoomBRRow`
- 4.1.1.44 `DWORD CHwHeaderInfo::dwZoomBRRow`
- 4.1.1.45 `unsigned long CHwHeaderInfo::dwZoomULColumn`
- 4.1.1.46 `DWORD CHwHeaderInfo::dwZoomULColumn`
- 4.1.1.47 `unsigned long CHwHeaderInfo::dwZoomULRow`
- 4.1.1.48 `DWORD CHwHeaderInfo::dwZoomULRow`

4.2 CHwHeaderInfoEx Struct Reference

Public Attributes

- unsigned short [wAccess](#)
- WORD [wAccess](#)
- unsigned short [wBias](#)
- WORD [wBias](#)
- unsigned short [wBinningMode](#)
- WORD [wBinningMode](#)
- unsigned short [wCameratype](#)
- WORD [wCameratype](#)
- unsigned short [wClock](#)
- WORD [wClock](#)
- unsigned short [wCommand1](#)
- WORD [wCommand1](#)
- unsigned short [wCommand2](#)
- WORD [wCommand2](#)
- unsigned short [wCommand3](#)
- WORD [wCommand3](#)
- unsigned short [wCommand4](#)

- WORD [wCommand4](#)
- unsigned short [wDataSorting](#)
- WORD [wDataSorting](#)
- unsigned short [wDummy](#)
- WORD [wDummy](#)
- unsigned short [wFrameCnt](#)
- WORD [wFrameCnt](#)
- unsigned short [wFrmNrRows](#)
- WORD [wFrmNrRows](#)
- unsigned short [wFrmRowType](#)
- WORD [wFrmRowType](#)
- unsigned short [wGain](#)
- WORD [wGain](#)
- unsigned short [wHeaderID](#)
- WORD [wHeaderID](#)
- unsigned short [wLeakRows](#)
- WORD [wLeakRows](#)
- unsigned short [wNrColumns](#)
- WORD [wNrColumns](#)
- unsigned short [wNrRows](#)
- WORD [wNrRows](#)
- unsigned short [wPROMID](#)
- WORD [wPROMID](#)
- unsigned short [wRealInttime_microSec](#)
- WORD [wRealInttime_microSec](#)
- unsigned short [wRealInttime_milliSec](#)
- WORD [wRealInttime_milliSec](#)
- unsigned short [wResolutionX](#)
- WORD [wResolutionX](#)
- unsigned short [wResolutionY](#)
- WORD [wResolutionY](#)
- unsigned short [wRowTime](#)
- WORD [wRowTime](#)
- unsigned short [wStatus](#)
- WORD [wStatus](#)
- unsigned short [wTiming](#)
- WORD [wTiming](#)
- unsigned short [wUgComp](#)
- WORD [wUgComp](#)
- unsigned short [wZoomBRColumn](#)
- WORD [wZoomBRColumn](#)
- unsigned short [wZoomBRRow](#)
- WORD [wZoomBRRow](#)
- unsigned short [wZoomULColumn](#)
- WORD [wZoomULColumn](#)
- unsigned short [wZoomULRow](#)
- WORD [wZoomULRow](#)

4.2.1 Member Data Documentation

4.2.1.1 unsigned short CHwHeaderInfoEx::wAccess

4.2.1.2 WORD CHwHeaderInfoEx::wAccess

4.2.1.3 unsigned short CHwHeaderInfoEx::wBias

4.2.1.4 WORD CHwHeaderInfoEx::wBias

4.2.1.5 unsigned short CHwHeaderInfoEx::wBinningMode

4.2.1.6 WORD CHwHeaderInfoEx::wBinningMode

4.2.1.7 unsigned short CHwHeaderInfoEx::wCameratype

4.2.1.8 WORD CHwHeaderInfoEx::wCameratype

4.2.1.9 unsigned short CHwHeaderInfoEx::wClock

4.2.1.10 WORD CHwHeaderInfoEx::wClock

4.2.1.11 unsigned short CHwHeaderInfoEx::wCommand1

4.2.1.12 WORD CHwHeaderInfoEx::wCommand1

4.2.1.13 unsigned short CHwHeaderInfoEx::wCommand2

4.2.1.14 WORD CHwHeaderInfoEx::wCommand2

4.2.1.15 unsigned short CHwHeaderInfoEx::wCommand3

4.2.1.16 WORD CHwHeaderInfoEx::wCommand3

4.2.1.17 unsigned short CHwHeaderInfoEx::wCommand4

4.2.1.18 WORD CHwHeaderInfoEx::wCommand4

4.2.1.19 unsigned short CHwHeaderInfoEx::wDataSorting

4.2.1.20 WORD CHwHeaderInfoEx::wDataSorting

4.2.1.21 unsigned short CHwHeaderInfoEx::wDummy

4.2.1.22 WORD CHwHeaderInfoEx::wDummy

4.2.1.23 unsigned short CHwHeaderInfoEx::wFrameCnt

- 4.2.1.24 WORD CHwHeaderInfoEx::wFrameCnt
- 4.2.1.25 unsigned short CHwHeaderInfoEx::wFrmNrRows
- 4.2.1.26 WORD CHwHeaderInfoEx::wFrmNrRows
- 4.2.1.27 unsigned short CHwHeaderInfoEx::wFrmRowType
- 4.2.1.28 WORD CHwHeaderInfoEx::wFrmRowType
- 4.2.1.29 unsigned short CHwHeaderInfoEx::wGain
- 4.2.1.30 WORD CHwHeaderInfoEx::wGain
- 4.2.1.31 unsigned short CHwHeaderInfoEx::wHeaderID
- 4.2.1.32 WORD CHwHeaderInfoEx::wHeaderID
- 4.2.1.33 unsigned short CHwHeaderInfoEx::wLeakRows
- 4.2.1.34 WORD CHwHeaderInfoEx::wLeakRows
- 4.2.1.35 unsigned short CHwHeaderInfoEx::wNrColumns
- 4.2.1.36 WORD CHwHeaderInfoEx::wNrColumns
- 4.2.1.37 unsigned short CHwHeaderInfoEx::wNrRows
- 4.2.1.38 WORD CHwHeaderInfoEx::wNrRows
- 4.2.1.39 unsigned short CHwHeaderInfoEx::wPROMID
- 4.2.1.40 WORD CHwHeaderInfoEx::wPROMID
- 4.2.1.41 unsigned short CHwHeaderInfoEx::wRealInttime_microSec
- 4.2.1.42 WORD CHwHeaderInfoEx::wRealInttime_microSec
- 4.2.1.43 unsigned short CHwHeaderInfoEx::wRealInttime_milliSec
- 4.2.1.44 WORD CHwHeaderInfoEx::wRealInttime_milliSec
- 4.2.1.45 unsigned short CHwHeaderInfoEx::wResolutionX
- 4.2.1.46 WORD CHwHeaderInfoEx::wResolutionX
- 4.2.1.47 unsigned short CHwHeaderInfoEx::wResolutionY

- 4.2.1.48 WORD CHwHeaderInfoEx::wResolutionY
- 4.2.1.49 unsigned short CHwHeaderInfoEx::wRowTime
- 4.2.1.50 WORD CHwHeaderInfoEx::wRowTime
- 4.2.1.51 unsigned short CHwHeaderInfoEx::wStatus
- 4.2.1.52 WORD CHwHeaderInfoEx::wStatus
- 4.2.1.53 unsigned short CHwHeaderInfoEx::wTiming
- 4.2.1.54 WORD CHwHeaderInfoEx::wTiming
- 4.2.1.55 unsigned short CHwHeaderInfoEx::wUgComp
- 4.2.1.56 WORD CHwHeaderInfoEx::wUgComp
- 4.2.1.57 unsigned short CHwHeaderInfoEx::wZoomBRColumn
- 4.2.1.58 WORD CHwHeaderInfoEx::wZoomBRColumn
- 4.2.1.59 unsigned short CHwHeaderInfoEx::wZoomBRRow
- 4.2.1.60 WORD CHwHeaderInfoEx::wZoomBRRow
- 4.2.1.61 unsigned short CHwHeaderInfoEx::wZoomULColumn
- 4.2.1.62 WORD CHwHeaderInfoEx::wZoomULColumn
- 4.2.1.63 unsigned short CHwHeaderInfoEx::wZoomULRow
- 4.2.1.64 WORD CHwHeaderInfoEx::wZoomULRow

4.3 DETECTOR_BATTERY Struct Reference

Public Attributes

- DWORD [capacity](#)
- DWORD [charge](#)
- DWORD [current](#)
- DWORD [cycle_count](#)
- DWORD [energy](#)
- DWORD [serial_no](#)
- DWORD [status](#)
- DWORD [temperature](#)
- DWORD [voltage](#)

4.3.1 Member Data Documentation

4.3.1.1 `DWORD DETECTOR_BATTERY::capacity`

4.3.1.2 `DWORD DETECTOR_BATTERY::charge`

4.3.1.3 `DWORD DETECTOR_BATTERY::current`

4.3.1.4 `DWORD DETECTOR_BATTERY::cycle_count`

4.3.1.5 `DWORD DETECTOR_BATTERY::energy`

4.3.1.6 `DWORD DETECTOR_BATTERY::serial_no`

4.3.1.7 `DWORD DETECTOR_BATTERY::status`

4.3.1.8 `DWORD DETECTOR_BATTERY::temperature`

4.3.1.9 `DWORD DETECTOR_BATTERY::voltage`

4.4 DETECTOR_CURRENT_VOLTAGE Struct Reference

Public Attributes

- `int imA1`
- `int imA2`
- `int imA3`
- `int iV1`
- `int iV2`
- `int iV3`

4.4.1 Member Data Documentation

4.4.1.1 `int DETECTOR_CURRENT_VOLTAGE::imA1`

4.4.1.2 `int DETECTOR_CURRENT_VOLTAGE::imA2`

4.4.1.3 `int DETECTOR_CURRENT_VOLTAGE::imA3`

4.4.1.4 `int DETECTOR_CURRENT_VOLTAGE::iV1`

4.4.1.5 `int DETECTOR_CURRENT_VOLTAGE::iV2`

4.4.1.6 `int DETECTOR_CURRENT_VOLTAGE::iV3`

4.5 deviceInfo Struct Reference

Public Attributes

- char [device_version](#) [16]
The device version.
- char [manufacturer_name](#) [32]
The manufactures name.
- char [manufacturer_specific](#) [48]
Some manufacture info.
- char [model_name](#) [32]
The model name.
- char [serial_number](#) [16]
The serial number.
- char [spec_version](#) [16]
GigE vision spec version.
- char [user_name](#) [16]
Device specific.

4.5.1 Detailed Description

The device info

4.5.2 Member Data Documentation

4.5.2.1 char deviceInfo::device_version[16]

The device version.

4.5.2.2 char deviceInfo::manufacturer_name[32]

The manufactures name.

4.5.2.3 char deviceInfo::manufacturer_specific[48]

Some manufacture info.

4.5.2.4 char deviceInfo::model_name[32]

The model name.

4.5.2.5 char `deviceInfo::serial_number`[16]

The serial number.

4.5.2.6 char `deviceInfo::spec_version`[16]

GigE vision spec version.

4.5.2.7 char `deviceInfo::user_name`[16]

Device specific.

4.6 `discoveryReply` Struct Reference

Public Attributes

- struct [deviceInfo](#) `deviceInfo`
The device information.
- char [gvcp_ip](#) [16]
Which IP address is used for image transfer.
- struct [networkInfo](#) `lanInfo`
The LAN network setup.
- struct [networkInfo](#) `wlanInfo`
The WLAN network setup.

4.6.1 Detailed Description

The original discovery reply

4.6.2 Member Data Documentation

4.6.2.1 struct `deviceInfo` `discoveryReply::deviceInfo`

The device information.

4.6.2.2 char `discoveryReply::gvcp_ip`[16]

Which IP address is used for image transfer.

4.6.2.3 struct `networkInfo` `discoveryReply::lanInfo`

The LAN network setup.

4.6.2.4 struct `networkInfo` `discoveryReply::wlanInfo`

The WLAN network setup.

4.7 `discoveryReplyEx` Struct Reference

Public Attributes

- struct `deviceInfo` `deviceInfo`
The device information.
- char `gvcp_ip` [16]
Which IP address is used for image transfer.
- struct `networkInfo` `lanInfo`
The LAN network setup.
- unsigned `messageCount`
How many messages this reply carries.
- struct `discoveryReplyMsg` `messages` [32]
Info about the received reply packets.
- struct `networkInfo` `wlanInfo`
The WLAN network setup.

4.7.1 Detailed Description

The extended discovery reply Can carry additional information

4.7.2 Member Data Documentation

4.7.2.1 struct `deviceInfo` `discoveryReplyEx::deviceInfo`

The device information.

4.7.2.2 char `discoveryReplyEx::gvcp_ip`[16]

Which IP address is used for image transfer.

4.7.2.3 struct `networkInfo` `discoveryReplyEx::lanInfo`

The LAN network setup.

4.7.2.4 unsigned `discoveryReplyEx::messageCount`

How many messages this reply carries.

4.7.2.5 struct `discoveryReplyMsg` `discoveryReplyEx::messages`[32]

Info about the received reply packets.

4.7.2.6 struct `networkInfo` `discoveryReplyEx::wlanInfo`

The WLAN network setup.

4.8 `discoveryReplyMsg` Struct Reference

Public Attributes

- char `local_ip` [16]
- char `remote_ip` [16]

4.8.1 Detailed Description

Information about a single discovery Reply network packet

4.8.2 Member Data Documentation

4.8.2.1 char `discoveryReplyMsg::local_ip`[16]

4.8.2.2 char `discoveryReplyMsg::remote_ip`[16]

4.9 `EPC_REGISTER` Struct Reference

Public Attributes

- DWORD `active_network_config`
- `DETECTOR_BATTERY` battery
- DWORD `channel`
- DWORD `exam_flag`
- DWORD `lan_status_register`
- DWORD `power_state`
- `RTC_STRUCT` `rtc_value`
- DWORD `sdcard_state`
- DWORD `sdcard_usage`
- DWORD `signal_strength`
- DWORD `spartan_id`
- `CHwHeaderInfoEx` `spartan_register`
- DWORD `temperature_error_level` [8]
- DWORD `temperature_value` [8]

- DWORD [temperature_warning_level](#) [8]
- DWORD [version](#)
- DWORD [wlan_status_register](#)

4.9.1 Member Data Documentation

4.9.1.1 DWORD EPC_REGISTER::active_network_config

4.9.1.2 DETECTOR_BATTERY EPC_REGISTER::battery

4.9.1.3 DWORD EPC_REGISTER::channel

4.9.1.4 DWORD EPC_REGISTER::exam_flag

4.9.1.5 DWORD EPC_REGISTER::lan_status_register

4.9.1.6 DWORD EPC_REGISTER::power_state

4.9.1.7 RTC_STRUCT EPC_REGISTER::rtc_value

4.9.1.8 DWORD EPC_REGISTER::sdcard_state

4.9.1.9 DWORD EPC_REGISTER::sdcard_usage

4.9.1.10 DWORD EPC_REGISTER::signal_strength

4.9.1.11 DWORD EPC_REGISTER::spartan_id

4.9.1.12 CHwHeaderInfoEx EPC_REGISTER::spartan_register

4.9.1.13 DWORD EPC_REGISTER::temperature_error_level[8]

4.9.1.14 DWORD EPC_REGISTER::temperature_value[8]

4.9.1.15 DWORD EPC_REGISTER::temperature_warning_level[8]

4.9.1.16 DWORD EPC_REGISTER::version

4.9.1.17 DWORD EPC_REGISTER::wlan_status_register

4.10 FPGAType Struct Reference

Public Attributes

- unsigned char [wTiming](#)
- unsigned char [wValue0](#)

- unsigned char [wValue1](#)
- unsigned char [wValue2](#)
- unsigned char [wValue3](#)
- unsigned char [wValue4](#)
- unsigned char [wValue5](#)
- unsigned char [wValue6](#)

4.10.1 Member Data Documentation

4.10.1.1 unsigned char [FPGAType::wTiming](#)

4.10.1.2 unsigned char [FPGAType::wValue0](#)

4.10.1.3 unsigned char [FPGAType::wValue1](#)

4.10.1.4 unsigned char [FPGAType::wValue2](#)

4.10.1.5 unsigned char [FPGAType::wValue3](#)

4.10.1.6 unsigned char [FPGAType::wValue4](#)

4.10.1.7 unsigned char [FPGAType::wValue5](#)

4.10.1.8 unsigned char [FPGAType::wValue6](#)

4.11 GBIF_Detector_Properties Struct Reference

Public Attributes

- char [cDetectorType](#) [32]
- char [cDeviceIdentifier](#) [16]
- char [cDummy](#) [48]
- char [cManufacturingDate](#) [8]
- char [cPlaceOfManufacture](#) [8]
- char [cUniqueDeviceIdentifier](#) [16]

4.11.1 Member Data Documentation

4.11.1.1 char [GBIF_Detector_Properties::cDetectorType](#)

4.11.1.2 char [GBIF_Detector_Properties::cDeviceIdentifier](#)[16]

4.11.1.3 char [GBIF_Detector_Properties::cDummy](#)

4.11.1.4 char [GBIF_Detector_Properties::cManufacturingDate](#)

4.11.1.5 char GBIF_Detector_Properties::cPlaceOfManufacture

4.11.1.6 char GBIF_Detector_Properties::cUniqueDeviceIdentifier[16]

4.12 GBIF_DEVICE_PARAM Struct Reference

Public Attributes

- char [cDeviceName](#) [16]
- CHAR [cDeviceName](#) [16]
- char [cGBIFFirmwareVersion](#) [32]
- CHAR [cGBIFFirmwareVersion](#) [32]
- char [cManufacturerName](#) [32]
- CHAR [cManufacturerName](#) [32]
- char [cModelName](#) [32]
- CHAR [cModelName](#) [32]
- unsigned long [dwIPCurrentBootOptions](#)
- DWORD [dwIPCurrentBootOptions](#)
- GBIF_STRING_DATATYPE [ucAdapterIP](#) [GBIF_IP_MAC_NAME_CHAR_ARRAY_LENGTH]
- GBIF_STRING_DATATYPE [ucAdapterMask](#) [GBIF_IP_MAC_NAME_CHAR_ARRAY_LENGTH]
- GBIF_STRING_DATATYPE [ucGateway](#) [GBIF_IP_MAC_NAME_CHAR_ARRAY_LENGTH]
- GBIF_STRING_DATATYPE [ucIP](#) [GBIF_IP_MAC_NAME_CHAR_ARRAY_LENGTH]
- GBIF_STRING_DATATYPE [ucMacAddress](#) [GBIF_IP_MAC_NAME_CHAR_ARRAY_LENGTH]
- GBIF_STRING_DATATYPE [ucSubnetMask](#) [GBIF_IP_MAC_NAME_CHAR_ARRAY_LENGTH]

4.12.1 Member Data Documentation

4.12.1.1 char GBIF_DEVICE_PARAM::cDeviceName[16]

4.12.1.2 CHAR GBIF_DEVICE_PARAM::cDeviceName[16]

4.12.1.3 char GBIF_DEVICE_PARAM::cGBIFFirmwareVersion[32]

4.12.1.4 CHAR GBIF_DEVICE_PARAM::cGBIFFirmwareVersion[32]

4.12.1.5 char GBIF_DEVICE_PARAM::cManufacturerName[32]

4.12.1.6 CHAR GBIF_DEVICE_PARAM::cManufacturerName[32]

4.12.1.7 char GBIF_DEVICE_PARAM::cModelName[32]

- 4.12.1.8 CHAR GBIF_DEVICE_PARAM::cModelName[32]
- 4.12.1.9 unsigned long GBIF_DEVICE_PARAM::dwIPCurrentBootOptions
- 4.12.1.10 DWORD GBIF_DEVICE_PARAM::dwIPCurrentBootOptions
- 4.12.1.11 GBIF_STRING_DATATYPE GBIF_DEVICE_PARAM::ucAdapterIP
- 4.12.1.12 GBIF_STRING_DATATYPE GBIF_DEVICE_PARAM::ucAdapterMask
- 4.12.1.13 GBIF_STRING_DATATYPE GBIF_DEVICE_PARAM::ucGateway
- 4.12.1.14 GBIF_STRING_DATATYPE GBIF_DEVICE_PARAM::ucIP
- 4.12.1.15 GBIF_STRING_DATATYPE GBIF_DEVICE_PARAM::ucMacAddress
- 4.12.1.16 GBIF_STRING_DATATYPE GBIF_DEVICE_PARAM::ucSubnetMask

4.13 networkAdapterConfiguration Struct Reference

Public Attributes

- int [bridged](#)
Is the device in a bridge, this is not changeable by API.
- char [dns](#) [16]
DNS server.
- int [enabled](#)
Enabled flag, this is not changeable for LAN by API, but changeable for WLAN.
- char [gateway](#) [16]
Gateway.
- int [hw_accel](#)
Used image transfer, this is not changeable by API.
- char [ifname](#) [16]
Interface name (eth0), this is not changeable by API.
- char [ipaddr](#) [16]
IP address.
- char [macaddr](#) [18]
MAC address, this is not changeable by API.
- char [netmask](#) [16]
Netmask.
- char [not_used](#) [110]
To fill up the struct to 320 byte.
- char [proto](#) [16]
"static" or "dhcp", only these two options are available

4.13.1 Detailed Description

Structure for holding the adapter part of a configuration

4.13.2 Member Data Documentation

4.13.2.1 `int networkAdapterConfiguration::bridged`

Is the device in a bridge, this is not changeable by API.

4.13.2.2 `char networkAdapterConfiguration::dns[16]`

DNS server.

4.13.2.3 `int networkAdapterConfiguration::enabled`

Enabled flag, this is not changeable for LAN by API, but changeable for WLAN.

4.13.2.4 `char networkAdapterConfiguration::gateway[16]`

Gateway.

4.13.2.5 `int networkAdapterConfiguration::hw_accel`

Used image transfer, this is not changeable by API.

4.13.2.6 `char networkAdapterConfiguration::ifname[16]`

Interface name (eth0), this is not changeable by API.

4.13.2.7 `char networkAdapterConfiguration::ipaddr[16]`

IP address.

4.13.2.8 `char networkAdapterConfiguration::macaddr[18]`

MAC address, this is not changeable by API.

4.13.2.9 `char networkAdapterConfiguration::netmask[16]`

Netmask.

4.13.2.10 char networkAdapterConfiguration::not_used[110]

To fill up the struct to 320 byte.

4.13.2.11 char networkAdapterConfiguration::proto[16]

"static" or "dhcp", only these two options are available

4.14 networkConfiguration Struct Reference

Public Attributes

- int [gbif_enabled](#)
Is the GBif enabled.
- char [hostname](#) [80]
The hostname.
- struct [networkAdapterConfiguration lan](#)
LAN.
- char [name](#) [80]
The configuration name.
- char [notUsed](#) [184]
For later extensions.
- char [path](#) [128]
The configurations path, this is not changeable by API.
- int [readonly](#)
Is the configuration readonly, this is not changeable by API.
- int [sshd_enabled](#)
SSH daemon enabled.
- struct [wifiConfigurationEx wifi](#)
Wifi configuration, extended version.
- struct [networkAdapterConfiguration wlan](#)
WLAN.

4.14.1 Detailed Description

Structure for holding the complete network configuration.

4.14.2 Member Data Documentation

4.14.2.1 int networkConfiguration::gbif_enabled

Is the GBif enabled.

4.14.2.2 char networkConfiguration::hostname[80]

The hostname.

4.14.2.3 struct networkAdapterConfiguration networkConfiguration::lan

LAN.

4.14.2.4 char networkConfiguration::name[80]

The configuration name.

4.14.2.5 char networkConfiguration::notUsed[184]

For later extensions.

4.14.2.6 char networkConfiguration::path[128]

The configurations path, this is not changeable by API.

4.14.2.7 int networkConfiguration::readonly

Is the configuration readonly, this is not changeable by API.

4.14.2.8 int networkConfiguration::sshd_enabled

SSH daemon enabled.

4.14.2.9 struct wifiConfigurationEx networkConfiguration::wifi

Wifi configuration, extended version.

4.14.2.10 struct networkAdapterConfiguration networkConfiguration::wlan

WLAN.

4.15 networkInfo Struct Reference

Public Attributes

- char [broadcast](#) [16]

The IP broadcast as string.

- char [ip](#) [16]

The IP (v4) address as string.

- char [mac](#) [18]

The MAC address as string.

- char [mask](#) [16]

The IP mask as string.

4.15.1 Detailed Description

The network information

4.15.2 Member Data Documentation

4.15.2.1 char `networkInfo::broadcast`[16]

The IP broadcast as string.

4.15.2.2 char `networkInfo::ip`[16]

The IP (v4) address as string.

4.15.2.3 char `networkInfo::mac`[18]

The MAC address as string.

4.15.2.4 char `networkInfo::mask`[16]

The IP mask as string.

4.16 RTC_STRUCT Struct Reference

Public Attributes

- DWORD [day](#)
- DWORD [hour](#)
- DWORD [minute](#)
- DWORD [month](#)
- DWORD [second](#)
- DWORD [year](#)

4.16.1 Member Data Documentation

4.16.1.1 DWORD RTC_STRUCT::day

4.16.1.2 DWORD RTC_STRUCT::hour

4.16.1.3 DWORD RTC_STRUCT::minute

4.16.1.4 DWORD RTC_STRUCT::month

4.16.1.5 DWORD RTC_STRUCT::second

4.16.1.6 DWORD RTC_STRUCT::year

4.17 wifiConfiguration Struct Reference

Public Attributes

- char [agmode](#) [32]
agmode
- int [channel](#)
Channel.
- char [description](#) [64]
Contains the description in case of a station.
- char [mode](#) [32]
Accesspoint or client.
- char [ssid](#) [64]
Own SSID if mode == "ap" or the accesspoints ssid.

4.17.1 Detailed Description

Wifi configurations

4.17.2 Member Data Documentation

4.17.2.1 char wifiConfiguration::agmode[32]

agmode

4.17.2.2 int wifiConfiguration::channel

Channel.

4.17.2.3 char `wifiConfiguration::description`[64]

Contains the description in case of a station.

4.17.2.4 char `wifiConfiguration::mode`[32]

Accesspoint or client.

4.17.2.5 char `wifiConfiguration::ssid`[64]

Own SSID if mode == "ap" or the accesspoints ssid.

4.18 `wifiConfigurationEx` Struct Reference

Public Attributes

- char `agmode` [32]
agmode
- int `channel`
Channel, only valid options will be accepted, otherwise will return with error.
- char `description` [64]
Contains the description in case of a station.
- char `mode` [32]
Accesspoint or client.
- char `passphrase` [68]
new Passphrase (station or accesspoint) (may be 64 byte)
- int `scan_ssid`
only of station mode: scan the SSID
- char `ssid` [64]
Own SSID if mode == "ap" or the accesspoints ssid.

4.18.1 Detailed Description

Wifi configurations, extended version

4.18.2 Member Data Documentation

4.18.2.1 char `wifiConfigurationEx::agmode`[32]

`agmode`

4.18.2.2 `int wifiConfigurationEx::channel`

Channel, only valid options will be accepted, otherwise will return with error.

4.18.2.3 `char wifiConfigurationEx::description[64]`

Contains the description in case of a station.

4.18.2.4 `char wifiConfigurationEx::mode[32]`

Accesspoint or client.

4.18.2.5 `char wifiConfigurationEx::passphrase[68]`

new Passphrase (station or accesspoint) (may be 64 byte)

4.18.2.6 `int wifiConfigurationEx::scan_ssid`

only of station mode: scan the SSID

4.18.2.7 `char wifiConfigurationEx::ssid[64]`

Own SSID if mode == "ap" or the accesspoints ssid.

4.19 WinHeaderType Struct Reference

Public Attributes

- unsigned short [BRX](#)
- WORD [BRX](#)
- unsigned short [BRY](#)
- WORD [BRY](#)
- unsigned short [Correction](#)
- WORD [Correction](#)
- unsigned long [FileSize](#)
- UINT [FileSize](#)
- unsigned short [FileType](#)
- WORD [FileType](#)
- unsigned short [HeaderSize](#)
- WORD [HeaderSize](#)
- unsigned short [HeaderVersion](#)
- WORD [HeaderVersion](#)
- unsigned short [ImageHeaderSize](#)

- WORD [ImageHeaderSize](#)
- double [IntegrationTime](#)
- unsigned short [NrOfFrames](#)
- WORD [NrOfFrames](#)
- unsigned short [TypeOfNumbers](#)
- WORD [TypeOfNumbers](#)
- unsigned short [ULX](#)
- WORD [ULX](#)
- unsigned short [ULY](#)
- WORD [ULY](#)
- unsigned char [x](#) [WINRESTSIZE]
- BYTE [x](#) [WINRESTSIZE]

4.19.1 Member Data Documentation

4.19.1.1 unsigned short WinHeaderType::BRX

4.19.1.2 WORD WinHeaderType::BRX

4.19.1.3 unsigned short WinHeaderType::BRY

4.19.1.4 WORD WinHeaderType::BRY

4.19.1.5 unsigned short WinHeaderType::Correction

4.19.1.6 WORD WinHeaderType::Correction

4.19.1.7 unsigned long WinHeaderType::FileSize

4.19.1.8 UINT WinHeaderType::FileSize

4.19.1.9 unsigned short WinHeaderType::FileType

4.19.1.10 WORD WinHeaderType::FileType

4.19.1.11 unsigned short WinHeaderType::HeaderSize

4.19.1.12 WORD WinHeaderType::HeaderSize

4.19.1.13 unsigned short WinHeaderType::HeaderVersion

4.19.1.14 WORD WinHeaderType::HeaderVersion

4.19.1.15 unsigned short WinHeaderType::ImageHeaderSize

4.19.1.16 WORD WinHeaderType::ImageHeaderSize

- 4.19.1.17 double WinHeaderType::IntegrationTime
- 4.19.1.18 unsigned short WinHeaderType::NrOfFrames
- 4.19.1.19 WORD WinHeaderType::NrOfFrames
- 4.19.1.20 unsigned short WinHeaderType::TypeOfNumbers
- 4.19.1.21 WORD WinHeaderType::TypeOfNumbers
- 4.19.1.22 unsigned short WinHeaderType::ULX
- 4.19.1.23 WORD WinHeaderType::ULX
- 4.19.1.24 unsigned short WinHeaderType::ULY
- 4.19.1.25 WORD WinHeaderType::ULY
- 4.19.1.26 unsigned char WinHeaderType::x[WINRESTSIZE]
- 4.19.1.27 BYTE WinHeaderType::x[WINRESTSIZE]

4.20 WinHeaderType101 Struct Reference

Public Attributes

- WORD [BRX](#)
- WORD [BRY](#)
- WORD [Correction](#)
- UINT [FileSize](#)
- WORD [FileType](#)
- WORD [HeaderSize](#)
- WORD [HeaderVersion](#)
- WORD [ImageHeaderSize](#)
- double [IntegrationTime](#)
- WORD [NrOfFrames](#)
- WORD [TypeOfNumbers](#)
- WORD [ULX](#)
- WORD [ULY](#)
- WORD [wMedianValue](#)
- BYTE [x](#) [WINRESTSIZE101]

4.20.1 Member Data Documentation

4.20.1.1 WORD WinHeaderType101::BRX

- 4.20.1.2 WORD WinHeaderType101::BRY
- 4.20.1.3 WORD WinHeaderType101::Correction
- 4.20.1.4 UINT WinHeaderType101::FileSize
- 4.20.1.5 WORD WinHeaderType101::FileType
- 4.20.1.6 WORD WinHeaderType101::HeaderSize
- 4.20.1.7 WORD WinHeaderType101::HeaderVersion
- 4.20.1.8 WORD WinHeaderType101::ImageHeaderSize
- 4.20.1.9 double WinHeaderType101::IntegrationTime
- 4.20.1.10 WORD WinHeaderType101::NrOfFrames
- 4.20.1.11 WORD WinHeaderType101::TypeOfNumbers
- 4.20.1.12 WORD WinHeaderType101::ULX
- 4.20.1.13 WORD WinHeaderType101::ULY
- 4.20.1.14 WORD WinHeaderType101::wMedianValue
- 4.20.1.15 BYTE WinHeaderType101::x[WINRESTSIZE101]

4.21 WinImageHeaderType Struct Reference

Public Attributes

- unsigned long [dwPROMID](#)
- DWORD [dwPROMID](#)
- float [fAmpere](#)
- float [fKVolt](#)
- unsigned short [n_avframes](#)
- WORD [n_avframes](#)
- char [strPrefilter](#) [9]
- char [strProject](#) [6]
- char [strSystemused](#) [3]

4.21.1 Member Data Documentation

- 4.21.1.1 unsigned long WinImageHeaderType::dwPROMID

4.21.1.2 `DWORD WinImageHeaderType::dwPROMID`

4.21.1.3 `float WinImageHeaderType::fAmpere`

4.21.1.4 `float WinImageHeaderType::fKVolt`

4.21.1.5 `unsigned short WinImageHeaderType::n_avframes`

4.21.1.6 `WORD WinImageHeaderType::n_avframes`

4.21.1.7 `char WinImageHeaderType::strPrefilter`

4.21.1.8 `char WinImageHeaderType::strProject`

4.21.1.9 `char WinImageHeaderType::strSystemused`

4.22 XislLoggingErrorHandler Class Reference

Public Member Functions

- virtual void [error](#) (const log4cplus::tstring &err)
- virtual void [reset](#) ()
- [XislLoggingErrorHandler](#) ()
- virtual [~XislLoggingErrorHandler](#) ()

4.22.1 Constructor & Destructor Documentation

4.22.1.1 `XislLoggingErrorHandler::XislLoggingErrorHandler ()`

4.22.1.2 `XislLoggingErrorHandler::~~XislLoggingErrorHandler ()` [virtual]

4.22.2 Member Function Documentation

4.22.2.1 `void XislLoggingErrorHandler::error (const log4cplus::tstring & err)`
[virtual]

4.22.2.2 `void XislLoggingErrorHandler::reset ()` [virtual]

4.23 XRpad_BatteryStatus Struct Reference

Public Attributes

- int [authenticated](#)
- [XRpad_ChargeMode](#) [charge_mode](#)
- int [charge_state](#)

- int [cycle_count](#)
- int [design_capacity](#)
- int [health](#)
- [XRpad_BatteryPresence](#) [presence](#)
- int [remaining_capacity](#)
- int [temperature](#)

4.23.1 Member Data Documentation

4.23.1.1 int [XRpad_BatteryStatus::authenticated](#)

4.23.1.2 [XRpad_ChargeMode](#) [XRpad_BatteryStatus::charge_mode](#)

4.23.1.3 int [XRpad_BatteryStatus::charge_state](#)

4.23.1.4 int [XRpad_BatteryStatus::cycle_count](#)

4.23.1.5 int [XRpad_BatteryStatus::design_capacity](#)

4.23.1.6 int [XRpad_BatteryStatus::health](#)

4.23.1.7 [XRpad_BatteryPresence](#) [XRpad_BatteryStatus::presence](#)

4.23.1.8 int [XRpad_BatteryStatus::remaining_capacity](#)

4.23.1.9 int [XRpad_BatteryStatus::temperature](#)

4.24 XRpad_ShockEvent Struct Reference

Public Attributes

- unsigned int [critical_sensor1](#)
- unsigned int [critical_sensor2](#)
- unsigned int [critical_sensor3](#)
- unsigned int [timestamp](#)
- unsigned int [warning_sensor1](#)
- unsigned int [warning_sensor2](#)
- unsigned int [warning_sensor3](#)

4.24.1 Member Data Documentation

4.24.1.1 unsigned int [XRpad_ShockEvent::critical_sensor1](#)

4.24.1.2 unsigned int [XRpad_ShockEvent::critical_sensor2](#)

4.24.1.3 unsigned int `XRpad_ShockEvent::critical_sensor3`

4.24.1.4 unsigned int `XRpad_ShockEvent::timestamp`

4.24.1.5 unsigned int `XRpad_ShockEvent::warning_sensor1`

4.24.1.6 unsigned int `XRpad_ShockEvent::warning_sensor2`

4.24.1.7 unsigned int `XRpad_ShockEvent::warning_sensor3`

4.25 XRpad_ShockSensorReport Struct Reference

Public Attributes

- [XRpad_ShockEvent largest](#)
- [XRpad_ShockEvent latest](#)

4.25.1 Member Data Documentation

4.25.1.1 `XRpad_ShockEvent XRpad_ShockSensorReport::largest`

4.25.1.2 `XRpad_ShockEvent XRpad_ShockSensorReport::latest`

4.26 XRpad_TempSensor Struct Reference

Public Attributes

- char [index](#)
- BOOL [is_virtual](#)
- char [name](#) [32]
- double [temperature](#)
- unsigned char [warn_level](#)

4.26.1 Member Data Documentation

4.26.1.1 `char XRpad_TempSensor::index`

4.26.1.2 `BOOL XRpad_TempSensor::is_virtual`

4.26.1.3 `char XRpad_TempSensor::name[32]`

4.26.1.4 `double XRpad_TempSensor::temperature`

4.26.1.5 `unsigned char XRpad_TempSensor::warn_level`

4.27 XRpad_TempSensorReport Struct Reference

Public Attributes

- unsigned char [sensor_count](#)
- [XRpad_TempSensor](#) [sensors](#) [16]
- unsigned int [shutdown_time](#)
- unsigned char [system_warn_level](#)

4.27.1 Member Data Documentation

4.27.1.1 unsigned char XRpad_TempSensorReport::sensor_count

4.27.1.2 XRpad_TempSensor XRpad_TempSensorReport::sensors[16]

4.27.1.3 unsigned int XRpad_TempSensorReport::shutdown_time

4.27.1.4 unsigned char XRpad_TempSensorReport::system_warn_level

4.28 XRpad_VersionInfo Struct Reference

Public Attributes

- char [hwdriver](#) [32]
- char [linux_kernel](#) [32]
- char [msp_firmware](#) [32]
- char [pld_firmware](#) [32]
- char [software](#) [32]
- char [spartan_firmware](#) [32]
- char [subversion](#) [256]
- char [wlan](#) [32]
- char [xrpd](#) [32]
- char [zynq_firmware](#) [32]

4.28.1 Member Data Documentation

4.28.1.1 char XRpad_VersionInfo::hwdriver[32]

4.28.1.2 char XRpad_VersionInfo::linux_kernel[32]

4.28.1.3 char XRpad_VersionInfo::msp_firmware[32]

4.28.1.4 char XRpad_VersionInfo::pld_firmware[32]

4.28.1.5 char XRpad_VersionInfo::software[32]

4.28.1.6 char XRpad_VersionInfo::spartan_firmware[32]

4.28.1.7 char XRpad_VersionInfo::subversion[256]

4.28.1.8 char XRpad_VersionInfo::wlan[32]

4.28.1.9 char XRpad_VersionInfo::xrpd[32]

4.28.1.10 char XRpad_VersionInfo::zynq_firmware[32]